INSTRUCTIONS FOR USE

# Kapitex<sup>®</sup>

## Laryngectomy Tube

**KAPITEX<sup>®</sup>**

# Kapitex®

## Laryngectomy Tube

<b>EN</b>	English.....	3
<b>FR</b>	Français .....	5
<b>DE</b>	Deutsch.....	7
<b>ES</b>	Español.....	10
<b>IT</b>	Italiano .....	12
<b>PT</b>	Português.....	14
<b>NL</b>	Nederlands .....	16
<b>DA</b>	Dansk .....	18
<b>SV</b>	Svenska .....	20
<b>EL</b>	Ελληνικά.....	22
<b>RU</b>	Русский.....	25
<b>PL</b>	Polski .....	28
<b>SL</b>	Slovensko .....	31
<b>ET</b>	Eesti .....	33
<b>LT</b>	Lietuvių .....	35
<b>TR</b>	Türkçe .....	37

## Kapitex<sup>®</sup> Laryngectomy Tube

**FEDERAL (US) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.**

### DESCRIPTION AND INTENDED USE

The Kapitex Laryngectomy Tube is a transparent, flanged, cuffless, curved silicone tube which is intended to be placed within a tracheostoma, in order to maintain stoma patency and a clear airway.

The exposed flange features 2 eyelet slots for attachment of a Twill Tape neck strap (supplied with the Laryngectomy Tube). The neck strap should be worn at all times when the Laryngectomy Tube is fitted, in order to help maintain the Laryngectomy Tube in position within the tracheostoma. The Kapitex Laryngectomy Tube is available in a range of different diameters and lengths, to account for variation in patient physiology.

The Kapitex Laryngectomy Tube is a **SINGLE PATIENT USE** device, indicated for use in adult and paediatric patients who have had a laryngectomy, who now breathe via a surgically created tracheostoma and who are not mechanically ventilated.

The maximum lifetime of a Kapitex Laryngectomy Tube is 28 days, at which time the device should be replaced with a new Laryngectomy Tube and neck strap and the old device disposed of.

Replacement of the Kapitex Laryngectomy Tube earlier than 28 days is recommended where the device shows signs of damage.

The Kapitex Laryngectomy Tube is supplied non-sterile, individually packed (together with a Twill Tape neck strap) in a polyethylene grip seal bag, which features a tamper-proof heat-seal.

### LARYNGECTOMY TUBE SIZING AND INITIAL FITTING

**NOTE: Initial fitting of the Kapitex Laryngectomy Tube and guidance on routine insertion and removal must be provided by an appropriately qualified clinician or healthcare professional.**

It is important that the Kapitex Laryngectomy Tube fits snugly within the patient's tracheostoma, without

exerting pressure on the tracheostoma wall, in order to provide adequate tracheal support. Correct sizing is based on the length, angle and diameter of the patient's stoma and, to obtain a proper fit, it is sometimes necessary to try a number of tube sizes.

### LARYNGECTOMY TUBE INSERTION

Before you insert the laryngectomy tube, you may wish to moisten the curved tube with a little clean water - this can help insertion of the device. Insert the Kapitex Laryngectomy Tube into the tracheostomal opening, curved end down, leaving the exposed flange slightly forward of the surface of the skin surrounding the tracheostoma opening. Check proper tube placement and secure the tube in position by attaching the Twill Tape neck strap to the eyelet slots on the front of the flange and placing the neck strap around the patient's neck.

### CLEANING

The Kapitex Laryngectomy Tube should be removed and cleaned daily to remove mucus and other secretions from its surfaces.

**NOTE: Assistance should be provided with removal, cleaning and re-fitting the Kapitex Laryngectomy Tube for patients who are unable to apply the device themselves.**

**Removal:** To remove the Kapitex Laryngectomy Tube, untie the neck strap, then grip the exposed flange of the Laryngectomy Tube between thumb and forefinger and gently pull the device out of the tracheostoma.

**Cleaning:** Soak the Laryngectomy Tube in a container of warm water for 10 minutes. Transfer the Laryngectomy Tube to a container of fresh warm water and physically clean with a specifically designed cleaning swab (Trachi-Swab - supplied separately) for 2 minutes. Clean inside and out with the swab, making physical contact with all surfaces using a backwards/ forwards/twisting motion. Rinse inside and out under cold running tap water for 1 minute. Shake off excess water from the device. Allow to air dry.

**Inspection and Refitting:** Prior to refitting the Kapitex Laryngectomy Tube, inspect it to ensure it is clean,

## ENGLISH

free from blockages and any signs of damage or deterioration. If inspection shows the Laryngectomy Tube is damaged or worn in any way, dispose of it and fit a new Laryngectomy Tube of the same size instead. Refit the neck strap to the front flange and re-insert the Laryngectomy Tube, following the instructions above for insertion of the device into the tracheostoma.

### DISPOSAL

In a healthcare environment, dispose of as clinical waste. In a home environment, dispose of as general household waste.

### WARNINGS AND PRECAUTIONS

- Avoid using a Laryngectomy Tube that is too small for the size of the patient's stoma, since this increases the risk of the device falling out of the stoma.
  - Do not apply excessive force during Laryngectomy Tube insertion, as this could result in trauma to the tissue of the stoma and/or the trachea.
  - Do not try to fit a Laryngectomy Tube that is too large for the stoma opening, since this may result in trauma to the tissue surrounding the stoma.
  - An incorrectly fitted or oversized Laryngectomy Tube may cause skin irritation, in which case the device should be removed immediately and the stoma allowed to recover before another Laryngectomy Tube is fitted.
  - Do not sterilize the Laryngectomy Tube using gamma irradiation, steam autoclave, boiling water, alcohol or other chemical agents.
  - Do not clean using any other method than that recommended for routine cleaning of the Laryngectomy Tube as described in this document. Use of other cleaning agents may result in damage to the Laryngectomy Tube cannula and flange, which could adversely impact the fit and/or the performance of the device.
  - Avoid aggressive or abrasive cleaning since this can lead to damage to the Laryngectomy Tube cannula and flange, which could adversely impact the fit and/or the performance of the device.
- Do not use an individual Laryngectomy Tube for more than 28 days.
  - Do not use the Laryngectomy Tube beyond the expiry date shown on the device packaging.
  - Do not use a Laryngectomy Tube if the polyethylene grip seal bag in which it is supplied is open or damaged, or the tamper-proof heat seal is not intact.
  - The Laryngectomy Tube must be inspected for signs of damage after cleaning and before reinsertion. Do not continue to use a Laryngectomy Tube that shows signs of damage, or wear, or build-up of mucus or secretions that cannot be removed by cleaning using the recommended cleaning method. Instead, dispose of the device and replace with a new Laryngectomy Tube.
  - Do not re-use a Kapitex Laryngectomy Tube on more than one patient, since this presents a risk of infection.
  - Do not use the Kapitex Laryngectomy Tube for any other purpose than the intended use specified in this document.

### CONTRAINDICATIONS

- Use of the Kapitex Laryngectomy Tube is contraindicated in patients who require mechanical ventilation.

**For further information or clarification, please contact your Kapitex distributor.**

Kapitex Healthcare Ltd,  
Unit 1, Erivan Park, Sandbeck Way, Wetherby,  
West Yorkshire, LS22 7DN, UK

**Tel:** +44 (0)1937 580211  
**Email:** sales@kapitex.com  
**Website:** www.kapitex.com

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Kapitex Healthcare Ltd, Unit 1, Erivan Park,