INHEALTH™

TECHNOLOGIES

Katz® Extractor

Oto-Rhino Foreign Body Remover

MEDICAL PROFESSIONAL Instructions For Use

R1

37593-01F

37593-01F Effective January 2016 / Gültig ab Januar 2016 /
Efectivo a partir de enero de 2016 / Date d'entrée en vigueur :
janvier 2016 / Valido da gennajo 2016 / Met ingang van januari
2016 / Vigente em janeiro de 2016 / V platnosti od ledna 2016 /
Wchodzi w życie w styczniu 2016 r.

Blom-Singer and **InHealth Technologies** are registered
trademarks of **Freudenberg Medical, LLC**

Patent(s): www.inhealthpatents.com

TABLE OF CONTENTS

Inhaltsverzeichnis / Índice / Sommaire / Sommario / Inhoud / Obsah / Spis treści

ENG	Katz® Extractor	4
DEU	Katz® Extraktor	8
SPA	Katz® Extractor	13
FRA	Extracteur Katz®	17
ITA	Katz® Extractor	21
NLD	Katz® Extractor	25
POR	Extrator Katz®	29
CES	Extraktor Katz®	33
POL	Ekstraktor Katz®	37


MEDICAL PROFESSIONAL Instructions For Use

Gebrauchsanweisung für MEDIZINISCHE FACHKRÄFTE
Instrucciones de uso para PROFESIONALES MÉDICOS
Instructions d'utilisation destinées aux PROFESSIONNELS DE SANTÉ
Istruzioni per l'uso destinate all'PERSONALE SANITARIO
Gebruiksaanwijzingen voor de MEDISCHE PROFESSIONAL
Instruções de uso para PROFISSIONAIS MÉDICOS
Návod k použití pro ZDRAVOTNICKÝ PERSONÁL
Instrukcja użytkowania dla PERSONELU MEDYCZNEGO

KATZ® EXTRACTOR

Oto-Rhino Foreign Body Remover

Developed in conjunction with David L. Katz, MD, MPH

Introduction

The Katz Extractor oto-rhino foreign body remover is a small, lightweight, single-use device for quick and easy extraction of foreign bodies from the nose or external ear canal.

This product is intended for single patient use only. Not for use between patients. See Warnings and Precautions.

Do not use the device if it becomes contaminated.

INDICATIONS

The Katz Extractor oto-rhino foreign body remover is designed and intended principally for extraction of impacted foreign bodies in the nasal passages. This device may also be used in the external ear canal, based upon clinical judgment.

CONTRAINDICATIONS

Do not use if the device appears to be damaged.

Do not use the device if bleeding is evident in the nasal or ear canal.

Do not use the device if you suspect there is not adequate space for the balloon to inflate. Increased pain would be an indicator of this.

HOW SUPPLIED

The Katz Extractor is supplied nonsterile; **do not sterilize as damage to the device may occur.** This device is designed as a single unit. Do not remove the catheter from the syringe. The Katz Extractor is available in a convenient, ready-to-use, five-pack dispenser box.

WARNINGS AND PRECAUTIONS

This device is intended for single use only. Use in the nasal canal or external ear canal only. Do not apply traction to a foreign body with sharp edges, as it may abrade or lacerate as it is displaced.

If, after inflation of the balloon, the impacted foreign body does not freely move out of the canal with minimal resistance (pull force), then deflate the balloon and remove the device. Repositioning of the catheter should be considered, as dictated by clinical judgment.

Insertion of this or any device into the ear canal may cause injury to the tympanic membrane; use of other methods, such as flushing, should be considered and medical judgment applied.

This device is designed as a single unit. Do not remove the catheter from the syringe.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

INSTRUCTIONS FOR USE

Note: The maximum inflation level of 1cc is monitored by a stop. Do not attempt to exceed the stop by use of force.

The following procedural instructions are provided by David L Katz, M.D., MPH. Please refer to the diagrams located at the back of this instruction manual.

1. Insert the catheter, with the distal balloon in the deflated state, along the canal until the catheter tip (balloon) has been positioned behind the impacted foreign body. Use the attached finger grip for additional stability, diagram 1a, 1b.

-
2. Inflate the balloon by depressing the syringe plunger until limited resistance from the canal wall is perceived, diagram 2a, 2b.
 3. Retract the catheter out of the canal with the balloon in the inflated state (plunger depressed), dislodging the foreign body in the process, diagram 3a, 3b.

Discard the device after use in accordance with CDC or institutional guidelines.

ORDERING INFORMATION

USA

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337, Monday–Friday, 9:30am–7:00pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: order@inhealth.com ORDER ON-LINE: www.inhealth.com POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

International

Please contact our customer service department for a distributor referral.

Consumer Affairs

If you have any questions or dissatisfaction with a product please contact our customer service department via telephone, fax, post, or Email: productcomplaints@inhealth.com

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