Blom-Singer®

voice restoration systems

Blom-Singer® Manometer with manometer adaptor

MEDICAL PROFESSIONAL Instructions For Use

INHEALTH®
TECHNOLOGIES

R1

37469-01E

37469-01E Effective March 2014 / Gültig ab März 2014 / Efectivo a partir de marzo de 2014 / Date d'entrée en vigueur : mars 2014 / Valido da marzo 2014 / Met ingang van maart 2014 / Vigente em de março de 2014 / V platnosti od března 2014 / Wchodzi w życie w marcu 2014 r.

Blom-Singer is a registered trademark of Hansa Medical Products, Inc.
InHealth Technologies is a registered trademark of Helix Medical, LLC.
© 2014 Helix Medical, LLC (InHealth Technologies)

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MANOMETER

with Manometer Adaptor

The manometer (diagram 1A), with manometer adaptor, is designed as a training aid for clinicians to teach new and current patients how to identify and self-regulate their own respiratory pressure requirements during tracheoesophageal speech. Proper airflow control results in better voice quality, less fatigue, and longer duration of the tracheostoma adhesive housing seal.

INDICATIONS

The manometer (with manometer adaptor) is for use by qualified clinicians as a training aid to teach laryngectomized patients how to recognize and self-regulate their intratracheal pressure during tracheoesophageal speech.

CONTRAINDICATIONS

The manometer (with manometer adaptor) is only indicated to assess and monitor intratracheal speech airflow pressure produced by laryngectomized patients who speak using a tracheoesophageal voice prosthesis.

HOW SUPPLIED

The manometer (with manometer adaptor) is supplied non-sterile in a kit. Replacement manometer adaptors are available separately. The manometer is reusable, however, the manometer adaptor is disposable and intended for **single patient use only**. **The manometer adaptor must not be reused by more than one patient to prevent cross contamination.**

Manometer with Adaptor	Manometer	Adaptor	TruSeal Housing
Manometer Kit	•	•	•
Replacement Manometer Adaptor Pack		•	•

INSTRUCTIONS FOR USE

Manometer Adaptor Attachment

Please refer to the diagrams located at the back of this instruction manual.

The clinician's hands should always be thoroughly cleaned to help avoid introducing contaminants.

The manometer adaptor (diagram 1B) should be connected to the manometer before attaching the manometer adaptor to the tracheostoma adhesive housing.

- 1.** Attach the manometer adaptor to the manometer by inserting the female luer lock connection (diagram 1C) of the manometer adaptor into the male luer lock connection (diagram 1D) of the manometer. Twist the two components together until securely connected.
- 2.** Attach the manometer adaptor to the patient by inserting the tracheostoma connector end to the tracheostoma adhesive housing (diagram 2).
- 3.** When the manometer adaptor is securely attached to the patient, and the manometer securely connected to the manometer adaptor, the clinician may proceed with the following assessment of intratracheal speech pressure.

Intratracheal Speech Pressure Assessment Procedure

The following procedure is provided by Eric D. Blom, Ph.D.

Failure to maintain the tracheostoma adhesive housing seal for a practical length of time is frequently the result of one or more of the following factors: (a) excessive intratracheal speech pressure; (b) excessive phlegm discharge; (c) ineffective valve seal application; or, (d) irregular tracheostomal architecture.

Excessive intratracheal speech pressure is a direct expression of resistance to airflow through a valved voice prosthesis, a tracheoesophageal vocal tract, or a combination of both. This "back pressure" pushes the tracheostoma adhesive housing seal loose from within, causing an air leak that can be durably repaired only by replacing the seal. Before initiation of therapy to teach a patient to attach and use a tracheostoma valve, a pressure measurement is recommended.

With the manometer adaptor inserted into the tracheostoma adhesive housing, which has been securely taped over the stoma, the patient takes a breath, lightly occludes the adaptor with a finger and counts to 15 at a normal conversational speech loudness level. Optimal levels are designated by color markings on the manometer: the green range indicates a pressure reading of 10 – 25 cm H₂O, which is excellent, and the blue range indicates a pressure reading of 26 – 40 cm H₂O, which is good. (See chart below.) Higher levels contribute to proportionally shorter duration of the tracheostoma adhesive housing seal.

Pressure Ranges (cm of H₂O)

10 — 25 cm H ₂ O	Excellent – Green
26 — 40 cm H ₂ O	Good – Blue
41 — 120 cm H ₂ O	High

For some patients, simply monitoring the manometer, as they vary their loudness from softest (usually registering 10 – 15 cm H₂O) to loudest (often measuring 90cm H₂O or more), teaches the importance of controlling their voice loudness. Many tracheoesophageal speakers attempt to speak too loudly resulting in excessive intratracheal “back pressure,” deteriorating voice quality, and shortened duration of the tracheostoma adhesive housing seal.

WARNINGS AND PRECAUTIONS

If there are tears, cracks, or structural damage, do not continue to use the product. The manometer adaptor is disposable and intended for single patient use only. The manometer adaptor must not be reused by more than one patient to prevent cross contamination.

The manometer is pre-calibrated at the factory. Tampering with the instrument, e.g., unscrewing the housing, may result in malfunction, inaccurate pressure readings and/or leaks of the gauge. Tampering with the connection fittings on the manometer adaptor may result in inaccurate pressure readings.

Caution: Federal (USA) law restricts this device to sale by or on the order of a qualified medical professional.

ORDERING INFORMATION

USA

Blom-Singer products may be ordered directly from InHealth Technologies. TELEPHONE: Toll-Free (800)477-5969 or (805)684-9337, Monday–Friday, 9:30am–7:00pm, Eastern Standard Time. FAX: Toll-Free (888)371-1530 or (805)684-8594. EMAIL: order@inhealth.com ORDER ON-LINE: www.inhealth.com POST: InHealth Technologies, 1110 Mark Avenue, Carpinteria, CA 93013-2918, USA, Attention: Customer Service

International

Please contact our customer service department for a distributor referral.

Consumer Affairs

If you have any questions or dissatisfaction with a product please contact our customer service department via telephone, fax, post, or Email: productcomplaints@inhealth.com

RETURNED GOODS POLICY

All returned merchandise must have a Return Merchandise Authorization number (RMA), and must be unopened and undamaged. RMA numbers are obtained from InHealth Customer Service. Returns without an RMA will not be accepted. Products must be returned in unopened packages, with manufacturer's tamper evident seals intact to be accepted for replacement or credit. Products will not be accepted for replacement or credit if they have been in possession of the customer for more than 14 days. Special order products are nonreturnable. There is a 20% restocking fee on all returned merchandise.

LIMITED WARRANTY, LIMITED REMEDY AND LIMITATION OF LIABILITY

InHealth Technologies, a Division of Helix Medical, LLC, warrants that the Product will be free from defects in workmanship and materials on the date of purchase. THE FOREGOING WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES OR REPRESENTATIONS, EXPRESS OR IMPLIED, ARISING BY LAW OR CUSTOM, INCLUDING IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR

PURPOSE, SAID IMPLIED WARRANTIES BEING HEREBY EXPRESSLY DISCLAIMED. Purchaser's sole and exclusive remedy for a valid warranty claim is replacement of the Product. This remedy does not include the cost of installation, removal, dismantling, delivery or reinstallation. InHealth shall have no liability to the extent Products are or have been modified or subject to misuse, abuse or improper storage, installation or maintenance or damaged in transit. No returns will be accepted without prior approval and authorization. INHEALTH SHALL NOT BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT, CONSEQUENTIAL, PUNITIVE, OR EXEMPLARY DAMAGES OF ANY KIND.

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