



Maximum
User Weight

Please make sure that you read and follow these instructions carefully. Failure to do so can result in injury. Keep these instructions in a safe place for future reference.




1.0 INTENDED USE

This medical device is designed to maximize mobility during your recovery from foot/ankle surgery or injury. Place the knee of your injured leg on the knee platform and stand as straight as possible. Adjust the handlebar to approximately waist high or the top of your hips. This will allow you to find the most comfortable position for your hands while holding the grips. Adjust the knee platform height so that your injured leg is supported at a 90 degree angle when standing.



Step 1

Put the knee walker as the picture shows, fold the steering column switch as shown in the figure.



Step 2

Unfold the knee walker, and insert the fold switch into the center of the connection, tighten and lock the fold switch.



Step 3

Install the front wheel: Connect the two iron pieces of the front wheel, insert the release switch, gasket and lock the screw.




Step 4

Insert the handle bars into the steering column. The screw and the release switch rotate at the same time until the release switch can be pressed tightly. (Note: The switch need to tighten, otherwise the steering column may still move.)




Step 5

Insert the knee pad into the tube and insert the screw into the hole, then tighten it. Insert the basket into the steering column and push down.



Brake and Parking Brake Operation

To engage the brake, squeeze the brake lever on the handle bar. To apply the parking brake, squeeze brake handle firmly and while squeezing the brake handle press the silver button to engage the lock (Top Right). While still holding the button down, release the brake handle and your parking brake will be set. To release the lock, squeeze the brake handle again and the silver button will pop up.



Parking Brake Unlocked

Parking Brake Locked



WARNINGS:



BEFORE OPERATING KNEE WALKER:

- ⚠ DO NOT exceed the maximum weight capacity of 300 pounds for this knee walker. Exceeding the maximum weight capacity may result in injury to the user or damage to the knee walker.
- ⚠ After fully assembling the knee walker according to this user manual, ensure that the brakes are working correctly, meaning they engage the rear wheel(s) when the brake handle is squeezed. Failure to do this may result in injury to the user.
- ⚠ Inspect the knee walker periodically to ensure the brakes and wheels are functioning properly. Be sure steering column clamp and height adjustment knob are tightened and fully engaged before each use. Also, ensure that the knee platform height adjustment clamp is tightened and fully engaged and the locking pin is locked in place prior to use.