

care and cleaning of equipment

Thank you for purchasing the DMA Clinical Reformer. Designed for years of safe trouble free use they will need only some basic care to ensure all parts continue to work smoothly.

Over time dust and dirt may build up on tracks in affecting the smoothness of the carriage running.

Once / twice a month wipe tracks and running wheels on carriage with methylated spirits.

All wheel bearings are sealed, do not spray with lubricants.

DO NOT use methylated spirits on vinyl top, damp soapy rag will clean best.

To ensure optimum operation of the footbar mechanism keep bolts between frame and footbar firm **(A)**.

Chrome locking bar and upright bolts **(B)** should be slightly loose so that the locking mechanism locks reliably. A slight spray of lubricant may help occasionally.

contents

Reformer Frame
Carriage
Shoulder Pads
5 Springs
Instruction Booklet

warranty

Clinical Pilates P/L designed and manufactured exercise equipment is guaranteed against defects in workmanship and materials under normal use for 24 months. This includes ropes and springs.

If any part thereof proves to be defective within the warranty period it will be repaired or replaced, free of charge, and at the discretion of Clinical Pilates Pty. Ltd. This includes parts and labour.

This warranty is extended to the original purchaser from date of purchase.

- Damage to structural or other integral parts of equipment due to negligence is not covered by warranty, in such cases warranty coverage will be subject to the manufacturers' discretion.
- The warranty will be voided if the equipment has been modified, altered or contains components that are not part of the Clinical Pilates Pty Ltd line.

The benefits conferred by this warranty are in addition to all other rights and remedies which the consumer has under the Trade Practices Act.

to obtain warranty service

Call Clinical Pilates P/L and fax / email your written claim with details of the problem experienced, Appropriate action will be taken as quickly as possible.



1 Yarra Street, South Yarra VIC 3141 Australia
p 61 3 9827 4511 / 61 3 9525 1566
e info@clinicalpilates.com
w www.clinicalpilates.com



reformer assembly

reformer assembly instructions

1. Attach springs under carriage with light spring (gold end) in centre, feed through webbing slings.
2. Insert shoulder pads into brackets, offset to outside.
3. Wipe tracks and wheels down thoroughly with damp rag to remove any residue from packaging.
4. Place carriage on tracks and run to check for smooth running. Wipe / check wheels if any bumps.
5. To attach ropes place foot loops around shoulder pads feed through rear pulleys and forward through cleats. Drop remainder through hole in carriage.

