

FOLDING WALKER

QUICK START GUIDE

WHAT'S IN THE PACKAGE

- 1x Folding Walker
- 1x Removable Walker Bag
- 4x Adjustable Leg Sections (pre-installed)
- 2x Adjustable Wheel Legs

OVERVIEW

We are constantly answering questions and recording helpful videos to make using your Vive Folding Walker as easy as possible. Throughout the manual we will include QR codes and links to important videos that you can visit for help.



Check out our list of Frequently Asked Questions at vhealth.link/jpr for helpful answers.

OPERATING INSTRUCTIONS

Product is fully assembled and ready to use.

(This item will require no tool to assemble.)



For a video tutorial on folding the walker, check out vhealth.link/z3e

To begin folding the walker, simply press down on each rubber tab location on the center support bar to unlock each side. Then pull the back of each side in towards the center of the walker to fold the walker.

To adjust the length of the legs, press down on the rounded metal button on the leg to unlock the leg section. Then push/pull the leg to the desired length, and align the depressed button with the adjustment hole corresponding to the desired length. The button will click back into place and lock the leg in the desired height. Repeat with all 4 legs and adjust to an even height that is stable and comfortable to use.

Walking with the Walker

- Align the middle of your foot with the back legs of the walker.
- Lift the walker up and place it a step's distance in front of you, making sure all 4 legs are stable on the ground.
- Carefully walk forward one step, and repeat.

SPECIFICATIONS

- Width - 20"
- Depth - 16'
- Adjustable height - 32" - 38"
- Product Weight - 6.5 LBS
- Max User Weight - 250 LBS

WARNINGS

- Do not push one side of the walker ahead and then other side without lifting the walker first.
- Use Vive accessories only. Use of other manufacturers' accessories can result in improper fit and risk fall and/or injury.
- Do not use on stairs or escalators.

vive®