

EXTENDED TRIGGER FINGER SPLINT

QUICK START GUIDE

WHAT'S IN THE PACKAGE

- 1x Full Finger Splint
- 1x Long Palm Strap
- 1x Short Finger Strap
- 1x Wrist Strap



PRODUCT OVERVIEW

We are constantly answering questions and recording helpful videos to make using your Finger Splint as easy as possible.



To see all of the FAQ's in one place visit vhealth.link/8ec14

HOW TO PUT IT ON

To get the Splint safely secured onto an injured finger, it is recommended that you have the assistance of another person. With that, just complete the following steps:



Want to watch us do this instead? Just go to the following link to see how to properly put on your Finger Splint. vhealth.link/501b2

1. Place the hand over the splint so that the injured finger lays over the splint. Secure in place with the top finger strap.



2. Adjust the wrist strap using the small velcro tab. Position it so that it is able to comfortably wrap around your wrist.



3. Loop the wrist strap through the D-Ring and adjust so that it fits comfortably.



4. Adjust the double-sided velcro tab on the Long Hand Strap along the Splint around the center of the palm. Wrap the Hand Strap around the center of the hand and secure the opposite end to the other side of the double-sided tab.



5. If extra support is needed, use the Short Finger Strap to further secure the injured finger to the Splint wherever you feel is comfortable.



CARE INSTRUCTIONS

- To clean the Splint and the various Straps, hand wash in cold water with a mild soap, rinse thoroughly, and let them air dry. Make sure to keep the Splint dry while wearing.

- The Splint is NOT machine-washable and should not be machine-dried either.
- The fabric is latex-free, and the splint inside is aluminum to resist corrosion.

vive®