

# WALKER WHEELS + SKI GLIDES WITH COVERS

## QUICK START GUIDE

### OVERVIEW

We are constantly answering questions and recording helpful videos to make using your Vive Walker Wheels + Ski Glides with Covers as easy as possible.




To see all FAQs in one place visit  
[vhealth.link/cac19](https://vhealth.link/cac19)

### USER INSTRUCTIONS


#### Note:

- For walker without front wheels, install walker glides on front 2 legs.
  - For walker with front wheels, install walker glides in the rear 2 legs.
1. Flip walker over and remove the front or rear two rubber feet. See note above.

2. Insert glides on each of the front or back legs. The arrows on the glides should be pointing towards the direction of forward movement. See Figure 1.



3. Use Phillips screw driver to expand or retract diameter of insert. Tighten screw until glides are secured to the walker. See Figure 2.



**Covers:**

Place the included covers over each glide if using on wooden, laminated, or tiled floors. Covers may be hand washed in cold water and air dry.

**To Remove:**

To remove glides simply loosen screw until glides are loose enough to come off. Remove and store away.

**Warning:**

**These glides are designed to reduce friction, movement may feel slightly different for a first time user.**

**Walker Wheels:**

1. Remove the front legs of the walker.
2. Slide new wheels on the front legs of the walker.

The logo for 'vive' is displayed in a white, lowercase, sans-serif font against a teal background. A registered trademark symbol (®) is located to the upper right of the letter 'e'.