

# MOBILITY SIDE BAG

## QUICK START GUIDE

### OVERVIEW

We are constantly answering questions and recording helpful videos to make using your Mobility Side Bag as easy as possible. Check out the included link and QR code to help you through the process.



For answers to all your FAQs  
check out [vhealth.link/h85](https://vhealth.link/h85)

The Vive mobility side bag comes equipped with 2 top straps, 1 bottom strap, 1 hand strap, and 1 back column strap to make it convenient and adaptable to different types of mobility devices such as mobility scooters, wheelchairs, and walkers.

### ATTACHMENT

Use the top 2 straps to secure onto the armrest handles of your mobility scooter, wheelchair, transport chair, or any other mobility device which uses an armrest.



In addition to the top 2 straps, you may use the bottom strap to further secure your side bag onto your mobility device if there are vertical supports to strap them onto, such as the side of a folding walker. This strap is optional and most often will not be needed.



Use the back column strap to secure your side bag onto a single support such as a rollator or the front of a scooter.

vive®