- EN** Instructions for use of **Aniball INCO** medical device
- DE** Gebrauchsanweisung für das medizinische Produkt **Aniball INCO**
- FR** Instructions pour l'utilisation d'**Aniball INCO** medical device
- IT** Istruzioni per l'uso di **Aniball INCO**
- ES** Instrucciones de uso del dispositivo médico **Aniball INCO**
- PT** Instruções de utilização do dispositivo médico **Aniball INCO**
- NL** Instructies voor het gebruik van **Aniball INCO** medisch hulpmiddel
- FI** Käyttöohjeet lääkinälliselle **Aniball INCO - laitteelle**
- SE** Bruksanvisning för den medicintekniska produkten **Aniball INCO**
- DK** Brugsanvisning til det medicinske udstyr **Aniball INCO**
- RO** Instrucțiuni de utilizare al dispozitivului medical **Aniball INCO**
- BG** Указания за използване медицинско устройство **Aniball INCO**
- GR** Οδηγίες χρήσης της ιατρικής συσκευής **Aniball INCO**
- HR** Upute za korištenje medicinskog proizvoda **Aniball INCO**
- SI** Navodilo za uporabo medicinskega pripomočka **Aniball INCO**
- HU** Útmutató az **Aniball INCO** gyógyászati segédeszköz használatához
- PL** Instrukcja użycia urządzenia medycznego **Aniball INCO**
- SK** Návod na používanie zdravotníckej pomôcky **Aniball INCO**
- CZ** Návod na používání zdravotnického prostředku **Aniball INCO**

AniBall®
INCO

EN CONTENT

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
The page features a decorative background consisting of numerous small red dots scattered across the white surface. Additionally, there are four registration marks, each consisting of a circle with a crosshair inside, positioned at the top, bottom, left, and right centers of the page. The text is centered and reads:

Instructions
for use of **Aniball INCO**
medical device

EN 1. INTENDED PURPOSE

DE

Aniball INCO is a medical device intended for the prevention and treatment of stress urinary incontinence.

FR

IT

Properties

ES

Aniball INCO is a medical device (hereinafter also referred to as "balloon") for the training of conscious relaxation and activation (contraction) of the pelvic floor muscles of an adult woman. Balloon will help to discover the pelvic floor and find out how to work with it effectively in the context of the breath.

PT

NL

FI

Aniball INCO is used to train the proper function (activation/relaxation) of the pelvic floor muscles for the prevention and treatment of stress incontinence. The balloon inserted into the vagina is easily inflated with air to the desired comfortable size. Thanks to its variable shape and size, it adapts optimally to the individual anatomy of each woman. The balloon inserted in this way must not cause pain or other discomfort. In case the balloon pressure is perceived unpleasantly or painfully, it is possible to reduce it at any time by releasing the air valve on the pump body. However, light balloon pressure is desirable to feel its position and resistance during exercise. By activating and, conversely, releasing the pelvic floor muscles, there is a slight movement of the balloon inside the vagina (inwards of the vagina and outwards), which provides very important feedback on the quality of work of the pelvic floor muscles.

SE

DK

RO

BG

GR

HR

SI

Aniball INCO consists of:


1. Inflatable balloon made of silicone
2. Mechanical pump
3. Tube
4. Connecting adapter
5. Two-way valve

HU

PL

SK

CZ

Please read the following instructions carefully!

2. CONTRAINDICATIONS OF THE USE

The use of the Aniball INCO balloon is limited by the woman's health. Consult a physician, physiotherapist, or midwife to assess the suitability of the exercise and to consider contraindications.

Do not use Aniball INCO in the following cases:

- pregnancy
- any vaginal inflammation (infectious or non-infectious)
- any vaginal bleeding
- presence of condyloma or other infectious pathology in the vagina or vaginal entrance area
- vaginal injury
- any bladder inflammation
- uterine prolapse
- genital herpes
- precancers and tumour pathologies of the cervix, vagina or external genitalia
- during the postnatal puerperium and within the first six weeks following gynaecological interventions and operations (including genital laser treatment)

In the following cases, you need to be extra cautious and consult your doctor:

- presence of varicose veins in the vagina and external genitalia
- reduced sensitivity of the vagina or external genitalia (neurological disease, painkillers administration, alcohol and drug abuse)
- after treated vaginal inflammation (vaginal fragility persists)
- after vaginal surgery (scarring may negatively affect the vaginal elasticity)
- lichen sclerosus of the external genitalia (chronic skin disease)
- cancer in the pelvic area

Any serious incident that has occurred in relation to the Aniball INCO should be reported to the manufacturer and the competent authority of the Member State in which the user is established.

EN

DE

FR

IT

ES

PT

NL

FI

SE

DK

RO

BG

GR

HR

SI

HU

PL

SK

CZ

EN 3. SAFETY PRECAUTIONS FOR THE USE

DE

FR

IT

ES

PT

NL

FI

SE

DK

RO

BG

GR

HR

SI

HU

PL

SK

CZ

- Before the first, and after each subsequent use, wash the silicone balloon thoroughly using oil-free pH-neutral soap and warm water – at least 40 °C / 104 °F – for 1–2 minutes and rinse under running water. Subsequently, it is advisable to disinfect the balloon using a product suitable for silicone. Perform the disinfection according to the manufacturer’s instructions. Or boil the balloon in boiling water for at least 5 minutes. Then, allow it to dry on a clean and smooth surface. You may then place the balloon onto the plastic adapter. The balloon fits to the adapter relatively tightly, but it can be disconnected and reconnected at any time.

- After each use, the balloon must be thoroughly washed with pH neutral oil-free soap, rinsed under warm running water, disinfected with a suitable product and stored in a (the original) bag.
- If pauses between exercise sessions exceed 5 days, wash and disinfect or boil the balloon again before the next use.
- If in doubt about choosing a suitable soap and disinfectant, visit our website: www.myaniball.com
- Before inserting the balloon into the vagina, it is necessary to check its integrity and correct functionality by inflating it using the pump to 20 pumps. **Perform this check before each use! Never use the product in the case of any damage or signs of material degradation!**
- For comfortable and safe insertion, it is advisable to lubricate the surface of the balloon using **exclusively water-based lubricants suitable for internal use (contact with mucous membranes). Do not use products or condoms containing oils.** Non-water-based lubricants may damage the silicone, causing the balloon to burst. Please be aware that such damage will not be considered a product defect in the case of a warranty claim.
- If you practice any kind of pelvic massage before exercising your pelvic floor muscles, **use a suitable water-based lubricant instead of massage oil.**

- Prevent contact between the pump and water when washing the balloon, as water could prevent its proper functioning.
- Always let the balloon dry thoroughly after each wash.
- Store the device separately in the original bag in a dry place and do not expose it to radiant heat or sunlight!
- Prevent the device from making contact with sharp objects.
- Do not use in water.
- Keep the device out of reach of children.
- **For safety and hygienic reasons, the device is intended for use by one person only!**

Recommendations for Using Aniball INCO

- In order to use this device effectively and reap the rewards of your efforts, you need to exercise regularly and with persistence. **Regular exercise** (every day for 5–15 minutes) alleviates the difficulties within approximately three months.
- **During exercise, never hold your breath! (Except before the intentional coughing during the pressure resistance training).**
- During exercise, it is a common mistake **to pull the belly button** into the abdomen or **to pull the gluteus muscles** towards you. This is an inappropriate substitute that prevents the true deep activation of your pelvic muscles.
- Sometimes, you may feel fatigue in your pelvic muscles. That's perfectly all right; it indicates you are exercising correctly. If the fatigue is so great that the muscles cannot be activated, stop the exercise.
- **Do not try to work through the pain nor intentionally induce pain while exercising. Pain is not a sign of correct exercise. What you should feel is a pleasant pull from the stretching or the muscle activation.**
- If you feel a sharp or disproportionate pressure and tension in your lower abdomen, take a break or stop exercising.
- **The balloon should never slip out of the vagina too quickly, so hold it by the plastic neck.** A quick slip indicates a loss of control and could result in injury.
- **Never allow an inserted and simultaneously inflated balloon to be expelled directly out of the vagina during exercise. We always recommend that you release it to the minimum size using the valve before removing it.** Along with the breath, relax the pelvic floor muscles, grasp the balloon by the plastic neck and slowly pull it out of the vagina, or slide it outwards so that it slides out spontaneously.

EN

DE

FR

IT

ES

PT

NL

FI

SE

DK

RO

BG

GR

HR

SI

HU

PL

SK

CZ

EN

DE

FR

IT

ES

PT

NL

FI

SE

DK

RO

BG

GR

HR

SI

HU

PL

SK

CZ

- If you have decided to use **Aniball INCO** after giving birth, keep in mind that you can **exercise at the earliest after the end of the sixth week and any birth injuries have healed**.
- Learn to control your pelvic floor muscles during the day, even when not exercising with the balloon. Before coughing or sneezing, imagine you are squeezing the balloon and preventing it from sliding out of the vagina. This activation will help you to regulate urine leakage in a strenuous situation. **Never “strengthen” your pelvic floor by interrupting urination! This incorrect procedure increases the risk of urinary tract infection, and may lead to disruption of the micturition reflex which contributes to incontinence!**

If you have any questions, problems or suggestions, please do not hesitate to contact us at: hello@myaniball.com. We are always here for you and ready to help.

4. EXERCISE INSTRUCTIONS

You should exercise at least once a day for 5–15 minutes. All exercises should be performed calmly, slowly and as closely to these instructions as possible. Exercise your pelvic floor muscles with an empty bladder.

1. Body Position During Exercise

Start with a position lying down on your back, bend your knees and anchor your feet on the ground. Place your feet at hipwidth distance apart. If you feel your head is leaning back, place a pillow or a folded towel under it, so the cervical spine is straight, without a backward bend. Your back should rest comfortably on the mat without a bend or arch in your lower back. Leave your hands loosely alongside your body or place them comfortably on your lower abdomen.


2. Inserting the Balloon

Before the balloon is inserted, relax and do not rush. The vagina may be dry and sensitive; therefore, we recommend **applying appropriate lubricant** (see **Safety Precautions for Use**). You may leave the balloon deflated or slightly inflate it for 1–2 pumps. Then, grab the balloon behind the plastic neck and insert it into the vagina. If you have inserted the balloon inflated and now feel uncomfortable pressure inside your vagina, you can reduce its size using the valve at the end of the pump to achieve a comfortable volume. If you

have inserted the balloon but do not feel its slight resistance, we recommend inflating it to a size which is comfortable for you but also makes you feel slight resistance of the balloon against the vaginal walls. Once you reach optimal feeling without uncomfortable pressure or pain, but at the same time a slight resistance of the balloon, you can start with the exercise. Perform the exercises below with respect to your own body.

If you are tired or have a headache, take a break at any time during the exercise or stop with the exercise for that day. You must not feel pain during any of these exercises.

3. Breathing Exercise

In this exercise, you will learn to activate and release the pelvic floor muscles in coordination with your breathing. First, slow down your breath and be aware of individual inhales and exhales.

While **inhaling**, imagine your breath flowing through your chest down to your abdomen and lower abdomen. Visualize the inhaled air reaching the pelvic floor muscles, all the way to the place of the inserted balloon. If your palms are placed on your abdomen or groin, you can feel the growing inner volume in this area which aims to gently push out the balloon. While inhaling, teach your pelvic floor muscles to consciously relax in order to feel the balloon moving under the gentle pressure in the outward direction. However, you should not loosen your muscles too much to avoid expelling the balloon completely out of the vagina (the goal is a conscious and controlled release).

EN

DE

FR

IT

ES

PT

NL

FI

SE

DK

RO

BG

GR

HR

SI

HU

PL

SK

CZ

**EN**

Upon each **exhale**, learn to firmly grip the balloon and move it upward deeper inside the vagina (a conscious activation of the pelvic muscles). At the same time, the internal volume of the abdomen created by the previous inhale may shrink again. It is very useful to visualize the balloon squeezed by the vaginal walls, not by the external urethral or rectal sphincters. Try not to squeeze your buttock muscles when trying to hug the balloon, nor to draw your belly button into your belly. **While exhaling, leave the abdominal area relaxed without clenching it** (learn to use only your pelvic muscles).

DE

FR

IT

ES

PT

NL

FI

SE

DK

4. Prolonging the Grip Time

During this exercise, you will increase the strength and endurance of your pelvic muscles.

RO

Perform the exercise in the same way as the previous one (see **Breathing Exercise**) with one difference in the activation stage (grip of the balloon), where you consciously prolong the grip of the balloon to hold it for at least ten seconds (you may count to ten). While holding the balloon, try either to do long and continuous exhales, or to inhale again and exhale (but, if possible, try it without losing the grip of the balloon). Repeat the entire cycle of continuous activation 5–10 times.

BG

GR

HR

SI

5. Resistance Against Increased Pressure

The aim of this exercise is to manage situations of excessive abdominal pressure, which can lead to unwanted release of the pelvic floor or urine leakage.

HU

PL

SK

CZ

Take the position of lying on your back as in previous exercises. While inhaling, lead your breath to the lower abdomen and groin to increase the volume of this area. Now try squeezing the balloon upon inhaling. Exhale and try to not lose the obtained volume and consciously keep the balloon squeezed using your pelvic muscles. Repeat 5–10 times.

Continue with exhaling while trying to cough at the same time, without letting the balloon slide down or all the way out of the vagina. Repeat 3–5 times.

6. Relaxation

Get ready for the most popular and calming part of your exercise dedicated to relaxing and loosening your pelvic muscles after their previous strengthening. You can deflate the balloon first and remove it carefully from the vagina (see Removing the Balloon). Alternatively, you can leave the balloon inside the vagina; this is not wrong. With the balloon still inserted, you will feel the relaxing of your pelvic muscles more intensely; however, do not let the inserted and inflated balloon slip all the way out of the vagina during the exercise. You should remove the balloon, or let it slip out freely, always while deflated!

Stay on your back, pull your knees up to your chest with your hands, so your feet no longer touch the floor, and assume the position of a "curled ball". Let your head rest on a thin pillow or folded towel. Now focus on the inhaling phase and imagine the entire area of your pelvic floor expanding and extending in all directions, like a blooming flower. While exhaling, stop squeezing the balloon, and simply relax. Repeat the inhale and let your pelvic floor once again open up like a flower with no signs of a grip or forced squeezing. Dedicate a lot of time for the relaxing exercise in order to give your muscles enough time to relax. Therefore, repeat the entire inhale-exhale cycle 15–30 times.

7. Removing the Balloon

Before removing the balloon from the vagina, we recommend deflating it to a minimum size using the valve. Inhale and relax your pelvic floor muscles, then take the balloon by the plastic neck and slowly pull it out of the vagina. Alternatively, slide it outward, so it slips out on its own. Continue with washing the balloon, see **Safety Precautions for Use.**

EN

DE

FR

IT

ES

PT

NL

FI

SE

DK

RO

BG

GR

HR

SI

HU

PL

SK

CZ

EN

You can find more information at: www.myaniball.com

DE

Shelf life: indicated on the medical device label and packaging.

FR


Prior to disposal, clean and disinfect in accordance with the instructions for use!

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