



## Introduction




The low dead volume SilFlow® GC capillary column 3- or 4-port splitter system enables you to install and remove a column without using a wrench; giving you flexibility in the GC irrespective of your configuration.

## Packing list



	P/N 123720 SilFlow 3 port splitter GC kit (1.1)	P/N 123721 SilFlow 3 port splitter GC kit (0.53)	P/N 123722 SilFlow 3 port splitter GC kit (0.25)	P/N 123730 SilFlow 4 port splitter GC kit (1.1)	P/N 123731 SilFlow 4 port splitter GC kit (0.53)	P/N 123732 SilFlow 4 port splitter GC kit (0.25)
	Qty	Qty	Qty	Qty	Qty	Qty
MultiChannel Device (MCD)	1	1	1	1	1	1
FingerTite tool	1	1	1	1	1	1
Mounting bracket & screw	1	1	1	1	1	1
SilFlow FingerTite nuts	4	4	4	5	5	5
SilFlow FingerTite ferrules (1.1 mm)	10	-	-	10	-	-
SilFlow FingerTite ferrules (0.4 mm)	10	10	10	10	10	10
SilFlow FingerTite ferrules (0.5 mm)	10	10	10	10	10	10
SilFlow FingerTite ferrules (0.7 mm)	-	10	-	-	10	-
Blanking pins	2	2	2	2	2	2
SilFlow FingerTite jig (0.4 mm)	1	1	1	1	1	1
SilFlow FingerTite jig (0.5 mm)	1	1	1	1	1	1
SilFlow FingerTite jig (0.7 mm)	-	1	-	-	1	-

## Instructions

### Column splitter installation



1. Attach the mounting bracket (1) to a suitable point in the GC Oven – ensure the locator pins (2) can be easily accessed to install the SiFlow MCD.
2. Locate the SiFlow MCD (3) and identify the appropriate port location suitable for the GC capillary column connection (see Figure 2a and 2b – the port dimensions are located on the back of the MCD).




3. Locate the appropriate SiFlow FingerTite ferrule (4), SiFlow FingerTite nut (5), capillary or blanking pin (6) and FingerTite tool (7)

### Assembly with FingerTite jig

The FingerTite jig is designed to pre-swage the ferrule for SiFlow devices outside of the oven. Pre-swaging the ferrule enables you to use the connection




onto any SiFlow port suitable for the OD of the fused silica capillary and ensures the capillary will not be crushed into the SiFlow channels.



There are 3 pre-swaging jigs; 0.4, 0.5 and 0.7.

When the ferrule is tightened into the jig using the FingerTite tool, (never use a wrench to tighten SiFlow FingerTite ferrules), the ferrule is fixed on a position of the fused silica capillary at preset offset (see Figure 4).



4. Once all the capillaries are connected, mount the SiFlow MCD on the mounting bracket (if it is not already mounted).

For application configuration suggestions visit [www.trajanscimed.com/silflow](http://www.trajanscimed.com/silflow)

### Information and support

Visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)

Specifications are subject to change without notice.