# Biopsy pads

<b>Intended use</b>	For in vitro diagnostic use. Biopsy pads are support material intended for holding tissue biopsies in place.
<b>Product</b>	Tissue biopsy pads (blue or black). Provided in bags – each containing 500 or 1000 pads.
<b>Precautions</b>	<ul style="list-style-type: none"> <li>• Single use – do not reuse.</li> <li>• Store in a dry place.</li> <li>• Do not store in direct sunlight</li> <li>• For laboratory use only by trained personnel and professional use.</li> <li>• Some leaching of dye may occur.</li> </ul>
<b>Storage instructions</b>	<ul style="list-style-type: none"> <li>• Store the product in original packaging, in a dry environment, preferably protected from direct sources of heat and light.</li> <li>• Store at room temperature – Recommended storage temperature between 2°C and 30°C.</li> <li>• Avoid large variations in temperature conditions during storage.</li> <li>• Do not use after use-before date.</li> <li>• Shelf-life is indicated by the use-before date, printed on outer packaging.</li> <li>• If tearing of biopsy pads occurs, or if suboptimal performance is noticed and cannot be attributed to variations in laboratory procedures, please contact your regional Trajan representative. Some leaching of dye may occur during processing.</li> <li>• Trajan biopsy pads resist temperatures from -40°C to 121°C, and are autoclavable.</li> </ul>
<b>Usage</b>	<ul style="list-style-type: none"> <li>• Trajan biopsy pads should be gently separated prior to use. They are compatible with all routine sized tissue cassettes used in histopathology laboratories. The pads prevent tissue from being lost during processing.</li> <li>• To avoid artefacts, biopsy pads should be pre-soaked with 10% neutral buffered formalin or a fixative of choice before use.</li> <li>• One or two biopsy pads should be used to line the tissue cassette and secure the specimen. The cassette can be secured with either a plastic or metal lid.</li> <li>• When the tissue processing is complete, carefully open the cassette and extract tissue whilst separating the pads gently.</li> <li>• After tissue is extracted, dispose of biopsy pads as per laboratory and local regulations.</li> </ul>
<b>Compatible products</b>	Trajan biopsy pads are compatible with most tissue cassettes used in histopathology laboratories.
<b>Compliance certification</b>	 

Technical data	
Dimensions (length x width)	31 mm x 26 mm
Thickness	2.5 mm
Material	100% reticulated polyester urethane

For more information about this product visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)

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## **Europe, Middle East, Africa and India**

Trajan Scientific Europe Ltd  
Crownhill Business Centre  
14 Vincent Avenue  
Crownhill, Milton Keynes  
MK8 0AB, United Kingdom  
Tel: +44 (0) 1908 568 844  
Fax: +44 (0) 1908 566 790  
europe@trajanscimed.com

## **Asia**

Trajan Scientific Japan Inc  
RK CUBE 3F, 6-85 Otamachi  
Naka-ku, Yokohama-shi  
Kanagawa 231-0011, Japan  
Tel: +81 (0) 45 222 2885  
Fax: +81 (0) 45 222 2887  
japan@trajanscimed.com

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Trajan Scientific Americas Inc