## Abstract

The demand for separation techniques for intact proteins is increasing with the introduction of a new generation of high resolution mass spectrometers which are able to measure the mass of small to medium size proteins very accurately. Liquid chromatography is a valuable tool for separating these proteins prior to the MS analysis. Intact protein chromatography is most commonly used in a top-down approach in proteomics and to determine expression levels during recombinant protein synthesis.

The size of the protein molecule results in very low diffusion coefficients and therefore slow mass transfer in and out of the pore system. A sufficiently large pore diameter is required to minimise the effects of restricted pore diffusion. We show examples of the separation of intact proteins on a column packed with 3 µm C8 silica with 1000 Å pore size. The molecular weight of the protein examples reported here cover ribosomal proteins (<40 kDa), monoclonal antibodies (~150 kDa) and intact membrane proteins derived from mouse liver.

Physico-Chemical Background **Diffusion Inside Pores** Diffusion in Free Solution |1-2.104| $D_n = D_f$ The diffusion rate is described by:  $D = \frac{n B^2}{2}$ **\_\_** \_\_ \_  $6\pi\mu r$ Renkin Equation (E.M. Renkin, J.Gen.Physio., 38 (1954) 225.) **Einstein-Stokes Equation**  $D_{f}$  = Diffusion rate constant in free solution  $D_n = Diffusion$  coefficient inside the pores k<sub>B</sub> = Boltzmann constant  $r_s = Stokes$  radius of the analyte T = Temperaturer<sub>∞</sub> = Pore radius r = Stokes radius  $\mu = Viscosity$ Fast diffusion rates but low capacities -0.9 Pores are larger than needed. 100% kinetics and high capacity. Size[Å] Compound | 3.89E-10 Caffeine 7.91E-11 Insulin **~** 0.5 2.75E-11 hSA 04 120Å - 5kDa 170 1.44E-11 IgG 🔵 300Å - 100kDa 🔰 100Å - 5kDa 20.3% 7 1 0 / /.1% 3.7% 0.2 0.1 Caffeine

Figure 1: Normalized diffusion rates of small, medium and large molecules in free solution

# Dedicated Protein Columns - C8, 3 µm, 1000Å

The SGE column range for biological analytes focuses on an optimized pore size for both peptides (200 Å and 300 Å) and proteins (1000 Å) and small particle size (3 µm) to enable lowest possible mass transfer restrictions. Exposure to the analyte to metal surfaces is minimized through the use of glass-lined or PEEK<sup>TM</sup>-lined stainless steel tubing. All tubing is PEEK<sup>TM</sup> coated fused silica. Great emphasis is put





Figure 3: Design of the ProteCol<sup>™</sup> Range Columns



AUSTRALIA & PACIFIC REGION

# Analysis of Intact Proteins Using Liquid Chromatography



"Sweet spot", where diffusion rates are between 50% and 80% of the free diffusion. Good compromise between fast diffusion kinetics and high capacity

> Slow diffusion rates - high capacity -Pores are too small for effective mass transfer.

#### 💊 200Å - 100kDa

120Å - 100kDa 100Å - 100kDa

0.4 0.5 0.6 0.7 0.8 0.9 Analyte diameter/pore diameter

Figure 2: Relative diffusion rates in free solution and inside pores as function of molecule size to pore

# **Accurate Mass Analysis of Intact Ribosomal Proteins**

## Introduction

Since ribosomal proteins are relatively small (Mw=6,000 to 40,000), they can be rapidly identified by accurate LC MS analysis of the intact proteins. In the present application ribosomal proteins isolated from rat liver were separated on a ProteCol<sup>™</sup> C8 HQ1003 column.

### **Sample Preparation**

80S ribosomal proteins were isolated from a rat liver microsomal preparation (Williamson et al; 1997, Eur. J. Biochem. 246: 786-793). One optical density unit at 260 nm of 80S ribosomal proteins was mixed with two volumes of 6M Guanidine HCl to denature the proteins. 1 % (v/v) formic acid was subsequently added to precipitate the nucleic acids. The mixture was centrifuged for 15 min at 13,000 rpm and the supernatant was collected into a sample vial ready for LCMS analysis.



Figure 4: Base peak chromatogram of ribosomal proteins

#### **Data Analysis and Results**

All data were acquired and reference mass corrected via a dual-spray electrospray ionisation (ESI) source. Each scan or data point on the Total Ion Chromatogram (TIC) is an average of 15,000 transients, producing a spectrum every second. Mass spectra were created by averaging the scans across each peak and background subtracted against the first 10 seconds of the TIC. The resulting base peak chromatogram shows very high peak capacity - 119 discrete protein masses were identified; 46 of which were identified as 80S ribosomal proteins. In some cases several different masses of the same protein were identified which correlated with known N- and/or C-terminal processing.

## **Accurate Mass Determination of Intact Monoclonal Antibodies**

A 10 mm x 2 mm ID trap column was used for desalting and sample focussing prior to the MS analysis. The trap column was packed with 3 µm - 1000 Å C8 silica. Again, the large pore size of the stationary phase facilitates a narrow elution profile.



Figure 5: Base peak chromatogram of mAb and MS traces

No.	RT [min]	Mass	Protein	No.	RT [min]	Mass		Protein
1	22.4	10943	L37	24	54.5	23191	23647	L14 (native & with mod)
1	24.2	10943	L37	25	55.3	23345		L13a
2	26.5	6648	S30	26	55.3	16503		L27a
3	28.7	6276	L39	27	56.2	24015		L15
4	31.9	12321	L36a(L44)	28	56.7	9170		S21, N-acetylmethionine
5	36.7	17279	L26	29	56.7	14776		L23 N-acetylserine
6	36.7	9399	S27a	30	57.8	15954		S19
7	37.9	9270	S27a; cleaved C-term.	31	60	13284		S20
8	42	18449	L21 NG to KR	32	60.8	21527		L18
9	42.3	17623 17779	L24 & L24 (cterm)	33	61.4	29464		S4
10	42.3	18448	L21	34	63.2	29862		L7a
11	43.8	15667	L27	35	63.2	16314		S16
12	44.6	12122	L36	35	63.6	16314		S16
13	45.6	15644	L28	36	64.5	29466		S4
14	46.6	8087	L38	37	64.5	15379		S17
15	47.3	12465	L35a	38	65.8	17091		S13
16	47.3	23922	L13 (terminal KK)	39	66.8	17629		S18 N-acetylserine
16	47.7	23922	L13 (terminal KK)	40	66.8	22169		S7 N-acetylmeth
17	49.1	27908	L8	40	67.4	22169		S7 N-acetylmeth
17	49.5	27908	L8	41	69.5	22460		S9
18	49.5	14164	L31	42	71.6	14708		S15a
19	50	18343	S11	43	72.5	21893		L9
20	51	14421	L35	44	75.3	11772		P2
21	51.6	15465	S24	44	77	11852		P2
21	52.2	15466	S24	45	77	30355		L7
22	53.4	15727	L32	45	78.5	30356		L7
23	53.9	28680	S6	46	78.5	26585		S3 N-acetylalanine



# **Analysis of Membrane Proteins**

The experiment was performed in three steps: Separation of intact membrane protein sample derived from mouse liver

- using a 3 µm C8 column with 1000 Å pore size.
- 2. Fractions collected from 1. were digested with trypsin off-line.
- 300 Å pore size.

## **Chromatographic conditions for analytical separation**

System:	Agilent 1100 CapLC with Agilent MSD-
Column:	ProteCol™-C8 HQ1003 3 µm; 1000 Å ´
Sample:	3µl AOHUPO-MPI standard
Flow rate:	5.0 µl/min
Temperature:	80°C
Mobile Phase A:	0.1% formic acid in water
Mobile Phase B:	0.09% formic acid in acetonitrile

1-minute fractions of the intact protein separation were collected and digested with trypsin off-line. The dried fractions were reconstituted in 100 mM  $NH_4HCO_3$  and digested with trypsin overnight, acidified with 1% formic acid (FA), concentrated and re-diluted to 10  $\mu$ l with 1% FA.

Chromatograp	hic conditions
System:	TSP4000 pump, Surveyor autosampler,
Column:	ProteCol <sup>™</sup> -C18 HQ303 150 µm ID x 100
Sample volume:	10 µl
Mobile phase A:	0.1% FA in 5% acetonitrile
Mobile phase B:	0.1% FA in 90% acetonitrile
MS:	400-1500 mass range; top 6 ions fragm

#### **Data analysis**

Raw data files were converted to mzXML and searched against the Ensembl mouse database using Xtandem algorithm (GPM-XE software). The found proteins were statistically filtered either by number of identified peptides (>4) and/or by a cutoff score (10<sup>-10</sup>). As a result 542 proteins were identified with high confidence.

## Conclusions

By combining stationary phases with the right surface chemistry and the appropriate pore size for the analyte it is possible to separate even "difficult" samples such as membrane proteins. Our investigations into the effect of particle porosity led to the following recommendation for analyte size:

## Acknowledgements

The ribosomal protein results were obtained from Nicholas Williamson and Paul O'Donnell from Bio21 Mass Spectrometry and Peptide Institute, Bio21 Molecular Science and Biotechnology, Melbourne. Monoclonal antibody data was provided by Matthias Pelzing, Bruker Biosciences, Melbourne.

The membrane protein analysis was done in collaboration with Mibel Aguilar, Monash University and Paul Haynes, APAF, Australia.

#### R.P. Freeman, H.J. Wirth, A.A. Gooley SGE Incorporated, 2007 Kramer Lane, Austin Texas. 78758, USA

3. Digested fractions were analyzed by LC-MS/MS using 3 µm C18 column with

-iontrap MS Å 150 mm x 300 µm ID

Thermo linear ion trap ) mm

fragmented with 39% collision energy



Figure 6: Base peak chromatogram of membrane proteins on a capillary column







Figure 8: Relative pore and protein sizes for 1000 Å and 300 Å pores

