ACCESS DIAGNOSTIC TESTS UK

www.ukdrugtesting.co.uk

7/10/11/13/16/17 DRUG TEST CUP - INSTRUCTIONS

Open the foil-wrapped kit when you are ready to proceed with the test.


- Remove the cup device from the foil.
- Unscrew the cover and fill it with a fresh urine sample up to the level indicator.
- If the urine is at body temperature a green temp will be displayed on the temp strip.
- Peel the label & read the results at Five minutes (Do not read after more than 10 minutes.)

READ THE RESULTS

- Each test kit has a label for the Test line and a Control line. These are marked T, and C. before the test is performed no lines show in the result area. After the test is done the control line should show in all tests as a coloured line. The intensity of the test line may be less than the control line. Where the control line appears, any line in the test line area indicates a negative test. In positive tests the control line will show but no line will appear in the test zone. If no lines appear after 10 minutes the test is invalid. If only the test line appears with no control line the test is invalid (this is the in-built quality control)
- (if included) Compare the specific gravity and pH for the best colour match using the comparison chart enclosed with the test cup. If no match shows the sample may have been adulterated or diluted and the test should be repeated using a fresh sample

All Negative tests can be read at 1-2 minutes. Always allow a full five minutes development time for positive results. The line intensity may vary, but any line in the test zone should be read as a negative result, even if very feint. We recommend you confirm all positive tests. Prescribed medication may give positive results (e.g. opiate/codeine pain killers) Do not read beyond 10 minutes.

qqq

AMP = Amphetamine COC = Cocaine MET = Methamphetamine MOP = Heroin (& opiates) THC = Cannabis (marijuana) BZO = Benzodiazepines MTD = Methadone PPX = Propoxythene KET = Ketamine MDMA = Methylenedioxymethamphetamine (Ecstasy) PCP=PCP BAR = barbituates, TCA = Tri Cyclic Antidepressants BUP = Buprenorphine K2 = Synthetic Cannabis OXY= Oxycodone TML= Tramadol FYL= Fentanyl PGB= Pregabalin Adulteratio (Ox Oxidants SG specific gravity ph=pH NIT=Nitrite GLU= Glucose CRE= Creatinine)

FOR ADVICE ON RESULTS CALL 01263 731168 MON- FRI 8.30am – 5pm
Or email customerservices@accessdiagnostics.co.uk

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