


# 3Helix

Empowering collagen targeting for the diagnosis and treatment of human conditions



## What are CHPs?

Collagen is the most abundant protein in mammals. It is the major structural component of almost all organs and tissues. Excessive collagen remodeling is associated with a variety of chronic pathological conditions, such as cancer, arthritis, and fibrosis.


**3Helix's Collagen Hybridizing Peptide (CHP)** is a synthetic peptide that can specifically bind to denatured collagen by structural recognition of exposed alpha-strands. Our *in vivo* CHP is conjugated to a near-infrared dye to enable **accurate and reliable measurements** of collagen turnover in live animal models. CHPs can be used in sub-cutaneous, *in situ*, or intra-venous injections and a single animal can be dosed several times. Multiple images can be taken after a single injection and the fluorescent signal can be quantified using image analysis software, such as ImageJ.

## In vivo imaging of nAMD

in collaboration with Roche Ophthalmology


Neovascularized age-related macular degeneration (nAMD) is a condition that affects 17 million people worldwide and is accompanied by subretinal fibrosis. CHPs are able to directly detect fibrotic regions in a laser-induced model for nAMD in an *in vivo* setting and the signal can be easily quantified.

*In vivo* imaging and *ex vivo* fluorescent immunofluorescence of sulfo-cyanine 7.5 collagen hybridizing peptide (sCy7.5-CHP) binding to remodeling collagen in JR5558 mice. JR5558 mice received targeted or control sCy7.5-CHPs intravenously and were analyzed 5 days later using *in vivo* imaging followed by *ex vivo* IHC. (A) Representative images of JR5558 retinas showing infrared reflectance (IR), fluorescein angiography (FA), and sCy7.5-CHP *in vivo* binding. (B) MFI of control and targeted sCy7.5-CHPs in JR5558 retinas (control, n = 4; targeted, n = 7). (C) Representative images of retinal pigment epithelium (RPE)/choroid flat mounts from JR5558 mice (previously injected with sCy7.5-CHPs) co-stained for sCy3 collagen hybridizing peptide (R-CHP) binding, isolectin B4, and fibronectin.



## In vivo imaging of bone destruction in Multiple Myeloma

(A) Micro-CT scans of the spine from the multiple myeloma (MM) mouse showing abundant osteolytic lesions as "hollow spots", in comparison to a normal control. (B) Luciferase bioluminescence showing myeloma cell growth only in the MM mouse. (C) Near-IR CHP images of the MM mouse in (B) and a normal control mouse, demonstrating highly elevated signals in the MM mouse.




## Why Use CHPs

CHPs specifically visualize degraded collagen *in vivo*. Unlike competitor products, such as MMPsense, that detect enzyme activities as an indirect measure, CHP **binds directly to degraded collagen at the molecular level** allowing for accurate visualization of disease states. Our *in vivo* CHPs can be directly injected without pre-activation.

- Superb affinity and signal intensity that can be used for small animal imaging.
- CHP has a strong capability to hybridize with denatured collagen strands *in vivo* by reforming the triple helix.
- CHPs are able to detect actively remodeling collagen in fibrotic conditions over time.
- Our *in vivo* CHP maintains high affinity for denatured collagen without the need for any pre-workup.
- Applicable to all types of collagen in any species with essentially no nonspecific binding

### In vivo imaging of chronic kidney disease



**CHPs can be used to study a wide variety of diseases states**  
Subretinal fibrosis nAMD, Kidney fibrosis, Idiopathic pulmonary fibrosis (IPF), Pancreatic cancer, NASH

In Vivo CHPs	SKU	Size
Target-sCy7.5-CHP	INVIVOTGT7.5	3 Doses (8nmole)
Control-sCy7.5-scCHP	INVIVOCTL7.5	3 Doses (8nmole)
sCY7.5 In Vivo Kit	INVIVOKIT7.5	10 Doses Targeted (24nmole), 3 doses control (8nmole)

**3Helix.com**

orders@3helix.com (sales)

info@3helix.com (inquiries and technical support)

**3Helix**

675 S Arapeen Dr, STE 302  
Salt Lake City, UT 84108

(385) 722-4772