



### Safeguard Medical Technologies Ltd.

The Old Rectory Manufacturer Hope Under Dinmore Herefordshire, HR6 0PW, United Kingdom



Authorized in the EC

Willow Grove, Delgany,

Unique Device Identifier 5060483571515



REF Catalogue Number PDF166



Co Wicklow, A63 XY89, Ireland

Safeguard Technologies Ltd.

Swiss authorised

#### Procamed AG CH REP

Schulstrasse 5 CH-8355 Aadorf Switzerland

### Intended use

The SWAT-T <sup>™</sup> is a non-pneumatic tourniquet designed to stop severe traumatic bleeding before or during transport to a care facility.

#### User profile

The intended users are first responders, to include self-care. laypersons, emergency medical personnel and health care professionals. Special training is not required.

#### Patient profile

To be used on those patients with severe bleeding from a limb.

#### Clinical benefits

SWAT-T™ controls bleeding from an extremity by circumferentially compressing the injured limb, at least (~5cm above the bleed site) to an extent where arterial occlusive pressure is reached.

#### Performance characteristics

- Is available in two colours: black and orange
- · Can remain secure wet or dry
- Is highly portable

#### Product claims

Device has been designed to control bleeding from an extremity by circumferentially compressing the injured limb, above the bleed site, to an extent where arterial occlusive pressure is reached.

#### Contraindications

N/A

may cause necessary loss of extremity. Reassess for loosening. Dispose of as clinical waste. Do not re-use.

## RESIDUAL RISKS/WARNINGS: Safe for 2 hours. Longer appplication



















Batch code

# Instructions for use



1. Apply between wound and heart.



2. Wrap tightly while pulling tension.



3. To secure: Tuck end under itself.



4. Tourniquet if: Ovals/ rectangles turn to circles/square.

Reporting: Any serious incident that has occurred in relation to this device should be reported to manufacturer (vigilance@safequardmedical.com) and applicable national competent authority.