<p>Threshold™ PEP Positive Expiratory Pressure Device ENGLISH Instructions for Use Make sure there are no foreign objects in the mouthpiece and the device. For Single Patient Use Only. Ref. HS735</p>	<p>CONTENTS: Threshold™ PEP Positive Expiratory Pressure Device Instructions for Use Therapy Diary</p>	<p>1. Turn control knob (Fig.1,A) to align red edge of the pressure indicator (Fig. 1B) to the setting prescribed by your physician. A higher number equals greater effort. Attach mouthpiece firmly.</p>	<p>2. Put nose clip on your nose and breathe through your mouth.</p>	<p>3. Seal lips around mouthpiece, take a full breath in and breathe out 2 or 3 times longer than you breathe in. Continue this pattern for 10-20 breaths (Fig.3).</p>	<p>4. Remove mouthpiece from mouth and perform 2-3 "huff" coughs.</p>	<p>5. Repeat steps 3 and 4 approximately 4-6 times or until your therapy has lasted for about 10-20 minutes, or until secretions have been cleared. Therapy should be performed 2-4 times per day as needed or as instructed by your physician.</p>	<p>Important: Training should be consistent and on-going. Record readings in the Therapy Diary.</p>
--	---	---	--	--	---	---	--

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