REF	G-18232	G-32936	G-32997	G-36212	G-36248	G-36278	G-36306	G-36334	G-36366	G-36522	G-36546	G-36576
	G-18235	G-32937	G-32998	G-36214	G-36252	G-36282	G-36312	G-36336	G-36368	G-36526	G-36552	G-36578
	G-18913	G-32938	G-33438	G-36216	G-36256	G-36284	G-36314	G-36342	G-36502	G-36528	G-36554	G-36582
	G-18918	G-32942	G-33831	G-36218	G-36258	G-36286	G-36316	G-36344	G-36504	G-36532	G-36556	G-36586
	G-18919	G-32944	G-33960	G-36222	G-36262	G-36288	G-36318	G-36348	G-36506	G-36534	G-36558	G-36588
	G-18928	G-32947	G-36202	G-36228	G-36264	G-36292	G-36322	G-36352	G-36512	G-36536	G-36562	G-36592
	G-18944	G-32948	G-36204	G-36232	G-36266	G-36294	G-36324	G-36356	G-36514	G-36538	G-36564	G-38606
	G-31306	G-32949	G-36206	G-36236	G-36272	G-36296	G-36326	G-36362	G-36516	G-36540	G-36566	G-38634
	G-32935	G-32996	G-36208	G-36238	G-36274	G-36302	G-36328	G-36364	G-36518	G-36541	G-36568	

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TUBE GUIDED INSTRUMENTS

PRODUCT RELATED PROCESSING INSTRUCTIONS

These product related processing instructions for critical* medical devices (MD) are exclusively valid for the above listed items by GEUDER AG. All medical devices that are provided in a non-sterile way must be reprocessed before each application.


Fig. 1 example illustrations, not to scale

ACCESSORIES

Product list together with accessories and tools							Accessories	
							Cleaning adapter	
							 G-36608 / G-36609	
							 G-31013	
Forceps								
G-18913	G-32997	G-31306	G-32935	G-33438	G-38606			
G-32944	G-32947						G-36608	
G-18918	G-32937	G-36212 *)	G-36252	G-36274	G-36302	G-36326	G-36362	
G-18919	G-32938	G-36216	G-36256	G-36282	G-36306	G-36336	G-36366	
G-18928	G-33960	G-36222	G-36262	G-36286	G-36312	G-36342		
G-18944	G-36202	G-36232	G-36266	G-36292	G-36316	G-36352	G-36608	
G-32936	G-36206	G-36236	G-36272	G-36296	G-36322	G-36356		
G-33831							G-38608	
G-38634							G-36608	
G-36204	G-36218	G-36248	G-36278	G-36294	G-36324	G-36344	G-36368	
G-36208	G-36228	G-36258	G-36284	G-36314	G-36328	G-36348		
G-36214 *)	G-36238	G-36264	G-36288	G-36318	G-36334	G-36364	G-36609	
Scissors								
G-32942							G-36608	
G-32948							G-31013	
G-36502	G-36516	G-36532	G-36546	G-36562	G-36582			
G-36506	G-36522	G-36536	G-36552	G-36566	G-36586	G-36608		
G-36512	G-36526	G-36540	G-36556	G-36576	G-36592			
G-32996	G-36514	G-36534	G-36554	G-36568				
G-32998	G-36518	G-36538	G-36558	G-36578				
G-36504	G-36528	G-36541	G-36564	G-36588	G-36609			
Needle holder								
G-32949							G-31013	
Glaucoma punches								
G-18232	G-18235						G-36608	

Tab. 1: Accessories

OVERVIEW

Cycle limit	Assembly/Disassembly	Accessories
YES - *) G-36212, G-36214 25 X NO - or all the others listed above	YES	Cleaning adapter see list of accessories (Tab. 1)
Pre-cleaning	Ultrasound	Automated processing
Yes	YES	YES
Manual processing	Compressed air	Steam sterilization
NO	YEA	YES

Tab. 2: Overview

AUXILIARY EQUIPMENT

First treatment	Soft tissue (e.g. compress), single-use syringe (5 times with 10 ml each)
Pre-cleaning	Soft plastic brush or steamer (steam jet, max. 2 bar) single use syringe (5 times with 10 ml each) or cleaning gun (water pressure gun, max. 2 bar) optionally: ultrasound bath 35-40 kHz
Drying	Medical compressed air (oil-free compressed air, max. 2 bar)
Checking pH value	pH test strips
Maintenance	Instrument oil apt for sterilization, lint-free wipes

Tab. 3: Auxiliary equipment

Take care that the medical device is adequately and securely stored during processing; the instrument shaft and the functional parts of the tip could get damaged due to their delicate design.

ASSEMBLY/DISASSEMBLY

- Be careful not to touch the instrument shaft or the functional parts of the tip when assembling the cleaning adapter, as these can be damaged due to their delicate design.

FIRST TREATMENT AT PLACE OF USE

- Remove macroscopically visible impurities (e.g. pharmaceuticals or corrosive solutions, etc.) during/immediately after application with a compress soaked in sterile water.
- Place dry tube guided instrument into a suitable storage or transport system and close it.

PREPARATION BEFORE CLEANING

Pre-cleaning

- Rinse the tube guided instrument under running cold (< 40 °C) tap water.
- Remove visible, gross contamination using a soft plastic brush or a steamer.
- Attach the respective cleaning adapters (see table 1).
- Rinse the tube guided instrument with a cleaning gun with a suitable flushing adapter or a disposable syringe.

* Classification recommended by GEUDER AG, the responsibility lies with the operator.

Assembling cleaning adapter G-36608 und G-36609:

- Carefully slide the cleaning adapter onto the tube-guided instrument as shown in fig. 2a and 2b.


Fig. 2a

Fig. 2b

Assembling cleaning adapter G-31013:

- Squeezing the handle (instrument tip is closed), carefully insert the tip of the instrument in an upright vertical position into the cleaning adapter (fig. 4a).


Fig. 3a

- Lift the instrument into a horizontal position and continue introducing the instrument until the handle touches the cleaning adapter (fig. 4b).

Fig. 3b

- Please note that the tip must not touch the insides of the cleaning adapter (fig. 3b).

OR

Ultrasound cleaning as supporting pre-cleaning (optional)

Basically, we recommend pre-cleaning with ultrasound all complex or very dirty medical devices (e.g. instruments with a lumen, slide switch, joint, etc.). Observe the instructions or exposure time and concentration of the producer of the chemical as well as the instructions of the manufacturer of the ultrasound bath.

- Introduce the tube guided instrument and make sure that it is completely moistened:
 - Fill tube guided instruments with cleaning agent.
- Ultrasound treatment at 35-40 kHz, wait until the exposure time is over.
- Rinse the tube guided instrument afterwards with demineralized water for neutralization:
 - Move tube guided instrument when doing so.
- Check tube guided instrument for unobstructed flow:
 - Rinse the tube guided instrument with a cleaning gun, a steamer or a disposable syringe.

AUTOMATIC CLEANING AND DISINFECTION

As additional load only ophthalmic products are permitted.

Thermic disinfection A_0 value > 3000, 93 °C, 5 min; VARIO TD program.

The selection of suitable racks, charging trolleys for the optimal rinsing and flushing of the medical devices should be done after consulting with the cleaning device manufacturer.

Process

- Connect the Luer-Lock connector of the cleaning adapter to the mobile unit.
- Place the tube guided instrument in the rack (e.g. sieve baskets, trays):
 - When placing it in the mobile unit, make sure that the instrument tips do not touch the cleaning adapters.
 - Avoid rinsing shadows.
 - After cleaning and disinfection remove the cleaning adapter in reverse order.

DRYING

Dry the tube guided instrument with medical compressed air, if necessary.

pH-Wert Kontrolle

For ophthalmic instruments a batch-related routine control for alkaline residues is recommended.

- Apply remaining humidity from the tube onto a pH test strip.
- Compare the pH test strip with the color scale:
 - pH must be neutral.

CONTROL AND MAINTENANCE

Visual control

- Check if the tube guided instrument is clean:
 - In case of residues, the cleaning and disinfection process must be repeated.
- Check if the tube guided instrument is damaged:
 - Remove corroded or tube guided instruments.
- Check if the tube guided instrument is complete.
- Check if the tube guided instrument is dry.

Maintenance

- Care for tube guided instrument with a suitable instrument oil:
 - let the tube guided instrument cool down.
 - Apply the instrument oil selectively on the area that needs it (fig. 5a and fig. 5b).
 - Remove excess instrument oil, if any, with a lint-free wipe.
- Carry out a functional control:
 - Check if the tube guided instrument is complete.
 - Check if the moveable parts can be moved easily.


Fig. 4a


Fig. 4b

PACKAGING

Use packaging materials and systems complying with the following standards: DIN EN (ANSI AAMI) ISO 11607 and DIN EN 868-8.

STERILIZATION

Steam sterilization in fractionated prevacuum method:

SAL 10⁶, at least 3 min at 134 °C, 3 bar.

STORAGE

Storage is at room temperature (max. 25 °C) in closed storage systems (dust- and light-protected also dry), protected from damage or mechanical influences, in a controlled environment.

- Damaged or open sterile barrier systems (SBS) are deemed non-sterile.
- The storage time depends on the type of packaging and must be validated accordingly.

TRANSPORT

Transport must not have a negative effect on the characteristics of the sterile goods.

Returns to GEUDER AG: Please note that we only accept processed (cleaned, disinfected and sterilized) medical devices.

DISPOSAL

The laws and provisions of each country must be complied with. Medical devices suspected to contain Creutzfeldt-Jakob Disease (CJD, vCJD): observe the recommendations of the Robert-Koch-Institute (www.rki.de).

National standards, regulations and/or restrictions must also be considered!

GENERAL INFORMATION FOR PROCESSING

Detailed information for processing is provided in the manufacturer information for the processing of medical devices from GEUDER AG.

The life cycle of reusable medical devices ends, either due to material changes, material wear, damages, as well as functional changes or loss of function of the medical product, or after reaching the maximum number of reprocessing cycles in case of a cycle limitation as stated in the product related processing instructions for the product (see table overview).

DETAILS OF VALIDATION

Summary of processes validated by GEUDER AG:

The selection of the parameters refers to worst case conditions.

Manual pre-cleaning

Pre-cleaning

Using tap water (according TrinkwV [German Drinking Water Ordinance]), cleaning gun or syringe, instrument cleaning brush with soft nylon bristles (15 mm).

Cleaning

- In the ultrasound bath, 0.5 %, 10 min/30-40 °C, 35 kHz, Bandelin RK 510 H, after ultrasound treatment the instruments are placed for 2 min in demineralized water and rinsed and/or operated.

Process chemicals (manual)

- Cleaning: Dr. Weigert: neodisher MediClean forte, demineralized water.

Automatic cleaning and disinfection

Cleaning

- Program: shortened VARIO TD program, 0.2 %, 55 °C/10 min, final rinsing 3 min/demineralized water, Miele & Cie. KG: Washer-Disinfector type G7836, charging trolley: E450/1 accessories for MIC instruments, sieve basket with lid.

Disinfection

- Thermal high level disinfection A_0 -value > 3000, 93 °C/5 min, Miele & Cie. KG: Washer-Disinfector type G7836, charging trolley: E450/1 accessories for MIC instruments, sieve basket with lid.
- Thermal disinfection A_0 -value = 600, 90 °C/1 min, Miele & Cie. KG: Washer-Disinfector type G7836, charging trolley: E450/1 accessories for MIC instruments, sieve basket with lid.

Process chemicals (automatic)

- Cleaning: Dr. Weigert: neodisher MediClean forte.
- Final rinsing: demineralized water.
- Disinfection: demineralized water.

Steam sterilization

- Program: half-cycle SAL 10⁶, fractionated pre-vacuum procedure 3 times, 132 °C/1.5 min, MMM GmbH: type Selectomat PL/669-2 CL, sterilization pouches, drying time 20 min, cooling time: outside of the autoclave 30 min at room temperature (max. 25 °C).
- The maximum sterilization temperature is 138 °C.

An independent, accredited test laboratory proved the basic suitability of the medical devices for an effective processing.

The above mentioned instructions were validated by the manufacturer of medical devices to be suitable for the processing/reprocessing of a medical product to be reused. The processing/reprocessing person(s) have the responsibility, that the actually applied procedure with the utilized equipment, materials and staff in the processing facility achieves the required result. Therefore the verification and/or validation and a routine monitoring of the procedure are indispensable.

Actual version available at GEUDER AG.

Subject to changes.