REF G-24205, 28133, G-28137, G-28140, G-28142, G-30139, G 31070, G-32697, G-33475, G-33919

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SILICONE TUBING SYSTEMS
PRODUCT RELATED PROCESSING INSTRUCTIONS

These product related processing instructions for semi-critical* medical devices (MD) are exclusively valid for the above listed items by GEUDER AG. All medical devices that are provided in a non-sterile way must be reprocessed before each application.


Fig. 1: example illustrations, not to scale

OVERVIEW

Cycle limit	Assembly/Disassembly	Accessories
YES 25 X	NO	NO
Pre-cleaning	Ultrasound	Automated processing
YES	YES	YES
Manual processing	Compressed air	Steam sterilization
YES	YES	YES

Tab. 1: Overview

AUXILIARY EQUIPMENT

First treatment	Soft tissue (e.g. compress), single-use syringe
Pre-cleaning	Soft plastic brush or steamer (steam jet, max. 2 bar), single use syringe or cleaning gun (water pressure gun, max. 2 bar) optionally: ultrasound bath 35-40 kHz
Manual cleaning and disinfection	Medical compressed air (oil-free compressed air, max. 2 bar)
Drying	pH test strips

Tab. 2: Accessories

FIRST TREATMENT AT PLACE OF USE

- Remove macroscopically visible impurities (e.g. pharmaceuticals or corrosive solutions, etc.) during immediately after application with a compress soaked in sterile water.
- Place dry silicone tubing system into a suitable storage or transport system and close it.

PREPARATION BEFORE CLEANING

Pre-cleaning

- Rinse the silicone tubing system under running cold (< 40 °C) tap water.
- Remove visible, gross contamination using a soft plastic brush or a steamer.
- Rinse the silicone tubing system with a cleaning gun with a suitable flushing adapter or a single-use syringe.
- Rinse the silicone tubing system against the flow direction.

OR

Ultrasound cleaning as supporting pre-cleaning (optional)

Basically, we recommend pre-cleaning with ultrasound for all complex or very dirty medical devices (e.g. instruments with a lumen, slide switch, joint, etc.). Observe the instructions of exposure time and concentration of the producer of the chemical as well as the instructions of the manufacturer of the ultrasound bath.

- Introduce the silicone tubing system and make sure that it is completely moistened:
 - Fill the silicone tubing system with cleaning agent.
- Ultrasound treatment at 35-40 kHz, wait until the exposure time is over.
- Rinse the medical devices afterwards with demineralized water for neutralization.
- Check the silicone tubing system for unobstructed flow:
 - Rinse the silicone tubing system with a cleaning gun, a steamer or a single-use syringe.
 - Rinse the silicone tubing system against the flow direction.

MANUAL CLEANING AND DISINFECTION IN ULTRASOUND BATH

Manual cleaning in ultrasound bath

- Introduce the silicone tubing system and make sure that it is completely moistened:
 - Fill the silicone tubing system with cleaning agent.
- Ultrasound treatment at 35-40 kHz, wait until exposure time is over.

- Afterwards, place the silicone tubing system into demineralized water for at least 2 minutes for neutralization:
 - Remove residues of debris with a soft plastic brush (without pressure).
 - Rinse the complete tubing system at least five times with demineralized water.
- Rinse the silicone tubing system under running demineralized water:
 - Check the silicone tubing system for unobstructed flow and rinse it with a cleaning gun or a single-use syringe (rinse the silicone tubing system against flow direction).

Manual disinfection in ultrasound bath

- Introduce the silicone tubing system and make sure that it is completely moistened:
 - Fill lumen instruments with disinfection agent.
- Ultrasound treatment at 35-40 kHz, wait until exposure time is over.
- Afterwards, place the silicone tubing system into demineralized water for at least 2 minutes for neutralization:
 - Rinse the complete tubing system at least five times with demineralized water.
- Rinse the silicone tubing system under running demineralized water.
- Check the silicone tubing system for unobstructed flow and rinse it with a cleaning gun or a single-use syringe (rinse the silicone tubing system against flow direction).

AUTOMATIC CLEANING AND DISINFECTION

As additional load only ophthalmic products are permitted. Thermic disinfection A₀ value > 3000, 93 °C, 5 min; VARIO TD program. The selection of suitable racks, charging trolleys for optimal rinsing and flushing of the medical devices should be done after consulting with the cleaning device manufacturer.

Process

- Place the silicone tubing system in the rack (e.g. sieve baskets, trays):
 - Avoid rinsing shadows.
 - Avoid kinks.
- Connect the silicone tubing system to the charging trolley.
- Tubing systems with joint:
 - For cleaning/disinfection the joints stay closed (fig. 2).


Fig. 2

* Classification recommended by GEUDER AG, the responsibility lies with the operator.

DRYING

Dry the medical devices with medical compressed air, if necessary.

Control of pH value

For ophthalmic instruments a batch-related routine control for alkaline residues is recommended.

- Apply the remaining humidity of the silicone tubing system onto a pH test strip.
- Compare the pH test strip with the color scale:
 - pH must be neutral.

CONTROL AND MAINTENANCE

Visual control

- Check if the silicone tubing system is clean:
 - In case of residues, the cleaning and disinfection process must be repeated.
- Check if the silicone tubing system is damaged:
 - Remove corroded or damaged medical devices.
- Check if the silicone tubing system is complete.
- Check if the silicone tubing system is dry.

Maintenance

- Carry out a functional control:
 - Check if the silicone tubing system is complete.
 - Check the connectors.
- Tubing systems with joint:
 - Check the joints.

PACKAGING

Use packaging materials and systems complying with the following standards: DIN EN (ANSI AAMI) ISO 11607 and DIN EN 868-8.

STERILIZATION

Steam sterilization in fractionated prevacuum method:

- SAL 10⁻⁶, at least 3 min at 134 °C, 3 bar.
- Tubing systems with joint:
 - For sterilization the joints must be opened (fig. 3).


Fig. 3

STORAGE

Storage is at room temperature (max. 25 °C) in closed storage systems (dust- and light-protected also dry), protected from damage or mechanical influences, in a controlled environment.

- Damaged or open sterile barrier systems (SBS) are deemed non-sterile.
- The storage time depends on the type of packaging and must be validated accordingly.

TRANSPORT

Transport must not have a negative effect on the characteristics of the sterile goods.

Returns to GEUDER AG: Please note that we only accept processed (cleaned, disinfected and sterilized) medical devices.

DISPOSAL

The laws and provisions of each country must be complied with. Medical devices suspected to contain Creutzfeldt-Jakob Disease (CJD, vCJD): observe the recommendations of the Robert-Koch-Institute (www.rki.de).

National standards, regulations and/or restrictions must also be considered!

GENERAL INFORMATION FOR PROCESSING

Detailed information for processing is provided in the manufacturer information for the processing of medical devices from GEUDER AG.

The life cycle of reusable medical devices ends, either due to material changes, material wear, damages, as well as functional changes or loss of function of the medical product, or after reaching the maximum number of reprocessing cycles in case of a cycle limitation as stated in the product related processing instructions for the product (see table overview).

DETAILS OF VALIDATION

Summary of processes validated by GEUDER AG:

The selection of the parameters refers to worst case conditions.

Manual pre-cleaning, cleaning and disinfection

Pre-cleaning

Using tap water (according TrinkwV [German Drinking Water Ordinance]), cleaning gun or syringe, instrument cleaning brush with soft nylon bristles (15 mm).

Cleaning

- In the ultrasound bath, 0.5 %, 10 min/30-40 °C, 35 kHz, Bandelin RK 510 H, after ultrasound treatment the instruments are placed for 2 min in demineralized water and rinsed and/or operated.

Disinfection

- 3.0 %, 15 min/20 °C, immersion bath, lumen were filled with disinfectant agent. After disinfection time rinse for 1 min with demineralized water.

Process chemicals (manual)

- Cleaning: Dr. Weigert: neodisher MediClean forte, demineralized water.
- Disinfection: Bode Chemie: Korsorex plus, demineralized water.

Automatic cleaning and disinfection

Cleaning

- Program: shortened VARIO TD program, 0.2 %, 55 °C/10 min, final rinsing 3 min/demineralized water, Miele & Cie. KG: Washer-Disinfector type G7836, charging trolley: E450/1 accessories for MIC instruments, sieve basket with lid.

Disinfection

- Thermal high level disinfection A₀-value > 3000, 93 °C/5 min, Miele & Cie. KG: Washer-Disinfector type G7836, charging trolley: E450/1 accessories for MIC instruments, sieve basket with lid.
- Thermal disinfection A₀-value = 600, 90 °C/1 min, Miele & Cie. KG: Washer-Disinfector type G7836, charging trolley: E450/1 accessories for MIC instruments, sieve basket with lid.

Process chemicals (automatic)

- Cleaning: Dr. Weigert: neodisher MediClean forte.
- Final rinsing: demineralized water.
- Disinfection: demineralized water.

Steam sterilization

- Program: half-cycle SAL 10⁻⁶, fractionated pre-vacuum procedure 3 times, 132 °C/1.5 min, MMM GmbH: type Selectomat PL/669-2 CL, sterilization pouches, drying time 20 min, cooling time: outside of the autoclave 30 min at room temperature (max. 25 °C).
- The maximum sterilization temperature is 138 °C.

An independent, accredited test laboratory proved the basic suitability of the medical devices for an effective reprocessing.

The above mentioned instructions were validated by the manufacturer of medical devices to be suitable for the processing/reprocessing of a medical product to be reused. The processing/reprocessing person(s) have the responsibility, that the actually applied procedure with the utilized equipment, materials and staff in the processing facility achieves the required result. Therefore the verification and/or validation and a routine monitoring of the procedure are indispensable.

Serious incidents related to the medical device must be reported to the manufacturer as well as to the competent authority.

Actual version available at GEUDER AG.

Subject to changes.