

# Medical Carriers User Guide



By Appointment to  
Her Majesty the Queen  
Supplier of Tamper Evident  
Bags and Security Seals  
Versapak International  
Limited, Erith, Kent

Tested to maintain temperature of 2-8 degrees

Product Name/Code	Carrier	Thermal Barriers	Freeze Boards	Ambient Temperature	*Temperature maintained for
Small Pathology Sample Carrier (PYTB1)		Thermal Barrier can be pre-conditioned to +4 degrees 	Freeze Boards should be pre-conditioned to -18 degrees 	 21°C 30°C	 6hrs 53min 4hrs 07min
Small Vaccine Carrier (PYTBVAC1)		Small (THBARRIER_SML)	(FREEZEB400)		
Medium Pathology Sample Carrier (PYTB2)				 21°C 30°C	 22hrs 11min 9hrs 41min
Medium Vaccine Carrier (PYTBVAC2)		Large (THBARRIER_LRG)	(FREEZEB750)		
Large Pathology Sample Carrier (PYTB3)				 21°C 30°C	 >24 hours 12hrs 39min
Large Vaccine Carrier (PYTBVAC3)		Small (THBARRIER_SML)	(FREEZEB750)		


## PLACEMENT



Place contents inside



Place Thermal Barriers



Place Freeze Boards



Close and Seal



\*Study by TuV GmbH. Freeze Boards preconditioned to -18°C for 24 hours. Thermal Barrier preconditioned to 4 degrees. For more information please contact us.