Indicadores Multi-variable

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ES

Multi-variable Indicators

For Steam sterilization processes

Type 4 Indicators according to ISO 11140-1:2014

Product description

Chemdye® Multi-variable Indicators have been designed to react to Steam sterilization processes within loads, ensuring an adequate control of the effectiveness of sterilization processes (temperature, time and steam quality). Yellow ink has been developed to turn to black when the process reaches the stated values of the critical process variables. Multi-variable indicators allow to rapidly detect whether the critical parameters of the sterilization process for which they have been designed to respond have been reached.

Characteristics

CD29: One-side laminated punched strip of 200 x 15 mm, this allows to use the entire strip or half strip depending on the size of the package to be sterilized.

CD29-AD: Self-adhesive one-side laminated punched double card strips 200 x 15 mm. This indicator version is adhesive, for an easy record results.

Indicators 100% free of toxic heavy metals.

Indications for use

Use Chemdye® Multi-variable Indicators for monitoring Steam sterilization processes. Conditions: 121 °C, 15 minutes // 134 °C, 3.5 minutes.

Precautions

Do not expose Chemdye® Multi-variable Indicators to Ethylene Oxide, Dry Heat or any sterilization processes other than Steam.

Instructions for use

• Place a Chemdye® Multi-variable Indicator within each package, bag or tray that will be sterilized by Steam.

• In rigid containers place the indicator in every corner, or at least two diagonally opposite corners, to ensure that the sterilization agent has reached the most inaccessible areas.

• After completing the process analyze the results. The indicator must change to the reference color for considering that the established conditions were reached.

Storage and shelf life

• Indicators should be stored in a dry place, protected from light at temperatures between 10-30 °C, 30-80 % Relative Humidity. Do not wet. Do not store close to sterilizing agents.

• Indicators have an expiration date of 5 years from the date of manufacture when stored in recommended conditions. Both the expiration date and batch number of the product are specified on the package label and in each unit of the product. Do not use indicators after their expiration date.

• After use, the indicator will not change within 1 year when stored at previously indicated conditions.

Disposal

Discard chemical indicators after use along with paper waste, according to your country's healthcare and safety regulations.

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Unprocessed

Failures

One or more sterilization parameters were not achieved. The materials cannot be considered sterile and must be processed again.

 The main failure causes are:

- Inadequate cycle duration
- Improper cycle selection
- Overloading / incorrect loading of the sterilizer

- Poor steam quality
- Sterilizer malfunction
- Blocking indicator with large metal mass

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NOTE: Reference colors exhibited in printed prospects as well as those displayed in our website and soft copies of other documents, do not necessarily represent the real color as shown in the actual indicators.