

Medical Carriers User Guide


Tested to maintain temperature of 2-8 degrees

Product Name/Code	Carrier	Thermal Barriers	Freeze Boards	Ambient Temperature	*Temperature maintained for
Small Pathology Sample Carrier (PYTB1)		 Small (THBARRIER_SML)	 (FREEZEB400)	 21°C	 6hrs 53min
Small Vaccine Carrier (PYTBVAC1)				 30°C	 4hrs 07min
Medium Pathology Sample Carrier (PYTB2)		 Large (THBARRIER_LRG)	 (FREEZEB750)	 21°C	 22hrs 11min
Medium Vaccine Carrier (PYTBVAC2)				 30°C	 9hrs 41min
Large Pathology Sample Carrier (PYTB3)		 Small (THBARRIER_SML)	 (FREEZEB750)	 21°C	 >24 hours
Large Vaccine Carrier (PYTBVAC3)				 30°C	 12hrs 39min


PLACEMENT



Place contents inside



Place Thermal Barriers



Place Freeze Boards



Close and Seal



*Study by TuV GmbH. Freeze Boards preconditioned to -18°C for 24 hours. Thermal Barrier preconditioned to 4 degrees. For more information please contact us.