





## Fitting Guide



## **LegStyle Fitting Guide – Detailed Instructions (FlexFit™ System)**

### **Step 1 – Clean**

Start with a clean and dry surface to ensure the adhesive tape sticks well.

- Use mild soap or detergent to clean your socket.
  - Rinse thoroughly with clean water.
  - Dry with a towel — avoid hairdryers or direct heat.
  - Make sure the surface is free from oils, lotions, or residue.
- 

### **Step 2 – Mark**

Position your LegStyle cover and mark the socket for tape placement.

- Place the cover around your socket and fasten the zipper fully.
  - Adjust the rotation as needed, and slide the cover upward until it sits flush with the socket.
  - Leave a gap of approximately 3 cm between the bottom of the cover and your foot shell.
  - Use a felt-tip pen or pencil to lightly mark the top edge of the cover on your socket.
- 

### **Step 3 – Adhere**

Apply the adhesive tape to secure the cover.

- Use the provided 5 cm wide tape.
  - Apply it about 1 cm above your marked line.
  - Stretch gently as you apply it to follow the socket's shape.
  - Overlapping is perfectly fine if you reach the end with extra tape.
  - If you ever need to reposition, each order includes a spare tape set.
  - You can layer new tape over old if the original shifts over time.
-

**Step 4 – Wear**

Attach and adjust the cover.

- Partially zip the cover (about half way)
- Align the top of the cover with the tape, adjusting for symmetry.
- Zip it closed like a well-fitted jacket.

Need help? Check out our [Support and FAQ](#) at for extra guidance