Descriptive Folder

Description

Exacta Number One is a single use infusion device, reliable, for infusion therapies requiring administration in continuos or with bolus. The combination of predefined flow and volume makes this infusion suitable for various applications. As a result of its elastomeric property, it does not require a power source such as a motor or batteries. EN1 permits continuous drug delivery. Exacta Number One Bolus is supplied with a control module (PCA) which additionally allows delivery of a predefined bolus of drug according to patient's needs. The module Multirate (Exacta Number One Multirate) enables the health professional to set and vary, tought the special key supplied, the flow rate also during the infusion. Exacta Number One Multirate with Bolus is supplied with both optional modules. All EN1 models have got an exclusive flow regulator that reduces the risk of any particlesand crystals occlusion thanks to its lumen which is larger than the glass capillary. All EN1 pumps are tested one by one.

Indications

EN1 is a single use medical device which can be used in complete safety for infusion therapies requiring continuous administration by intravenous, intra-arterial, peridural, subcutaneous routes and for infusion in soft tissue and body cavities.

EN1 is suitable for the continuos administration of cancer , anaesthetic , iron chelation (DFO) , antibiotic and analgesic therapies and for the treatment of acute and chronic pain, EN1 Bolus allows an additionally delivery (0,5 ml, 1 ml, 2 ml) in a predefined time (8, 15, 30, 60 min). EN1 Multirate enables to choose four different flows.

Precautions

- Carefully read this descriptive folder before using.
- Determine the integrity of the product packaging. If it has been opened or damaged, do not use.
- EN1 is for single use only. Do not re-sterilize or reuse.
- EN1 should be stored at a temperature less than 38° C and should be kept away from direct light.
- When preparing and administering the drug, follow the instructions supplied by the producer and the prescribing physician.
- When preparing and administering antiblastic chemotherapy drugs it is necessary to adopt the established prevention systems individual personal protection devices (DPI) as established by the norms in vigour.
- Keep out the special key of reach of patients (Multirate and Multirate with bolus models).

Collateral effects

The use of the device does not cause sensibilities or allergies.


Instruction for use

EN1 is a medical device that must be prepared using strict aseptic techniques.

- Before opening the package, check that the flow rate and the bolus shown on the package are the same as that required for use.
- Clamp the tube and remove the Luer lock on the filling port.
- Fill a sterile syringe with the drug.
- Connect the syringe to the fill port, and inject the medication into the membrane. This operation must be executed with high caution in order not to damage the valve.
- Do not fill beyond the allowed volume.
- After filling the reservoir, remove the syringe and close the fill port with the supplied plug..
- Remove the clamp and the plug on the flow regulator.
- Wait until the tube is filled and all the air has been expelled.
- When the first drops of the solution begin to flow freely, the air has been completely expelled and the system is ready to function.
- In case the liquid to be infused does not flow freely from the flow regulator, use tubing connected to a syringe with a Luer lock and aspirate.
- After starting the pump, connect the flow regulator on the device to the patient's access route.
- Remove the clamp.
- Remove the coloured slice of the bolus after exhausting the air (PCA Module).
- Use the special key to manage the flow-rate in the Multirate module. Ensure that it remains in the exclusive availability of the health care professional.

Single use

As indicated by the signs on the etiquette and in the recommendations contained in the present instructions for use, do not reuse the device because it is no more sterile and the elastomer performance after one usage is no longer guaranteed.

Do not reuse: biohazard

Technical information

- EN1 is a single use medical device with a predefined fixed flow rate
- EN1 Bolus has a control module to deliver a predetermined additional dose of drug .
- EN1 is a medical device that allows the exclusive delivery of the flow rate shown on the package.
- EN1 Multirate gives the opportunity to set and vary the flow rates during the infusion.
- The nominal flow rate at the level of the delivery flow is based on a temperature of 32°C/89.6° F.
- Tests of the accuracy of the tube flow rate have been performed using a 5% glucose solution (DSW). The flow rate may increase or decrease (± 10%) with decreasing or increasing viscosity of the solution to be injected.
- The flow rate has been tested by placing the reservoir and the flow delivery device at the same level. If the reservoir is positioned above the delivery device, the flow rate will tend to increase.

Specifications

Flows	From 0.5ml/h to 500ml/h	
Volumes and System Accuracy	Volume 60 ml, 100 ml, 150 ml 200 ml, 275 ml, 350ml	Accuracy ± 10% nominal rate
Residual volume	0.45 ml for 60 ml volume 0.70 ml for 100 ml, 150ml volume 1.20 ml for 200, 275, 350 ml volume	
Filter	Particles 1,2µ (on demand 5,0µ) Air venting 0.02µ	
Bolus	Volume 0.5 ml, 1 ml, 2 ml	Min. 8, 15, 30, 60
Multirate	From 0.5 ml/h to 12ml/h (four flow rates in each model)	
Materials	DEPH Free, Latex Free	

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