



# HOW TO PROPERLY PUT ON YOUR OMNA TOURNIQUET LEASH / LEG ROPE

**Step 1:** Visually inspect that all components are in good working order, and that the ladder strap is at full length. *(Addendum explains in detail)*





# HOW TO PROPERLY PUT ON YOUR OMNA TOURNIQUET LEASH / LEG ROPE

**Step 2:** Place the cuff so that the cord is on the back of your ankle, and wrap the neoprene around your ankle to engage a small amount of Velcro (*Just enough to hold it in place*).





# HOW TO PROPERLY PUT ON YOUR OMNA TOURNIQUET LEASH / LEG ROPE

**Step 3:** Grasp the D-ring (*and leash cord if needed*) and pull tension to remove any slack while attaching remainder of Velcro. (*This prevents bunching, which may create a poor Velcro connection while surfing*).





# HOW TO PROPERLY PUT ON YOUR OMNA TOURNIQUET LEASH / LEG ROPE

**Properly Applied – *No Bunching***



**Improperly Applied - *Bunching***





# HOW TO PROPERLY PUT ON YOUR OMNA TOURNIQUET LEASH / LEG ROPE

**Step 4:** Ensure the area between the D-ring and Ratchet is attached well, confirm by pulling on bite-strap, repeat steps 2-4 as needed.







# HOW TO PROPERLY PUT ON YOUR OMNA TOURNIQUET LEASH / LEG ROPE

**Step 5: Go Surf!**





# HOW TO PROPERLY PUT ON YOUR OMNA TOURNIQUET LEASH / LEG ROPE

## Addendum:

1. To verify the ladder strap is at full length press the Red release button and pull the ladder strap out of the ratchet, then reinsert the ladder strap back into the ratchet until you hear **1** audible click, pull to confirm that it is now locked in place.







# HOW TO PROPERLY PUT ON YOUR OMNA TOURNIQUET LEASH / LEG ROPE

## Addendum:

1. Remember the ratchet is only used when applying the device in its emergency tourniquet configuration.

