



Simplifying Protein Expression

Product Manual

Product name: GlycFree (λ DE3) Chemically Competent *E. coli* Cells.

Catalog# GF-201, **Size:** 10 X 0.1 ml/tube

Product details:

The protein expression strain GlycFree (λ DE3) is the λ DE3 lysogen of *E. coli* M226 strain (M226(λ DE3)). It has been engineered for expression of sugar-free recombinant proteins in *E. coli*. Some glycosyltransferases have ability to catalyze a self-glycosylation reaction using UDP-glucose (UDPG) as glucose donor. As such, these proteins can be aberrantly glycosylated during expression in common *E. coli* strains like BL21 and its derivatives. GlycFree (λ DE3) is an *E. coli* strain defective in the synthesis of UDP-glucose. GlycFree(λ DE3) allows to produce non-glycosylated forms of target proteins in native conformation and preserve their biological activity.

GlycFree (λ DE3) strain carries a copy of the T7 RNA polymerase under control of the IPTG inducible lacUV5 promoter, and as a result, is suitable for controlled expression of both *E. coli* and T7 promoter driven constructs.

Δ UDPG strains are not susceptible to phage P1 infection.

Transformation efficiency: $>1 \times 10^5$ cfu/ μ g pUC19

E. coli genotype: *e14- galU106 trp-50 relA1 spoT1 rpsL150(str^R) (λ DE3)*

[λ DE3 = λ sBamHIo Δ EcoRI-B int::(*lacI::PlacUV5::T7 gene1*) *i21* Δ nin5]

Suggested Transformation Procedure for Optimal Results:

1. Remove cells from -80°C and let thaw on ice.
2. Gently mix cells by lightly flicking tube. Aliquot ~ 50 - $100 \mu\text{l}$ of cells into chilled, 17×100 mm polypropylene tube(s). Unused cells may be refrozen, but a drop in efficiency may result. For optimal recovery, refreeze cells in a dry ice/ ethanol bath prior to -80°C storage.
3. Add DNA solution ($\leq 5 \mu\text{l}$ per $50 \mu\text{l}$ cells) to cell suspension and gently swirl tube(s) for a few seconds to mix.
4. Incubate on ice for 30 minutes.
5. Place tube(s) in 42°C water bath for ~ 30 to 45 seconds without shaking. For $50 \mu\text{l}$ aliquots, 30 seconds is recommended for maximum efficiency.
6. Place tube(s) on ice for ~ 2 minutes.
7. Add 200 - $250 \mu\text{l}$ SOC to transformation reaction(s).
SOC medium: 2% Tryptone, 0.5% Yeast Extract, 0.4% glucose, 10 mM NaCl, 2.5 mM KCl, 10 mM MgCl_2 & 10 mM MgSO_4 . Other media can be used to grow transformed cells, including standard LB or TB broths. However, SOC is the optimal choice for recovery of the cells and for obtaining maximum transformation efficiencies.
8. Shake tube(s) ~ 200 rpm for 1 hour at 37°C .
9. Plate by spreading 100 - $200 \mu\text{l}$ of cell transformation mixture on LB agar plates containing appropriate antibiotic and incubate overnight at 37°C .

Notes: If necessary, the entire amount of transformation mixture can be spread on the plates. For that, spin 1 ml of the mixture in microcentrifuge at $10,000$ rpm for 30 seconds to pellet cells. Remove 800 - $900 \mu\text{l}$ of supernatant. Resuspend cells in the remaining volume of media (100 - $200 \mu\text{l}$) and plate cells (using all volume) onto a pre-warmed selective plate.

Storage is recommended at -80°C .

Usage: This product is intended for LABORATORY RESEARCH USE ONLY. Not for diagnostic or therapeutic use.