



## Pinnacle BioLabs COVID-19 Novel Coronavirus Dual IgG/IgM Rapid Test Package Insert

### Intended Use

The Pinnacle BioLabs COVID-19 Novel Coronavirus (COVID-19)IgG/IgM Rapid Test (Whole Blood/Serum/Plasma) is a lateral flow immunoassay for the qualitative detection and differentiation of IgG and IgM of Novel Coronavirus(COVID-19) in human whole blood, serum or plasma. This test is intended to be used as an aid in the diagnosis of infection with Novel Coronavirus. Any reactive specimen with the Novel Coronavirus(COVID-19)IgG/IgM Rapid Test must be confirmed with alternative testing method(s) as this test is currently under an Emergency Use Authorization for professional use. **For professional use only. This test has not been reviewed by the FDA. Negative results do not rule out SARS-CoV-2 infection, particularly in those who have been in contact with the virus. Follow-up testing with a molecular diagnostic should be considered to rule out infection in these individuals. Results from antibody testing should not be used as the sole basis to diagnose or exclude SARS-CoV-2 infection or to inform infection status. Positive results may be due to past or present infection with non-SARS-CoV-2 coronavirus strains, such as coronavirus HKU1, NL63, OC43, or 229E. Not for the screening of donated blood.**

### Summary and Explanation of the Test

Coronavirus (CoV) belongs to the genus Nestovirus, Coronaviridae, and is divided into three genera:  $\alpha$ ,  $\beta$ , and  $\gamma$ . The  $\alpha$  and  $\beta$  gene are only pathogenic to mammals. The  $\gamma$  gene mainly causes bird infections. CoV is mainly transmitted through direct contact with secretions or through aerosols and droplets. There is also evidence that it can be transmitted through the fecal-oral route.

So far, there are 7 types of human coronaviruses (HCoV) that cause human respiratory diseases: HCoV-229E, HCoV-OC43, SARS-CoV, HCoV-NL63, HCoV-HKU1, MERS-CoV and Novel Coronavirus(COVID-19) (2019), it's an important pathogen of human respiratory infections. Among them, a COVID-19 was discovered in 2019 due to Wuhan virus pneumonia cases. The clinical manifestations are systemic symptoms such as fever and fatigue, accompanied by dry cough, dyspnea and so on. It can quickly develop into severe pneumonia, respiratory failure, acute respiratory distress syndrome, septic shock, multiple organ failure, severe acid-base metabolism disorders, etc., and even life-threatening.

### Test Principle

The Pinnacle BioLabs COVID-19 Novel Coronavirus (COVID-19)IgG/IgM Rapid Test is a lateral flow chromatographic immunoassay. The test cassette consists of: 1) a burgundy colored conjugate pad containing recombinant COVID-19 antigen conjugated with colloidal gold (COVID-19 conjugates) and quality control antibody gold conjugates, 2) a nitrocellulose membrane strip containing two test bands (T1 and T2 bands) and a control band (C band). The T1 band is pre-coated with monoclonal anti-human IgG for the detection of IgG anti-COVID-19, T2 band is pre-coated with reagents for the detection of IgM anti-COVID-19 and the C band is pre-coated with quality control antibody.

When an adequate volume of test specimen is dispensed into the sample well of the cassette, the specimen migrates by capillary action across the cassette. COVID-19 IgM antibodies if present in the specimen will bind to the COVID-19 conjugates. The immunocomplex is then captured on the membrane by the pre-coated anti-human IgM antibody, forming a burgundy colored T2 band, indicating COVID-19 IgM positive test result.

COVID-19 IgG antibodies if present in the specimen will bind to the COVID-19 conjugates. The immunocomplex is then captured by the pre-coated reagents on the membrane, forming a burgundy colored T1 band, indicating a COVID-19 IgG positive test result. Absence of any test bands (T1 and T2) suggests a negative result. The test card also contains a quality control band C. Regardless of the presence or absence of a detection band, the red quality control band C should appear. The quality control band is a color band of the quality control antibody immune complex. If the quality control band C does not appear, the test result is invalid, and the sample needs to be tested again with another test card.

### Kit Components

- **25 Individually wrapped test cassette device(s)**
- **Disposable pipette(s)**
- **5 mL buffer**
- **Package insert**
- **Results Interpretation Quick Guide**

### Materials Not Provided : Timer

### Storage and Stability

Do not use after expiration date indicated on the package. Do not use the test if its foil pouch is damaged. Do not reuse test. The kit should be stored at 2~30°C in cool and dry place, protected from light. After opening the aluminum foil pouch, the test card will become invalid due to moisture absorption. Therefore, it is important that the test is performed and resulted within one hour of opening the individually wrapped foil cassette.

### Specimen Collection and Preparation

The kit can be performed using a whole blood (finger-stick) specimen [recommended], or plasma or serum samples prepared from commonly used anticoagulants (EDTA, heparin, sodium citrate). Testing should be performed immediately after specimen collection. **For Serum and Plasma Only:** If the test cannot be performed immediately, the serum and plasma specimen to be tested can be stored at 2 ~ 8 °C for 5 days. For long-term storage, store at -20 °C. Avoid repeated freeze-thaw specimens. Anti-coagulated whole blood specimens should not be stored for more than 72 hours at room temperature; not more than 7 days at 2-8 °C. Before testing, slowly return the refrigerated or frozen specimens to room temperature and mix them carefully. When clearly visible particulate matter is present in the specimen, it should be centrifuged to remove sediment before testing. If the specimen contains a large amount of lipid, hemolysis or turbidity, please do not use it, so as not to affect the result judgment.

### Assay Procedure

1. Place the test device on a clean, flat surface. After washing your hands, choose the non-dominant hand and face it palm side up. Remove the cap from the finger-stick device and use the disposable finger-stick device to stick the ring finger. It is recommended to wipe off the first droplet of blood with the provided gauze pad. **For Serum and Plasma Only:** Bring the specimen and test components to room temperature if refrigerated or frozen.

2. Fill the pipette dropper with the blood specimen. Holding the dropper vertically, dispense 1 drop (about 10  $\mu$ L) of whole blood (include finger blood), serum, plasma into the sample well, making

sure that there are no air bubbles. Then add 2 drops (about 70-100 µL) of Sample Diluent immediately.

3. Set up timer for 15 minutes. Read and record results at the 14-15 minute mark. It is important not to read results after 15 minutes.

### Interpretation of Assay Results

#### Negative Result



If only the C band is present, the absence of any burgundy color in the both test bands (T1 and T2) indicates that no COVID-19 antibody is detected in the specimen. The result is negative.

#### Presumptive Positive Result (Must be verified by HHS approved facility)



In addition to the presence of the C band, if the T1 band is visible, the test indicates the presence of the COVID-19 IgG antibody. The result is a presumptive positive and additional confirmation testing is immediately needed.



In addition to the presence of the C band, if the T2 band is visible, the test indicates the presence of the COVID-19 IgM antibody. The result is a presumptive positive and additional confirmation testing is immediately needed.



In addition to the presence of the C band, if both the T1 and T2 bands are visible, the test indicates the presence of both COVID IgG and IgM antibodies. The result is a presumptive positive and additional confirmation testing is immediately needed.

#### Invalid Result



If no C band is developed, the assay is invalid regardless of any burgundy color in the test bands. Repeat the assay with a new device.



### Performance Characteristics

The Pinnacle BioLabs COVID-19 Novel Coronavirus (COVID-19) Dual IgG/IgM Rapid Test is CE marked in the European Union and available in the USA under an Emergency Use Authorization use only in the United States.

1. Positive Coincidence Rate: The test results of positive quality control are all positive.
2. Negative Coincidence Rate: The test results of negative quality control are all negative.
3. Analytical Specificity: The test results of specimen from non-infected by novel coronavirus should be negative.
4. Analytical Sensitivity: The detection result is positive when detection of a novel coronavirus IgG strongly positive serum 1:50 dilution sensitivity reference : The detection result is positive when
5. Intra-Assay: There is no different test results of the same quality control in the same batch.
6. Inter-Assay: There is no different test results of the same quality.

### Limitations of the Test

1. The Assay Procedure and the Assay Result Interpretation must be followed strictly when testing. Failure to follow the procedure may give inaccurate results.
2. This kit is only used for in vitro diagnosis and is only used for qualitative detection of Coronavirus IgG and/or IgM antibodies in blood samples.
3. Positive and negative results indicate the presence of IgG and/or IgM antibodies with/without detectable concentrations of Coronavirus in blood samples, but cannot be used as the sole criterion for the determination of Coronavirus infection. Other methods (such as nucleic acid testing) should be used for identification when necessary, and comprehensive judgment should be made based on the test results. All positive results should be deemed presumptive positives and appropriate follow-up testing should be immediately sought.



Pinnacle BioLabs  
315 Deaderick Street  
15th Floor  
Nashville, TN 37238  
United States of America  
[www.PBLabs.com](http://www.PBLabs.com)  
1.800.609.6419

