

# ENDOZIME AW EXTRA

**Multi –Tiered Low-Foaming Enzymatic Cleaner  
Specifically Designed for Manual Cleaning,  
Ultrasonic and Washer Disinfectors**

Endozime AW Extra is an innovative and powerful formulation of an extra concentrated multi – tiered enzymatic detergent that quickly removes the soil contamination from all surgical instruments and flexible/rigid scopes in as little as one minute. It is suitable for hematological and compounding labs, tracheal tubes, suction tubes, biopsy forceps and flexible fiber optic instruments (gastroscopy, colonoscopies, bronchoscopes, etc.) as it polishes optical fibers.

The cleaner penetrates into the hardest-to –reach channels of scopes and cannulated instruments, instantly and actively breaking down gross contaminants upon contact and speeds the process of liquification and solubilisation. It is designed to clean the most difficult instruments (i.e. orthopaedic, laparoscopic) yet is safe for use on the most delicate (i.e. ophthalmic, microsurgical). Endozime AW Extra was developed for universal applications, eliminating the need for all other cleaners or detergents. Endozime Aw Extra’s result is not affected by the quality of water.

Endozime AW Extra has proven to be the most extra concentrated and cost effective cleaner for manual cleaning and all washer sterilizers and ultrasonic cleaners. It is bacteriostatic, non-abrasive, biodegradable, low-foaming and has a neutral pH (6-7.5). It is a non-corrosive, non-oxidising and non-caustic product. It is free rinsing and safe for all instruments and flexible/rigid scopes and will not harm any metals, plastic, rubber or corrugated tubing. It removes effectively the biofilm.

## **MANUAL USE, ULTRASOUND DEVICES & AUTOMATIC WASHING MACHINES:**

From 0.5 – 4 ml per litre of water depending on the soil level and cleaning method used.

Soak for one (1) minute. Rinse with tap, distilled or sterile water.

Discard used solution. Works in all temperatures.

## **COMPOSITION**

Isopropyl alcohol, Enzymes (Lipase, Amylase, Protease, Carbohydrase, proteolytic Enzyme, Hyaluronidase), Non-ionic detergents.

This product is Comply with the criteria as laid down in Regulation (EC) No. 648/2004 on detergents and has tested according to the EN ISO 15883. The information contained in this bulletin is, to our best knowledge, true and accurate, but all recommendations or suggestions are made without guarantees since the conditions of use are beyond our control.



393 Sagamore Avenue, Mineola, NY 11501 – 1919  
1-800-537-8463 • (516) 294-5888 • Fax (516) 248-6456  
Web Site: <http://www.ruhof.com> • E-mail: [ruhof.com](mailto:ruhof.com)