

Endozime® Xtreme Power

Item No. 34530

SUPER-CONCENTRATE - USE LESS CLEAN MORE

**Multi-Tiered Enzymatic Detergent with
Advanced Proteolytic Action and Rust Inhibitors**

**Designed for use in all Endoscope Washers,
Washer Disinfectors, Ultrasonic Washers and
For Manual Cleaning**



Endozime® Xtreme Power combines the power of a super-concentrate with the performance of Ruhof's proven multi-enzyme formula to create the ultimate surgical instrument and scope cleaner. This unique combination gives Endozime® Xtreme Power it's superior cleaning abilities to deliver fast and thorough soil contamination removal. Enhanced with APA (Advanced Proteolytic Action) that greatly increases protease enzyme activity. Endozime® Xtreme Power penetrates into the hardest-to-reach channels of scopes and cannulated instruments, instantly and actively breaking down gross contaminants upon contact. Endozime® Xtreme Power contains biological additives that speed the process of liquification and solubilization, facilitating enzymatic action and contributing to the product's overall effectiveness.

Features

- **Super-Concentrate:** Dilution from 1/4 oz. per gallon (2ml/liter)
- Effective in all water temperatures (Including cold water)
- No foam formula
- Removes all soil contamination - blood, protein, fats, carbohydrates and starches
- Available in smaller bottles - easier to carry, store and pour
- Compatible with all medical device materials
- Neutral pH, Free rinsing and 100% Biodegradable
- Inhibits rust

Environmental and Cost Benefits:

- Less detergent usage
- 100% biodegradable
- Less plastic bottle and packaging waste
- Reduction in energy usage

Directions for Use

Endozime® Xtreme Power can be diluted as low as 1/8 oz. (1ml/liter) per (1) gallon of water up to 1/2 oz. (4ml/liter) per one (1) gallon of water, but for optimal results, the following is recommended.

Manual Cleaning

1. Add 1/4 oz. per gallon (2ml per liter) of Endozime® Xtreme Power (Works in all water temperatures)
2. Submerge instruments and scopes to be cleaned.
(For scopes: suction or flush through channels before soaking)
3. Soak for two minutes to remove all organic matter. (Heavy soil may require longer time)
4. Rinse thoroughly with tap, distilled or sterile water.
5. Discard used solution.

Automatic Cleaning Devices:

1. Pour 1/4 oz. per gallon (2 ml/liter) directly into detergent reservoir or use full strength.
 2. Calibrate the washer dilution at 1/4 oz. per gallon (2 ml/liter) of water
- Follow manufacturer's recommendations for correct dosing option.

Cautions: Avoid prolonged contact with skin. Do not swallow. Keep away from children.