

# CS/SPD PULL THRU™

## Single Use Cannulated Instrument Cleaning Device

ITEM #	DESCRIPTION	WAND COLOR
345PT1	MICRO: 1.2mm for lumen diameters of 1-1.2mm	Purple
345PT2	MINI: 2.6mm for lumen diameters of 1.4-2.6mm	Green


ITEM #	DESCRIPTION	WIPER COLOR
345PT5	5mm for lumen diameters of 3-5mm	White
345PT8	8mm for lumen diameters of 6-8mm	Green
345PT12	12mm for lumen diameters of 9-12mm	Blue
345PT15	15mm for lumen diameters of 13-15mm	Yellow

The CS/SPD PULL THRU™ is an instrument cleaning device, which can be effectively used to clean the hard-to-reach narrow lumen channels of laparoscopic and cannulated instruments in a single pass. The CS/SPD PULL THRU™ is designed to clean instruments with lumens from 1mm to 15mm. Use with Ruhof's multi-tiered enzymatic detergents for more effective cleaning.



The unique wiper design of the CS/SPD PULL THRU™ provides a complete circumferential seal in the lumen channel, thereby removing almost all residue in a single pass. **The 360° seal creates a vacuum which draws detergent through the lumen removing residue from crevices or other areas of minor damage in the lumen. The vacuum also ensures the lumen is completely filled with detergent to attack and remove bioburden.**

### Features

- Effectively cleans lumen channels from 1mm-15mm with a single pass.
- Unique wiper element design provides a complete circumferential seal in the lumen with a piston like action thereby removing almost all residue.
- Reduces the time required to manually clean cannulated instruments.
- Removes the residue bristle brushes leave behind, enabling detergents to attack and remove bioburden more effectively.
- CS/SPD PULL THRU™ action is consistent every time, which reduces any variation in cleaning method between users.
- Disposable/single use device, which removes the need for the complicated and expensive cleaning of traditional brushes.
- Significantly improves cleaning efficiency when compared to traditional brushes.
- Wipers are constructed of soft plastic, non-abrasive material and will not harm the lumen walls.



The multiple wiper design of the CS/SPD PULL THRU™ provides a complete circumferential seal in the lumen channel, thereby removing almost all residue in a single pass.



The CS/SPD PULL THRU™ is available in 6 sizes, each is color coded to clearly indicate the wiper circumference.



The CS/SPD PULL THRU™'s 61cm flexible wand allows for easy access through all lumen channels.


continued


# ScopeValet™ CS/SPD PULL THRU™

## Directions for Use

**Working Lumen: 1mm-15mm ID**

**Working Length: approximately 61cm**

1. Be sure to use the CS/SPD PULL THRU™  on wet lumen channels.
2. Open the pouch and remove the CS/SPD PULL THRU™ lumen cleaning device from the package.
3. Check that you are using the correct size CS/SPD PULL THRU™ by inserting the CS/SPD PULL THRU™ head ② into the distal end of the instrument to be cleaned; if too tight or too loose check channel inner diameter against instrument manufacturer's instructions.
4. Insert the short tip end ① of the CS/SPD PULL THRU™ into the proximal end of the instrument to be cleaned. Advance the device through the lumen and exit at the distal end, slowly pulling the device completely through.
5. Check for soil on wipers 3 and 4. If there is soil, rinse and repeat step 4.
6. Dispose of device after use.



393 Sagamore Avenue, Mineola, New York 11501  
1-800-537-8463 • [www.ruhof.com](http://www.ruhof.com)

