

ASSEMBLY INSTRUCTIONS FOR FOLDING PARALLEL BARS MODELS 3-3317 & 3-3310

MAXIMUM WEIGHT CAPACITY 350 LBS

INTENDED USE. Clinton Folding Parallel Bars, models 3-3317 and 3-3310 are specifically designed and engineered for physical therapy use only. Persons using the equipment should be observed by trained medical or physical therapy personnel. Never allow anyone to use the equipment unattended. Combination of weight and force exerted on the equipment is not to exceed a total of 350 lbs. (158 kg). Using this equipment for unintended purposes could result in damage to the equipment.



Expanded for Use



Folded for Storage



Unpack and expand parallel bars and place on a solid, level, flat surface. Insert a star knob in the hole on each of the four corners.





When fully expanded, insert the two, chained, screw pins as shown and tighten.

TO FOLD: Remove pins and compress unit, sliding the sides in opposite directions until the bases align.



Clinton Industries warrants its products to be free of defects in workmanship for five years from the date of purchase. This warranty does not include damage resulting from accident, abuse or misuse, alteration, or negligence of the product and expressly excludes normal wear of parts. Contact Clinton for full warranty details.

CLEANING

Clinton Folding Parallel Bars are built to be durable. With proper care, your Clinton equipment will provide years of superior looks and performance. Before using any cleaning agents, try them on a small inconspicuous area to insure that they will not harm the surface finish. Damage from cleaning is not covered under the warranty. **Do not use alcohol based cleaning agents.**

Periodically check that all parts are secure and operating properly.