## Instructions For Use

# DIOP8

## Single Use 8 MHz Surgical Doppler Probe

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Instructions for Use · Mode d'emploi · Gebrauchsanleitung · Istruzioni per l'uso · Instrucciones de uso  
Bruksanvisning · Gebruiksaanwijzing · Brugervejledning · Bruksanvisning · Käyttöohjeet ·

Single Use 8 MHz Surgical Doppler Probe · Sonde chirurgicale Doppler 8 MHz à usage unique · Chirurgische 8MHz  
Doppler-Sonden zum Einmalgebrauch · Sonda Doppler quirúrgica monouso, 8 MHz · Sonda Doppler quirúrgica de  
8 MHz para un solo uso · 8 MHz kirurgisk dopplersond för engångsbruk · 8MHz Chirurgische Dopplersonde voor  
eenmalig gebruik · Kirurgisk 8 MHz Doppler-probe til engangsbrug · 8 MHz kirurgisk dopplersonde for engangsbruk ·  
Kertakäyttöinen 8 MHz:n kirurginen doppler-anturi ·

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

DIPP10 - DIOP8 x 10

## Indications For Use

The DIOP8 surgical Doppler is indicated for use by qualified healthcare practitioners in a sterile condition in the operating theatre for the assessment of vascular blood flow by direct application to the vessel wall.

## Contraindications

- Do not use the DIOP8 on the eye.
- Do not use in the presence of flammable gases or oxygen rich environments.
- Do not apply the DIOP8 to the patient when using high frequency (HF) surgical equipment. Make sure when using such equipment that the DIOP8 is not in contact with the patient.

**Note. No pressure should be applied to the vessel wall.**

## Warnings/Cautions & Safety

- The sterile Intraoperative Probes are delicate and should be handled with care. Do not drop or strike against hard surfaces. Avoid excessive tension on the probe cable.**
- These probes are supplied STERILE. Always ensure that the packaging is intact and undamaged before opening.**
- The DIOP8 is a screening tool to aid the healthcare professional and can not provide a diagnosis. If there is doubt as to the status of the vessel after using the unit, further investigations should be undertaken immediately using alternative techniques.**
- We recommend that exposure to ultrasound should be kept As Low As Reasonably Achievable - (ALARA guidelines). This is considered to be good practice and should be observed at all times.**
- The DIOP8 is for single use only. DO NOT REUSE. DO NOT STERILISE. Failure to follow this advice may result in failure of the product, compromised sterility and potential harm to the patient.**
- This equipment is for use only by suitably qualified healthcare practitioners.**
- This equipment must not be modified.**
- Note. The control unit and adaptor are NOT sterile. Refer to the hand-held Doppler IFU for details of cleaning.**

## Operation

Connect the probe adaptor to the Dopplex cable, using the alignment mark on the barrel. Locate the probe adaptor into the probe holder at the side of the Dopplex control unit. The Dopplex control unit can then be mounted onto an IV pole using the clamp as shown in Fig.1. Using aseptic procedure, remove the DIOP8 probe from its packaging and insert plug into adaptor. Always maintain sterility of the probe.

**Note: The Dopplex control unit will automatically switch off at a predetermined time after switching the unit on. See control unit IFU for more information.**

The Dopplex control unit can be switched on by a member of staff positioned outside the sterile field.

Adjust the volume control accordingly.

Ensure that the probe connector is fully engaged in the socket on the probe adaptor.

Correct operation of the system should be confirmed by placing the DIOP8 tip wetted with patient's body fluid in light contact with an artery which is known to have blood flowing through it. A clearly audible pulsatile Doppler sound should be heard. The audibility and amplitude of the signal will be optimised when the probe is at an angle of approximately  $45^\circ \pm 15^\circ$  to the vessel.

The probe can now be used to assess blood flow in other vessels.

Aiming the probe proximally along the line of the vessel at an angle of  $45^\circ \pm 15^\circ$  and ensuring that the tip is fully wetted at all times (Fig.2.) will enable the optimal audio signal to be obtained. (If you are using a DMX, Multi Dopplex II or Super Dopplex II, an indication of the blood velocity and its direction will be shown on the LCD display.)

## Operation

The probe can then be moved along the length of the vessel, noting any change in pitch of the Doppler signal or height of the Doppler waveform. This may be indicative of a change in the lumen area.

The same procedure can then be carried out after a graft has been inserted to confirm that adequate blood flow has been restored. By placing the probe on the vessel distal to the anastomosis, confirmation of distal run off is provided.

## Product Labelling

	This symbol signifies that this product, including its accessories and consumables is subject to the WEEE (Waste Electrical and Electronic Equipment) regulations and should be disposed of responsibly in accordance with local procedures.
	This symbol signifies that this product complies with the essential requirements of the Medical Device Directive (93/42/EEC) - Medical Device Regulation (EU/2017/745)
	Legal Manufacturer in association with the CE mark in Europe ArjoHuntleigh AB, Hans Michelsensgatan 10 211 20 Malmö, Sweden
	3 year from date of manufacture.
<b>Manufactured by:</b>	<b>Huntleigh Healthcare Ltd.</b> 35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk www.huntleigh-diagnostics.com
	Sterile Barrier System (non sterile protective outer layer with an internal sterile barrier)
	Attention, consult accompanying documents / Instructions for Use
	Warning
	Medical Device
	Device Identification code
	Contents can be recycled
	Temperature Limitations
	Rx Only
	Single patient use
	Do not resternalise
	Batch Code
	Do not use if packaging is damaged
	Use By date
	Keep Dry
	Do Not Reuse
	Sterilised using ethylene oxide
	Date of Manufacture
	Atmospheric Pressure Limitation
	Fragile, handle with care

## Technical Specifications - IEC60601-1 Classification

Type of shock protection	Internally powered equipment (Dopplex control unit)
Degree of shock protection	Type CF equipment
Protection against ingress of liquids	<b>Probe Adaptor:</b> Ordinary Equipment <b>DIOP8:</b> Suitable for use in contact with body fluids.
Degree of safety in presence of flammable gases	Equipment not suitable for use in presence of flammable gases or oxygen rich environments
Mode of Operation	Continuous

## Ultrasound

$P_r$	< 1 MPa
$I_{ob}$	< 20 mW/cm <sup>2</sup>
$I_{sp1a,3}$	< 100 mW/cm <sup>2</sup>
Thermal Indices and Mechanical Index are 1.0 or less for all device settings.	

## Standards

<b>Complies With</b>	IEC60601-1 :2012
<b>EMC</b>	IEC60601-1-2: 2014
<b>CISPR11</b>	Group 1 Class A

## Environmental

<b>Operating</b>	
Temperature range	+10 °C to +30 °C
Relative Humidity	10% to 90% (non-condensing)
Pressure	860 to 1060 hPa
<b>Storage</b>	
Temperature range	-10 °C to 40 °C
Relative Humidity	90% maximum
Pressure	860 to 1060 hPa


Fig. 1


Fig. 2

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

Kuva 2

## Warranty and Service

Huntleigh Healthcare's standard terms and conditions apply to all sales. A copy is available on request. These contain full details of warranty terms and do not limit the statutory rights of the consumer.

### Service Returns

There are NO USER SERVICEABLE PARTS inside the probe adaptor or probe. If for any reason your DIOP8 or probe adaptor is being returned, please:

1. Contact Service Dept. to obtain instructions and authorisation for the product to be returned.
2. Failure to do this may result in the product being returned without investigation.

EN

### FOR UK

For further details refer to the NHS document HSG(93)26.

For service, maintenance and any questions regarding this, or any other Huntleigh Healthcare Dopplex product, please contact:

Huntleigh Healthcare Ltd, Diagnostic Products Division, 35 Portmanmoor Road, Cardiff, CF24 5HN UK, Tel : +44 (0)29 20485885, Fax: +44 (0)29 20492520  
Or your local distributor.

### FOR USA

For service, maintenance and any questions regarding this, or any other Huntleigh Healthcare Dopplex product, please contact:

Service Department, ArjoHuntleigh Inc, Suite 250, 2349 West Lake Street, Addison, IL60101, USA T: 1-800-323-1245 option 2 W: www.huntleigh-healthcare.us  
Or your local distributor.

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FI

This section is only applicable to United Kingdom (UK) market when UK marking is applied to the Arjo medical device labelling.

### UK Symbol:

UK marking indicating conformity with UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended)  
Figures indicate UK Approval Body supervision.

### UK Responsible Person:

Arjo (UK) Ltd., ArjoHuntleigh House, Houghton Regis. LU5 5XF

Is the appointed UK Responsible Person as defined in UK Medical Devices Regulations 2002 (SI 2002 No 618, as amended).

For Northern Ireland (NI) CE marking will still apply until further amendment to applicable regulations.

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If a serious incident occurs in relation to this medical device, affecting the user, or the patient then the user or patient should report the serious incident to the medical device manufacturer or the distributor.  
In the European Union, the user should also report the serious incident to the Competent Authority in the member state where they are located.

**Manufactured in the UK by Huntleigh Healthcare Ltd on behalf of;**

 ArjoHuntleigh AB  
Hans Michelsensgatan 10  
211 20 Malmö, Sweden

**Huntleigh Healthcare Ltd.**

35 Portmanmoor Road, Cardiff, CF24 5HN, United Kingdom  
T: +44 (0)29 20485885 sales@huntleigh-diagnostics.co.uk  
www.huntleigh-diagnostics.com

Distributed in the USA by:  
**ArjoHuntleigh, Inc.**  
2349 West Lake Street, Suite 250  
Addison, IL 60101  
T: 800-323-1245  
www.huntleigh-healthcare.us

1001057-2

Vertrieben in Deutschland Durch:  
Huntleigh Healthcare GmbH:  
Industriering Ost 66  
47906 Kempen  
Germany  
T: +49 02152 551110  
F: +49 02152 551120  
https://www.huntleigh.de  
1001058-2

www.huntleigh-diagnostics.com/

Registered No: 942245 England & Wales. Registered Office:  
ArjoHuntleigh House, Houghton Hall Business Park, Houghton Regis, Bedfordshire, LU5 5XF  
©Huntleigh Healthcare Limited 2003-2020

**A Member of the Arjo Family**

® and ™ are trademarks of Huntleigh Technology Limited

As our policy is one of continuous improvement, we reserve the right to modify designs without prior notice.

www.huntleigh-healthcare.us/

# HUNTLEIGH

1001048-3