Midmark 630 HUMANFORM[®] Procedures Chair

[English](#)
[Français](#)
[Español](#)

For Models:

630 (-010 /-011 /-012 /-013)
630 (-020 /-021 /-022 /-023)

**630 (-010 /-012 /-020 /-022)
Models Only**

User Guide

Style P

003-10431-99 Rev. AA2 (4/19/21)

TP202 20-42-FO-00014 Rev A2 C2169


Product Information

Dealer:

Date of Purchase:

Model / Serial Number:

Midmark Authorized Service Company:

Product Registration

To register your product warranty go to www.Midmark.com.

Table of Contents

Important Information

Safety Symbols.....	4
Intended Use.....	5
Electrical Requirements.....	5
Electromagnetic Interference.....	5
Transportation / Storage Conditions.....	6
Disposal of Equipment.....	6
Authorized Representatives.....	6
Safety Instructions.....	8
EMC - Manufacturer's Declaration and Guidance.....	9

Operation

Hand / Foot Control Connection.....	10
Base / Back / Foot / Tilt Functions.....	11
Active Sensing Technology™.....	15
Activating Multiple Functions.....	16
Home Function.....	17
Quick Chair® Function.....	18
Control Lockout.....	19
Association Procedure.....	20
Foot Extension.....	21
Headrest Positioning.....	22
Treatment Pan.....	23
Paper Roll.....	24
Stirrups.....	25
Duplex Receptacle.....	26
Programming Chair Positions.....	27
One-Touch Positioning.....	28
Rotational Base Operation (optional).....	29

Maintenance

Calling for Service.....	30
Cleaning.....	30
Preventative Maintenance.....	30
Troubleshooting.....	31

Specifications

Model Identification / Compliance Chart.....	33
Specifications.....	34
Range of Motion.....	35
Dimensions.....	36

Warranty Information

Limited Warranty.....	39
-----------------------	----

Important Information - Safety Symbols

WARNING

Indicates a potentially hazardous situation which could result in serious injury.

Caution

Indicates a potentially hazardous situation which may result in minor or moderate injury. It may also be used to alert against unsafe practices.

Equipment Alert

Indicates a potentially hazardous situation which could result in equipment damage.

Note

Amplifies a procedure, practice, or condition.

Manufacturer

Duty Cycle (motor run time)
30 seconds ON, 5 minutes OFF

Proper shipping orientation

Fragile

Type B, Applied Part
(Upholstery)

Protective earth ground

Patient Weight Limit

Catalogue Number

Serial Number


No Pushing

Keep dry

Maximum stacking height
(palletted units)

Fuse rating specification

Dangerous voltage / shock hazard

Refer to manual

Intended Use

This product is intended to be used as a chair / table to provide positioning and support of patients during general examinations and procedures conducted by qualified medical professionals in a medical office environment.

Electrical Requirements

Caution

To completely isolate the chair from electrical mains supply, power cord must be unplugged.

When using high frequency surgical devices or endocardial catheters, use non-conductive material to insulate patient from metal portions of the chair.

Failure to comply may result in electric shock or burns to patient.

Electromagnetic Interference

This product is designed and built to minimize electromagnetic interference with other devices. However, if interference is noticed between another device and this product:

- Remove interfering device from room
- Plug chair into isolated circuit
- Increase separation between chair and interfering device
- Contact Midmark if interference persists

Transportation / Storage Conditions

Ambient Temperature Range:.....-30°C to +60°C (-22°F to +140°F)
Relative Humidity.....10% to 90% (non-condensing)
Atmospheric Pressure500hPa to 1060hPa (0.49atm to 1.05atm)


Operating Conditions

Ambient Temperature Range:.....+10°C to +40°C (+50°F to +104°F)
Relative Humidity.....30% to 75% (non-condensing)
Atmospheric Pressure700hPa to 1060hPa (0.69atm to 1.05atm)

Disposal of Equipment

At the end of product life, the chair, accessories, and other consumable goods may have become contaminated due to normal medical use. Consult local codes and ordinances for proper disposal of equipment, accessories and other consumable goods.

Authorized Representatives

WARNING

In the event of any serious incident in relation to the device, please contact Midmark and the appropriate competent authority.

Customers in the EU should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

CEpartner4U
Esdoornlaan 13
3951 DB Maarn, The Netherlands
Phone: +31 343 442 524
Fax: +31 343 442 162

Customers in Australia should direct all questions, incidents, and complaints to Midmark's Sponsor listed below.

ICONA Pty Ltd.
1A/2A Westall Road
Clayton, Victoria 3168
Phone: 1300 442 662

Authorized Representatives - continued

Customers in India should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Midmark India
Kalpataru Point, Unit No. 12, 1st Floor Kamani road, Sion (E)
Mumbai, 400022
Tel: +91 22 4915 3000
Fax: +91 22 4915 3100

Customers in Saudi Arabia should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Attieh Medico Ltd
Abdullah M. Al Khalifi Street
Al Nakheel Dist. II
P.O. Box 116105
Jeddah (East)
21391
Saudi Arabia
Tel: + 966 2 286 4707
Fax: + 966 2 286 4744

Customers in the United Arab Emirates should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

GULF DRUG LLC.
Al Barsha,
Dubai, UAE
Phone: +971 4 501 4000
Fax: +971 4 501 4100

Customers in Hong Kong should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Associated Medical Supplies Company LTD. 沒有中文版使用說明
Room 1201-1202, 12/F, Fo Tan Industrial Center
26-28 Au Pui Wan Street, Fo Tan, Shatin
New Territories, Hong Kong
Phone: +852 2604 9389
fax: +852 2694 0866

Customers in Qatar should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Tadmur Trading
Midmac Roundabout
Opposite Al Mana Petrol Station
Doha, Qatar
Tel: +974.44337000
Fax: +974.44337100

Customers in Israel should direct all questions, incidents, and complaints to Midmark's Authorized Representative listed below.

Medtechnica Ltd.
HaTnufa St 7
Petah Tikva, Israel
Tel: +972 3.925.4040
Fax: +972 3.924.9977

Safety Instructions

WARNING

No modification of this equipment is allowed.

WARNING

Connecting equipment to the multiple socket outlet effectively leads to creating a Medical Electrical System and the result can be a reduced level of safety.

WARNING

Chair must be positioned in normal use such that the power cords can be easily accessed to unplug.

WARNING

To avoid risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

WARNING

Equipment is not suitable for use in the presence of a flammable anesthetic mixture.

Clarification: *Equipment is suitable for use in the presence of oxygen, air, or nitrous oxide.*

WARNING

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

WARNING

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Procedures Chair including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

EMC - Manufacturer's Declaration and Guidance

Note

The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals. The exam chair is not intended to be used in a residential environment.



Immunity Test	Immunity Test Level
Electrostatic Discharge	± 8 kV Contact, ±15 kV Air
Radiated RF EM Fields	3V/m 80MHz-2.7 GHz 80% AM at 1kHz
Rated Power Frequency Magnetic Fields	30 A/m
Electrical Fast Transients	±2 kV at 100kHz repetition frequency
Surges	Line to Line: ±1 kV Line to Ground: ±2 kV
Conducted RF	3V: 0.15MHz -80MHz 6V in ISM Bands between 0.15MHz - 80MHz 80% AM at 1kHz
Voltage Dips	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°
Voltage Interruptions	0 % UT; 250/300 cycle

Radiated Test	Radiated Test Level
Conducted and Radiated RF Emissions	CISPR 11 Group 1 Class A

Operation

Hand / Foot Control Connection

Although wireless controls are standard on some models, there are three hand / foot control ports for wired controls. One is located on the base column at the foot end of the chair, and two are located on the PC board cover at the head end. These control ports are interchangeable.

MA8830i

Base Up / Down



Attention!

The Active Sensing Technology™ will stop and retract the Base Up / Down function if the foot section contacts another object.

(You will hear a series of “beeps” if the Active Sensing Technology™ is activated.)


To resume Base Up / Down operation:
 Remove object from under the foot section.
 Check if foot extension is extended.
 Press Base Down button.

To move the Base Up or Down...
 Press and hold the appropriate button.

Movement stops when...
 Up / Down button is released.
 Maximum / minimum position is reached.
 Active Sensing Technology™ is activated.
 Stop button is pressed.

Models w/Rotational Base

Max.
 111 cm
 (43.5 in.)
Min.
 55 cm
 (21.5 in.)


Max.
 102 cm
 (40 in.)

Min.
 46 cm
 (18 in.)

MA8865i


Back Up / Down

To move the Back section Up or Down...
Press and hold the appropriate button.

Movement stops when...

*Up / Down button is released.
Maximum / minimum position is reached.
Active Sensing Technology™ is activated.
Stop button is pressed.*

MA8866i

Foot Up / Down

Attention!

The Foot Up and Down functions will be disabled when the treatment pan is not in the stowed position.

The Active Sensing Technology™ will stop and retract the Foot Up / Down function if the foot section contacts another object.

(You will hear a series of “beeps” if the Active Sensing Technology™ is activated.)

To resume Foot operation:

Push the treatment pan all the way in.



Remove object from under the foot section.

Check if foot extension is extended.

Press Foot Up / Down button.

To move the Foot section Up or Down...

Press and hold the appropriate button.

Movement stops when...

Up / Down button is released.

Maximum / minimum position is reached.

Treatment pan is not in the stowed position.

Active Sensing Technology™ is activated.

Stop button is pressed.

MA88671

Tilt Up / Down

Attention!

The Active Sensing Technology™ will stop and retract the Tilt Up / Down function if the top of the column or the foot section contacts another object.

(You will hear a series of “beeps” if the Active Sensing Technology™ is activated.)

To resume Tilt Up / Down operation:

Remove object from the top of the column.

Remove object from under the foot section.

Check if foot extension is extended.

Press Tilt Down button.


To Tilt the Seat section Up or Down...

Press and hold the appropriate button.

Down

Up

Up

Down

Stop


Movement stops when...

Up / Down button is released.

Maximum / minimum position is reached.

Active Sensing Technology™ is activated.

Stop button is pressed.

Max.
30°

Min.
0°

MA8868i


Active Sensing Technology™


Caution

Be sure that all personnel and equipment are clear of the chair before activating any function. Failure to do so could result in personal injury.


The Active Sensing Technology™ prevents the headrest from moving within four inches of the floor. If the Active Sensing Technology™ is activated, it will sound a series of “beeps” and disable the following functions: Base Down, Back Down, and Tilt Up.

The Active Sensing Technology™ prevents the foot section from moving within four inches of the floor when the foot extension is extended. If the Active Sensing Technology™ is activated, it will sound a series of “beeps” and disable the following functions: Base Down, Tilt Down, and Foot Down.

The Active Sensing Technology™ prevents the back section from moving within one half inch of the PC Board cover. If the Active Sensing Technology™ is activated, it will sound a series of “beeps” and disable the following functions: Base Down, Back Down, and Tilt Up.


Activating Multiple Functions

The *BASE / BACK / TILT / FOOT* functions may be activated at the same time.

Example:

Base Up / Back Down

MA8869i

To activate multiple functions ...

*Simultaneously press & hold the desired buttons.
(ex. Base Up & Back Down)*



Home Function


The Home function automatically lowers only the chair base so that your patient may get on or off of the chair.

To lower the base to 22 in. (55.9 cm)...
 Press, then immediately release the Home button.

To lower the base to 18 in. (45.7 cm)...
 Press and hold the Home button until you hear a single “beep”.

Movement stops when...
 “Home” position is reached.
 Active Sensing Technology™
 is activated.
 Stop button is pressed.

MA8870i

Quick Chair® Function (w/ Foot Control only)


The Quick Chair® function activates the Base Down, Tilt Down, Foot Down, & Back Up functions at the same time.

To activate the Quick Chair® function...
Press and release the Quick Chair® button.

Quick Chair®
Button

Movement stops when...
“Chair” position is reached. Active Sensing Technology™ is activated. Stop button is pressed.


MA8871i

Control Lockout

Control Lockout:

Press and hold the Stop and Back Down buttons simultaneously for two seconds. You will hear a single “Beep” to indicate control lockout is enabled / disabled.


Note: While the controls are locked out, you will hear a series of three “Beeps” while motion is attempted.

Association Procedure

Association Procedure:

- A) Disconnect power to the chair for at least 5 seconds.
- B) Reconnect power, wait for at least two seconds, then press & hold the **Stop & Home** buttons on the **wireless controller** until you hear three “beeps”.
- C) Check for proper operation of wireless controller.
- D) If associating a hand **and** foot control, repeat procedure using the other wireless controller.

Foot Extension

The foot section can be extended up to 3.5 in. (8.9 cm) to accommodate taller patients.


Attention!

The Active Sensing Technology™ prevents the foot section from moving within four inches of the floor when the foot extension is extended. If the Active Sensing Technology™ is activated, it will sound several “beeps” and disable the following functions: Base Down, Tilt Down, and Foot Down.

To resume Base Down, Tilt Down, and Foot Down operation:
Move foot extension in.
Press Base Down, Tilt Down or Foot Down buttons.

To extend / retract the foot section...

- A) Squeeze handle on either side of foot section.
- B) Move foot section in / out to desired position.
- C) Release handle to lock in position.

MA8872i

Headrest Positioning


Caution

Headrest must be securely locked in place before performing any exam or procedure.

To adjust the pivot angle...

- A) Loosen locking knob.
- B) Pivot headrest to desired position.
- C) Tighten locking knob to lock in place.

MA8877i

Treatment Pan

Attention!

The Foot Up and Down functions will be disabled when the treatment pan is not in the stowed position.

(You will hear a series of “beeps” if the treatment pan is not in the stowed position when foot Up / Down is activated.)


To resume Foot Up / Down operation:

Push treatment pan all the way in.

Press Foot Up or Foot Down button.

To access the treatment pan...

- A) Lower the foot section.
- B) Pull the treatment pan outward.


MA8875i

Paper Roll

Paper Roll Size (max.)..... 45.7 cm x 8.9 cm diameter
(18 in. long x 3.5 in.)

To install the paper roll...

- A) Pivot paper roll cover back to open.
- B) Install paper roll, then pull paper over chair.
- C) Pivot paper roll cover forward to close.

MA8873i

Stirrups


Caution

Be sure the stirrups are locked in place before using. The stirrups will not support the patient's entire weight / excessive force.

To position the stirrups...


- A) *Pull the stirrup out, then unfold.*
- B) *Lift the stirrup slightly, then move it left or right as desired.*
- C) *Release stirrup to lock in position.*

Duplex Receptacle (not available on models with rotational base.)

Duplex receptacles on each side of the chair provide power for accessories used during medical procedures. There is a circuit breaker located at the base of the chair. If the receptacle's maximum load is exceeded, the circuit breaker interrupts power to the receptacles.

Maximum Load [both duplex receptacles combined]..... 115 VAC, 3 amps, 50/60 Hz

MA8876i

Programming Chair Positions (Programmable Models only)

The Programming feature allows you to easily recall frequently used chair positions.

To program a chair position...


- A) Move the chair to desired position.
- B) Press the Program button. (you will hear a single “beep”)
- C) Press the desired Position button (1, 2, 3, or 4).
(you will hear three “beeps”)

Note: You must press the desired Position button within 5 seconds of pressing the Program button.

To recall a programmed position...

Press and release the desired Position button (1, 2, 3, or 4).

Note: The Position buttons (1, 2, 3, 4) can be changed to function as “press & hold”. Refer to: **One-Touch Positioning** for procedure to enable.

Movement stops when...

- Programmed position is reached.
 - Position button is released.
 - Active Sensing Technology™ is activated.
 - Stop button* is pressed.
- (* or any hand / foot control button)

One-Touch Positioning (Programmable Models only)


During normal operating mode, the One-Touch Positioning feature is activated. The Position buttons function as “press & release”. When disabled, Position buttons must be pressed & held until the programmed position is reached.

To enable (or disable) the One-Touch Positioning feature...
Press & hold the Stop & Program buttons simultaneously. (after two seconds, you will hear a single “beep” if successful)

Note: The One-Touch Position feature is preset (enabled) at factory.

To recall a programmed position...
Press and release the desired Position button (1, 2, 3, 4).

Note: You will hear a single “beep” and the chair will move to the programmed position.

Movement stops when...
Programmed position is reached.
Active Sensing Technology™ is activated.
Stop button* is pressed.
(* or any hand / foot control button)

Rotational Base Operation (optional)

Note

Chair motion is not allowed while rotation is activated.


To rotate the chair base...

- A) Press and release the brake pedal.
- B) Rotate the chair to desired position.

To lock in position...

- A) Press and release the brake pedal.

Note: The base will automatically lock in position after one minute.

Maintenance

Calling for Service

Note

Model / serial number information is required when calling for service.

If service is required, contact your authorized Midmark dealer.

To contact Midmark directly:

1.800.MIDMARK (1.800.643.6275)

8:00 am until 5:00 pm Monday - Friday (ET)

www.midmark.com

Cleaning

Upholstery

Equipment Alert

The upholstery is resistant to most medicinal-type stains, but may be damaged by solvents and dyes. Immediately remove any fluids spilled on the upholstery.

Wash your upholstery weekly with a mild liquid soap and water mixture, rinse with clear water and dry completely to remove disinfection cleaner build-up.

Disinfect your upholstery with a solution of standard bleach and water mixed 1 in 10 (10%) or chlorine based cleaners. Follow this with a clear water rinse and thorough drying of material. See current CDC Guideline for Disinfection & Sterilization in Healthcare Facilities.

To minimize damage caused by disinfectant cleaner residue build-up, do not allow disinfectants to pool on the upholstery surface. Once the approved contact time has been obtained, remove and dry excess liquid remaining on the surface.

Detailed care and maintenance instructions are included with your product. This information is also available on www.midmark.com in the Technical library under the User Information tab for your product.

Painted Metal / Plastic Surfaces

Clean the painted metal and plastic surfaces weekly using a clean soft cloth, and mild cleaner.

Preventative Maintenance

Periodically inspect the following areas:

- Power cord(s) should be free of cuts or other visible damage.
- All fasteners should be in place and tightened securely.
- All mechanical functions should operate properly.

Have a qualified service technician inspect your equipment every twelve months.

Troubleshooting

Symptom	Probable Cause	Correction
No functions will operate.	Facility supply voltage.	Secure power cord connections at wall outlet and at chair base.
		Reset facility circuit breaker.
No power at chair receptacles. All other functions work.	Chair circuit breaker tripped.	Reset chair circuit breaker located at the base of chair, press to reset.
When any control button is pressed, that controller sounds a single beep. (All functions operate)	Low battery. (Wireless controls only)	Replace batteries. (size: AA)
When any control button is pressed, that controller sounds a single beep. (No functions operate)	Controller not associated with chair. (Wireless controls only)	Perform association procedure.
When any control button is pressed, that controller sounds a series of three beeps. (No functions operate)	Control Lockout is enabled.	Disable Control Lockout.
Base Down function will not operate.	Active Sensing Technology™ activated.	Activate Tilt Up / Down and Back Up, then try Base Down again.
	Foot extension is extended.	Remove object from under foot section.
Base Up function will not operate.	Active Sensing Technology™ activated.	Move foot extension in or raise foot section.
Back Down function will not operate.	Active Sensing Technology™ activated.	Remove object from under foot section.
		Activate Base Up and Tilt Down, then try Back Down again.

Troubleshooting (continued)

Symptom	Probable Cause	Correction
Tilt Up function will not operate.	Active Sensing Technology™ activated.	Activate Base Up and Back Up, then try Tilt Up again.
		Remove object from the top of the column.
		Remove object from under foot section.
Tilt Down function will not operate.	Active Sensing Technology™ activated.	Activate Base Up and Foot Up, then try Tilt Down again.
		Remove object from the top of the column.
	Foot extension is extended.	Remove object from under foot section.
Foot Up function will not operate.	Treatment pan not in stowed position.	Move foot extension in or raise foot section.
	Active Sensing Technology™ activated.	Push treatment pan all the way in.
Foot Down function will not operate	Treatment pan not in stowed position.	Remove object from under foot section.
	Foot extension is extended.	Push treatment pan all the way in.
	Active Sensing Technology™ activated.	Move foot extension in.
"Home" function does not lower the chair to a height of 18 inches. (Rotation Models 21.5 inches)	"Home" button pressed, then released.	Remove object from under foot section.
	Foot extension is extended.	Activate Tilt Up / Down and Back Up, then try "Home" again.
	Active Sensing Technology™ activated.	Remove object from under foot section.
		Activate Tilt Up / Down and Back Up, then try "Home" again.

Model Identification / Compliance Chart

Model	Description	Complies To:				Electrical Ratings:		
		UL 60601-1	CAN/CSA 22.2, #601.1-M90	IEC 60601-1	EN 60601-1-2 (EMC)	CE	VAC +/-10%	Amps
630-010	Midmark 630 HUMANFORM® Procedures Chair (Non-Programmable) w/ Receptacle	•	•	•	•	•	Configurable 115 12 230 6	50 / 60
630-011	Midmark 630 HUMANFORM® Procedures Chair (Non-Programmable) w/ Receptacle and Wireless Controls	•	•	•	•	•	Configurable 115 12 230 6	50 / 60
630-012	Midmark 630 HUMANFORM® Procedures Chair (Non-Programmable) w/ Rotational Base	•	•	•	•	•	Configurable 115 9 230 4.5	50 / 60
630-013	Midmark 630 HUMANFORM® Procedures Chair (Non-Programmable) w/ Rotational Base and Wireless Controls	•	•	•	•	•	Configurable 115 9 230 4.5	50 / 60
630-020	Midmark 630 HUMANFORM® Procedures Chair (Programmable) w/ Receptacle	•	•	•	•	•	Configurable 115 12 230 6	50 / 60
630-021	Midmark 630 HUMANFORM® Procedures Chair (Programmable) w/ Receptacle and Wireless Controls	•	•	•	•	•	Configurable 115 12 230 6	50 / 60
630-022	Midmark 630 HUMANFORM® Procedures Chair (Programmable) w/ Rotational Base	•	•	•	•	•	Configurable 115 9 230 4.5	50 / 60
630-023	Midmark 630 HUMANFORM® Procedures Chair (Programmable) w/ Rotational Base and Wireless Controls	•	•	•	•	•	Configurable 115 9 230 4.5	50 / 60


Fire Code Ratings: All upholstery complies with California Bureau of Home Furnishing [Technical Bulletin 117](#) and California Code of Regulations, Sect. 93120-93120.12, Title 17. Optional upholstery is available that complies with California Bureau of Home Furnishing [Technical Bulletin 133](#).

Specifications



Specifications Chart	
Patient Weight (maximum): 630 (w/standard base) 630 (w/rotational base)	294 kg (650 lbs) 294 kg (650 lbs)
Paper Roll (maximum size):	45.7 cm x 8.9 cm diameter (18 in. long x 3.5 in.)
Range of Motion & Dimensions:	<i>[See Range of Motion & Dimensions pages]</i>
Weight of Chair: 630 w/upholstery (standard base) w/packaging & skid (no uph.) 630 w/upholstery (rotational base) w/packaging & skid (no uph.) Uph. w/packaging (shipped separately)	215 kg (475 lbs) 240 kg (530 lbs) 249 kg (550 lbs) 274 kg (605 lbs) 20 kg (45 lbs)
Power Cord Length:	243.8 cm (8 ft.) long
Electrical Requirements:	<i>[See Model Identification / Compliance Charts]</i>
Foot / Hand Control Voltage:	12 VDC, SELV (Safety Extra Low Voltage)
Duplex Receptacles Maximum Load: (not available on rotation models)	115 VAC, 3A, 50/60 Hz
Fuse Ratings: <i>F1, F2, & F3 (on Motor Control Hub board):</i> <i>F1 & F2 (on Power Supply board):</i> <i>F3, F4, & F5 (on Power Supply board):</i>	5A, 250V, Type-T, 5 x 20 mm 6.3A, 250V, Type-T, 5 x 20 mm 10A, 250V, Type-T, 5 x 20 mm
Duty Cycle (Motor Run Time):	Intermittent Operation <i>[Domestic: 30 seconds ON / 5 minutes OFF]</i> <i>[Export: 30 seconds ON / 5 minutes OFF]</i>
(Receptacles):	Continuous Operation
Classifications:	Class I, Type B Applied Part Intermittent Operation <i>[Domestic: 30 seconds ON / 5 minutes OFF]</i> <i>[Export: 30 seconds ON / 5 minutes OFF]</i>
Protection against ingress of fluids:	IPXO <i>[Foot Control only: IPX1]</i>

Range of Motion

Models w/Rotational Base


Dimensions

MA8847i



Dimensions - (continued)

Standard Base Casting

MA99021

Dimensions - (continued)

Warranty Information

SCOPE OF WARRANTY Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will, at Midmark's option, repair or replace components of the domestic and international medical products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. The sole remedy under this limited warranty is the repair or replacement, at Midmark's option, of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which are determined to exist upon examination by Midmark. This warranty extends only to the original retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

Midmark warrants to the original retail purchaser that during the applicable warranty period it will repair or replace software contained within the products manufactured by Midmark (except for those not warranted under "Exclusions") if: (1) the media on which the software is furnished exhibits defects in material or workmanship under normal use; or (2) the software does not substantially conform to its published specifications.

APPLICABLE WARRANTY PERIOD The applicable warranty period, measured from the date of invoice to the original retail purchaser of the product and shall be one (1) year for all warranted products and components.

OBTAINING WARRANTY SERVICE Warranty service must be obtained through either Midmark or an authorized dealer in the Midmark product line for which warranty service is requested. Midmark may be contacted for warranty service inquiries or issues via email at www.midmark.com; by phone at 1-800-MIDMARK or by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380. It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

EXCLUSIONS: This limited warranty does not cover and Midmark shall not be liable for the following:

- (1) Defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) Products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environmental conditions and electrical requirements;
- (3) Products considered to be of a consumable or sterile nature;
- (4) Accessories or parts not manufactured by Midmark;
- (5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (6) Costs and expenses of routine maintenance and cleaning; and
- (7) Representations and warranties made by any person or entity other than Midmark.
- (8) Matching of color, grain or texture except to commercially acceptable standards;
- (9) Changes in color caused by natural or artificial light;
- (10) Custom manufactured products;
- (11) Alterations or modifications to the product by any person or entity other than Midmark; and
- (12) Products that would otherwise be covered under this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical market.

SOFTWARE; WITH RESPECT TO SOFTWARE THAT IS A PRODUCT OR COMPONENT THEREOF, MIDMARK DOES NOT WARRANT THAT THE SOFTWARE: (1) IS ERROR FREE; (2) CAN BE USED WITHOUT PROBLEMS OR INTERRUPTIONS; OR (3) IS FREE FROM VULNERABILITY TO INTRUSION OR ATTACK BY VIRUSES OR OTHER METHODS.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, LOSS OF DATA, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. THIS DISCLAIMER SHALL SURVIVE ANY FAILURE OR ASSERTED FAILURE OF THE ESSENTIAL PURPOSE OF THIS LIMITED WARRANTY OR ITS REMEDIES SPECIFIED HEREIN. **WARRANTY DISCLAIMER** THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

STATUTE OF LIMITATIONS No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

NO AUTHORIZATION No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

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