

BICEP RESTRAINT

IMPORTANT PRODUCT INFORMATION—PLEASE READ BEFORE USING

INTENDED USE: To restrict motion of the bicep. Use only as specified in the instructions.



Leather Bicep Cuff

INSTRUCTIONS: Open restraint completely and lay flat. Put the bicep which is to be restrained on the flattened restraint and fold closed and secure restraint by securing the cuff strap. Continue by wrapping the long strap

around an anchor and feeding through the buckle; fasten securely.
NOTE: Do not use for prolonged periods of time. Improper adjustment may prohibit proper circulation.

- INVENTORY ID:**
- Leather Bicep Cuffs **LBIC-100**
 - Polyurethane Bicep Cuffs **PBIC-100**
 - Polypropylene Bicep Cuffs **PPBC-200**

WEIGHT: 2 lbs., 2 lbs., 1 lb.

DESCRIPTION: Cuffs are 20" long with 3" wide leather or polyurethane or 2" wide polypropylene straps. Strap to secure to bed or rail is 48" long.

SAFETY CHECK: Including, but not limited to stitching, rivets, leather, spring buckle, and hook and loop.

MATERIALS USED: Spring buckles, roller buckle, polypropylene, leather or polyurethane, and hook & loop.

CLEANING AND MAINTAINING: Wipe leather or polyurethane with a mild soap and warm water. After wiping leather, allow it to air dry and apply neat's-foot oil. Dry clean or launder polypropylene at less than 280°F. When decontaminating bloodborne pathogens, use bleach while laundering.

STORAGE: Store away from heat greater than 280°F or open flame.

Caution: Federal law restricts this device to sale by or on the order of a physician.

Rev 8/25/16