

## Welcome to the Joylux family!

From one eternally busy woman to another, I would like to share a few quick pointers with you:

**To get the best results, use the device every other day.** We all get busy and set things aside, especially when we can't see visible results. Trust us—your device is doing exactly what it is supposed to do even though you can't see it. Stay the course, and you will be thrilled with the benefits.

**Speak out.** Talk to your girlfriends, your doctor, your partner. You deserve the best care, and that comes with open and honest dialogue about intimate wellness.

**Live your life to the fullest.** We are thrilled to partner with you in your journey to optimal wellness. Go ahead—laugh hard, exercise, reconnect with your partner. We have you covered.

Here's to living our best lives,



Colette Courtion, XO


Founder & CEO



# QUICK START GUIDE

Models VS1100, VS1120, VS1150, VS1170 (Session times and vibration modes vary by model.)


Before you begin using your device, please read the enclosed User Guide. For your convenience, this Quick Start Guide is provided as a reference after **all instructions and warnings have been read**. Please note: If you are pregnant or think you may be pregnant or are using photosensitive medication, do not use the device.



**1** Charge 12 hours prior to first use. The charging cable is located inside the cardboard insert in the box.




**2** Apply desired amount of **Photonic Gel** externally to the genital area or to the tip of the device prior to vaginal insertion.



**3** Turn on your device by pressing and holding the **power button** until the center green light appears.



**4** Press and release the **light button** to select your preferred session time, 6, 8, 10 or 12 minutes.



**5** Once your preferred time is selected, press and hold the same **light button AGAIN** until the lights and vibration turn on.

**6** Insert device, lie back, and relax—knowing you are improving your vaginal wellness. Please note: At any point during use, if the device becomes too warm, stop using it. At next use, reduce the session time to fewer minutes until your body becomes accustomed to the heat.

**7** After use, the light function of the device will lock for 2.5 hours. During this time, you may continue using the device with vibration only. There are up to 10 different vibration modes. **Press and release the vibration button** to scroll through the modes until you find your favorite. You can also **turn vibration off by pressing and holding the vibration button**.

**8** When you are finished using your device, **clean it thoroughly** with warm water and antibacterial soap and ensure all residue is removed. Let it air dry, then store in a cool, dry space.

**Congratulations on investing in yourself. You deserve it!**