

Instructions: One Pass Ostomy Draining Device (OPODD)

<ol style="list-style-type: none"> 1. Take OPODD: un-clip, and open rollers to a V-shape. 2. Place a container near pouch opening. 	
<ol style="list-style-type: none"> 3. Place OPODD across pouch (just underneath stoma opening): Note: One roller (with clip) should be on front of pouch and other roller (no clip) at the rear of pouch. 	
<ol style="list-style-type: none"> 4. Re-clip OPODD and Secure OPODD across the pouch. 5. Prepare Pouch opening: remove Pouch-clip or undo Pouch-fastener; to open draining end of pouch. 	
<ol style="list-style-type: none"> 6. Fold the sealing end of Pouch back exposing the opening. This is to allow draining. 	
<ol style="list-style-type: none"> 7. Pass (roll) OPODD over entire length of pouch stopping before the end of pouch. To avoid contact with fecal matter, do not pass OPODD over pouch end opening 	
<ol style="list-style-type: none"> 8. Un-clip OPODD and remove OPODD from pouch. 9. Wipe excess material from end of pouch. 	
<ol style="list-style-type: none"> 10. Fold end of pouch back to original position. 11. Apply Pouch-clip or close Pouch-fastener. 	
<p>Thank you.</p>	





CAUTION



- FOR EXTERNAL USE ONLY
- FOR USE on EXTERNAL Ostomy POUCH surface ONLY
- NO ASSEMBLY REQUIRED
- DO NOT DIS-ASSEMBLE DEVICE

WARNINGS



- DO NOT USE IF PACKAGE IS DAMAGED
- DO NOT USE IF DEVICE (OPODD) IS BROKEN
- Cleaning agents with chlorine or chloride as the active ingredient are corrosive to stainless steel and must not be used. Enzymatic and cleaning agents with neutral pH are recommended.



- Contains or presence of Natural Rubber Latex in End Caps

CLEANING INSTRUCTIONS

Scope: Recommended for first use, daily, and higher frequency of cleaning/disinfecting to be based on daily use of product at discretion of user.

Warnings:

- Do **not** use cleaning agents with chlorine or chloride as the active ingredient are corrosive to stainless steel and must not be used.
- Do **not** use abrasive scouring powders or metal scouring pads or wire brushes.
- Enzymatic and cleaning agents with neutral pH are recommended
- Enzyme solutions should be used in accordance with manufacturer's instructions
- Disinfectants should be used in accordance with manufacturer's instructions
- A Manual cleaning process is recommended

Steps:

1. Cleaning: Clean surfaces using water and mild cleaning agent (following applicable Manufacturer's cleaning directions for that disinfectant). Required before disinfection. A soft brush may be used.
2. Disinfecting: Wipe all surfaces with a disinfectant or disinfectant wipes (following applicable Manufacturer's cleaning directions for that disinfectant).
3. In particular, focus on rollers and any area of Device that may have had contact with intact skin, or Ostomy pouch opening.
4. Allow to air dry after disinfection.
5. At a minimum the Device should be visibly clean.

