

Programming Instructions

- **Set CPAP Mode Pressure**

- Unlock Z2 – **press and hold start/stop button** until a number and *CPAP-LOC* appears.
- Press and hold **up and down arrow buttons** until a flashing number appears and *CPAP* appears.
- Use the **up and down arrow buttons** to set desired *CPAP* pressure.
- Once set, **press and hold start/stop button** until *ACCEPT* appears on screen. This will return you to home screen.

- **Set Ramp Time**

- While in standby mode, **press and hold the up arrow button** for 3 seconds until flashing *00 MINUTES* appears.
- While flashing, use **up and down arrow buttons** to increase or decrease the ramp time in 5-minute increments.
- When desired ramp time is reached, **press and hold start/stop button** until *ACCEPT* appears on screen. This will return you to home screen.
 - Ramp Time must be set to a minimum of 5 minutes prior to attempting to program Ramp Start Pressure or change to Ramp Start Pressure will not save, and you will receive *FAIL 2*.
 - To disable Ramp Mode, set ramp time to *00 MINTUES*.

- **Set Ramp Start Pressure**

- **Press and hold the down arrow button** for 3 seconds until ramp start pressure value flashes.
- While flashing, use **up and down arrow buttons** to increase or decrease the ramp start pressure by 0.5 cm H₂O increment. (Ramp start pressure must be lower than prescribed pressure.)
- When desired ramp start pressure is reached, **press and hold start/stop button** until *ACCEPT* appears on screen. This will return you to home screen.

Programming Instructions

- **Set Mode - CPAP or APAP and Z-Breathe® setting**

- Press and hold the **up and down arrow buttons** at the same time.
- Flashing *APAP-1, 2, 3* or flashing *CPAP-1, 2,3* will appear. Release buttons.
- Use the **up and down arrow buttons** to select desired mode and Z-Breathe® setting.
 - 1 – offers mild pressure swing relief
 - 2 – offers moderate pressure swing relief (this is the default setting)
 - 3 – most aggressive setting for pressure swing relief
- Press and hold the **start/stop button** to save selection – hold until *ACCEPT* appears on screen. This will return you to the home screen.

- **Set APAP Mode Pressures – Minimum and Maximum**

- Unlock Z2 – **press and hold start/stop button** until a number and *APAP-LOC* appears.
- Press and hold **up and down arrow buttons** until a flashing number appears and *APAP HI* appears.
- Use the **up and down arrow buttons** to set desired *APAP HI* pressure.
- Once set, press and hold **start/stop button** until *ACCEPT* appears on screen.
- A flashing number and *APAP LO* appears on screen. Set desired *APAP LO* pressure.
- Once set, press and hold **start/stop button** until *ACCEPT* appears on screen. This will return you to home screen.

Programming Instructions

- **Set Auto Start / Stop**

- **Press and hold the up arrow button and down arrow button** at the same time until the Z-Breathe setting flashes. Release the buttons.
- **Press and hold the up arrow button and down arrow button** at the same time again until flashing *ROFF-ON* appears on screen. Release the buttons.
- **Press down arrow button** to change the auto start/stop to off. *ROFF-OFF* will appear on the screen.
- **Press and hold the start/stop button** until *ACCEPT* appears on screen. This will return you to home screen.
- To change the auto start/stop back to on. Follow above steps, until flashing *ROFF-OFF* appears on screen. Release the buttons.
- **Press up arrow button** to change the auto start/stop to on. *ROFF-ON* will appear on the screen.
- **Press and hold the start/stop button** until *ACCEPT* appears on screen. This will return you to home screen.

- **Turn on Bluetooth (this will allow connection to Z2 Nitelog® Mobile App)**

- Immediately after plugging the power adapter into Z2, **hold the up arrow button** until *bt-ON* appears. This indicates your Bluetooth has been turned on. You will be returned to the home screen automatically.
- To turn the Bluetooth off, immediately after plugging the power adapter into Z2, *hold the down arrow button* until *bt-OFF* appears. This indicates your Bluetooth has been turned on. You will be returned to the home screen automatically.