

GHOST
STRYKE
FOR EXTREME HUNTING & TACTICAL ADVANTAGE™



*Ghost Stryke Sound Enhancement & Hearing
Protection Manual & Warranty*

*Sport***EAR**

GHOST
STRYKE

FOR EXTREME HUNTING & TACTICAL ADVANTAGE™



!!! Important Performance & Safety Tips for Your New Ghost Stryke

Congratulations on the purchase of your new Ghost Stryke personal sound amplification devices. Whether you are using this for general sound enhancement or hearing protection or both, Ghost Stryke offers you the best product for fit, performance and reliability.

To make sure you get the best performance from your new Ghost Stryke devices, please read the manual entirely, so you know how to best use your new devices.

Hearing loss happens gradually or quickly, depending on the levels of noise(s) exposed to your ears. Once hearing loss occurs, it does not reverse its course, THE DAMAGE IS DONE. Many who suffer with hearing loss wish they would have understood this better and sooner. Take action now and protect your ears from severe damage.


ALWAYS make sure you use the product as intended for best results. ALWAYS use your new Ghost Stryke devices when in damaging loud environments.

Make sure your product is working properly before use in extreme environments, so you get the performance and protection your ears need.

See the following charts on "How the Ear Works" and the decibel (dB) rating chart


How the Ear Works

1. Acoustic design of the outer ear captures the oncoming sound waves.
2. Sound waves enter the ear canal and penetrate the eardrum (tympanic membrane).
3. The eardrum taps the malleus causing a chain reaction with the stapes and the incus.
4. The incus taps the cochlea, causing a hydraulic effect on 35,000-40,000 hair-like nerves (cilia) in the cochlea.
5. As the cilia are pushed, they send a signal through the auditory nerve to the brain.
6. The brain registers the incoming signal as a specific sound.



How the Ear Gets Damaged

1. As the nerves continually bend over from oncoming sound, they eventually begin to wear out or even break off (especially from damaging loud sounds).
2. Therefore, they do not respond as sensitively as they used to and your ability to hear continues to diminish over time. Extreme loud noises can damage your hearing instantly and permanently.



How to Help the Ear Hear Better

A traditional piano has 88 strings that allow you to create many beautiful sounds for your enjoyment. A traditional guitar only has six strings. Your ear has 35,000-40,000 "strings on a piano". This is the most advanced sound system on the planet. But just as strings on a piano, they can be easily damaged if you are not careful. Once damaged, **they cannot be restored to normal**. If your ears have been damaged and you need to hear better, then modern technology offers digital hearing enhancement. This allows you to have personal hearing device(s) tuned to match the deficiency in your hearing. If you are missing more high tones in your hearing, then you can have your devices programmed to enhance those sounds only. Therefore, your hearing device(s) will give the added stimulation to those damaged nerves (cilia) without affecting the other sounds. This will help you hear things again more clearly. Please see the entire SportEAR product line and call one of our SportEAR product specialists to help determine any additional needs.

The Decibal Rating Chart

SOUND	NOISE LEVEL (dB)	EFFECT
Whisper	30	Very Quiet
General Household Electronics	70-90	Moderate to Loud & Invasive
Industrial Environments	85-140	For any noise over 85 db see OSHA note below
Music Concerts	110-140	See OSHA & Threshold of pain notes below
Power Tools	110-140	
Firearm(s) Shooting	110-185	Extreme Hearing Damage

OSHA Hearing Chart shows hearing damage and loss occur when exposed to sounds at 85 dB or louder for an 8 hour or greater period of time.

Threshold of pain occurs for most people around 125 dB.

General hearing enhancement & protection has become a widely discussed topic in recent years. More than ever, people and businesses realize the dangers of loud noise exposure on the human ear. Anyone suffering from hearing loss would reverse it immediately, if possible. Now with the Ghost Stryke, **you can protect your hearing while maintaining your ability to hear well**, even enhance your hearing above its normal capacity.

Once permanent hearing loss sets in, speech and communication become very difficult. Even with hearing aids, your hearing will never be the same. Whether you use your Ghost Stryke for work or recreation, please remember to use them anytime you are exposing your ears to damaging loud sounds. Ghost Stryke offers an incredible solution to help you work and play safely, by hearing important noises around you and protecting your hearing at the same time.

***Always Protect Your Hearing
in Every Noisy Environment, Maintain
Your Hearing Throughout Your Life!***


Ear Tip(s) Selection

Enhancement with Automatic Hearing Protection


For maximum hearing protection, choose the foam tip that best fits your ear size. Make sure that you have a complete seal of your ear before exposing your ears to loud sounds. This can be done by inserting the Ghost Stryke devices and then clapping next to your ear or presenting another semi-loud sound to your hearing, so you make sure you have them inserted correctly. **A leak in the seal from improper insertion can still lead to hearing damage.**



24 dB NRR




30 dB NRR



Always make sure you have pushed the tip completely on, so it will not come off in your ear canal. If it ever does come off deep in your ear canal, simply contact a hearing clinic or your physician to have it removed safely.

Battery Installation

1. Only use size 10 hearing aid batteries.
2. Remove a battery from the package.
3. Remove the yellow tab from the battery.
4. Gently open the battery door.
5. Place end of battery that has black circle around the edge into the battery, so the black circle lines up with the round lip in the battery door.
6. Gently close the battery door completely.



Ghost Stryke Product Insertion

1. Make sure you insert the correct product in each ear (Left and Right specific; watch for L & R markings).
2. Once you have the desired tip on your device, place the end with the tip on it in your ear. The smaller circle on the design should sit just outside your ear canal. The larger circle should rest comfortably in the lower bowl (concha) of your outer ear. The larger circle may also rotate upwards for smaller ears.




SecurFit™ Extenders

1. Simply place the SecurFit™ Extender on your device as shown.
2. Carefully place the device in your ear, as described above, and gently push the upper portion of the SecurFit™ Extender in the upper portion of your ear as shown. Make sure this part is not sticking out of your ear, as you will not have a secure fit for best comfort and performance.


SecurFit(TM) Extenders aren't necessary for all ears.

They are helpful to secure the device for certain ear shapes and sizes or when using your device(s) for extra rigorous activities.



Adjusting the Volume

- Your new Ghost Styke devices use a volume wheel to increase or decrease the volume. Once you have the devices in place in your ear, to turn the volume up, rotate the volume wheel towards your face. To reduce the volume, rotate the volume wheel away from your face.
- Choose the volume level that is most comfortable to you and best for your intended application.



Sound Check


- To make sure the device is working, place device in the palm of your hand and cup your hand around the device. If you can hear it whistling in the palm of your hand, it means the battery and the device are both working normally.
- To check sound quality, simply insert it in your ear and listen to yourself or some one else talking. Try the telephone or listen to the TV. If you have clear and enhanced sound, then it is working. If you have a muffled sound or plugged feeling, it may not be working correctly. If it's not working, then see the Troubleshooting section.

Hearing Protection Check

- To insure you are getting the hearing protection you need, insert the device(s) in your ear(s), and then clap your hands next to your ear(s). If the clap sounds sharp and clear still, you do not have it sealed properly for hearing protection. If the sound is flat and dull, then you have an accurate seal for hearing protection.

Cleaning

- For best results, use the Ghost Stryke Cleaning Brush provided. A dedicated medium to firm bristled toothbrush also works very well. Hold the device upside down (sound port pointing towards the ground) and brush firmly in the sound port to clear out any earwax or debris. For best results, allow earwax to dry prior to cleaning.
- Use a tissue to wipe down the outer shell periodically or dab a bit of rubbing alcohol on a tissue for a more thorough cleaning of the shell itself or the tips.



Care & Storage

- Always keep your Ghost Stryke clean for safe, hygienic use. Never wear if they appear dirty or broken.
- When not in use, always store your Ghost Stryke in the protective case provided.
- Keep case closed. Store case and batteries in dry, clean environment.
- Do not store other items in the case to avoid breakage or contamination of product.
- Avoid exposure to liquids and high humidity.
- Keep product and batteries out of the reach of children and pets.
- **DO NOT STORE BATTERIES OR PRODUCT IN A FREEZER OR REFRIGERATOR OR IN EXTREME HEAT. STORE IN A NORMAL ROOM TEMPERATURE ONLY FOR BEST RESULTS.**

Troubleshooting

- If your device is not working, make sure you clean out the sound port thoroughly, replace the battery with a battery you know to be good, clean anything else off the device that may be causing the problem and then hold it in your hand to see if it whistles. If problem still occurs call our customer service for additional help.

Customer Service

- At SportEAR, we appreciate your business and are here to help with all your questions. Please feel free to reach us Monday - Friday Mountain Standard Time (MST) from 9 am – 5 pm at 1-866-422-5502

LIMITED 1-YEAR CONSUMER WARRANTY

SportEAR, LLC, 12397 S. 300 E. Suite 100, Draper, UT 84020 USA
1-866-422-5502

Please read the following before submitting a claim. Also read below "How to Submit a Claim." SportEAR, LLC extends this warranty to You. Statutory or common law may provide You with additional rights or remedies, which shall not be affected by this Limited Warranty.

DEFINITIONS

"Adequate Use" means use of the Product (i) within home dwelling, (ii) for private (as opposed to commercial) purposes, (iii) in conformance with all applicable local, state and federal laws, codes or regulations (including without limitation building and/or electrical codes). (iv) in accordance with manufacturer recommendations and/or instructions in the materials and documentation that accompany the Product, and (v) if applicable, with proper electrical grounding.

“Authorized Dealer” means any distributor, reseller or retailer that (i) was duly authorized to do business in the jurisdiction where it sold the Product to You, (ii) was permitted to sell You the Product under the laws of the jurisdiction where You bought the Product, and (iii) sold You the product new and in its original packaging.

“Formal Warranty Claim” means a claim made in accordance with the section “Formal Warranty Claims” herein.

“Product” means a Product (i) that is an official SportEAR Product that You purchased from an Authorized Dealer new and in its original packaging, and (iii) whose serial number, if any, has not been removed, altered, or defaced.

“Product Defect” means an inadequacy of the Product that existed at the time when You purchased it from an Authorized Dealer and causes a failure of the Product to perform in accordance with SportEAR’s documentation accompanying the Product, unless such failure has been caused completely or partly by (a) any use other than Adequate Use, (b) transportation, neglect, misuse or abuse by anyone other than SportEAR employees, (c) alteration, tampering or modification of the product by anyone other than a SportEAR employee, (d) accident (other than a malfunction that would otherwise qualify as a Product Defect), (e) maintenance or service of the Product by anyone other than an SportEAR employee, (f) exposure of the Product to heat, bright light, sun, liquids, sand or other contaminants, or (g) acts outside the control of SportEAR, including without limitation acts of God, fire, storms, earthquake or flood.

“Warranty Period” means the time period during which SportEAR must have received your Formal Warranty Claim. SportEAR’s entire Ghost Strike line carries a One Year Limited Warranty. The One Year Period commences on the date of original purchase of Your product or when you actually received Your new product from the Authorized Dealer (whichever comes later) as evidenced by Your original receipt or the Authorized Dealer’s Invoice, sales receipt or packing slip. If You do not have written proof of the date of purchase or receipt, then the Warranty Period commences three (3) months after the date the product left SportEAR as evidenced by SportEAR’s shipping records.

PLEASE NOTE: YOU MUST CALL SPORTEAR AND OBTAIN A RETURN AUTHORIZATION NUMBER (as described in the “How to Make a Claim” section) within two (2) months after You discover a Product Defect (or should have discovered it, if Product Defect was obvious).

“You” means the first individual person that purchased the Product in its original packaging from an Authorized Dealer. This Limited Warranty does not apply to persons or entities that bought the Product (i) in used or unpackaged form, (ii) for resale, lease or other commercial use, or (iii) from someone other than an Authorized Dealer.

SCOPE OF THIS LIMITED WARRANTY

PRODUCTS If the Product contained a Product Defect when You bought it from an Authorized Dealer and SportEAR receives a Formal Warranty Claim from You within two (2) months after You discover such Product Defect (or should have discovered it, if such Product Defect was obvious) and before the end of the Warranty Period for Product Defects applicable to the affected Product, then SportEAR will provide You with one of the following remedies: (1) SportEAR will repair or, at SportEAR's sole discretion, replace the Product, or (2) refund to You the purchase price You paid to the Authorized Dealer for the affected Product if repair or replacement is not commercially practicable or cannot be timely made. NOTE: SPORTEAR DOES NOT ASSUME ANY LIABILITY FOR ANY INCIDENTAL, CONSEQUENTIAL OR INDIRECT DAMAGES UNDER THIS LIMITED WARRANTY.

GENERAL PROVISIONS

CHOICE OF LAW/JURISDICTION This Limited Warranty and any disputes arising out of or in connection with this Limited Warranty ("Disputes") shall be governed by the laws of the state of Utah, USA, excluding conflicts of law principles and excluding the Convention for the International Sale of Goods. The courts of Utah shall have exclusive jurisdiction over any Disputes.

OTHER RIGHTS

THIS LIMITED WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY ALSO HAVE OTHER RIGHTS, WHICH VARY FROM STATE TO STATE AND JURISDICTION TO JURISDICTION, AND WHICH SHALL NOT BE AFFECTED BY THIS LIMITED WARRANTY. THIS WARRANTY EXTENDS ONLY TO YOU AND CANNOT BE TRANSFERRED OR ASSIGNED. If any provision of this Limited Warranty is unlawful, void or unenforceable, that provision shall be deemed severed and shall not affect any remaining provisions. In case of any inconsistency between the English and other versions of this Limited Warranty, the English version shall prevail.

FORMAL WARRANTY CLAIM

HOW TO MAKE A CLAIM In the event of a product defect, You must follow these instructions: (1) Call SportEAR within two (2) months after You discover a Product Defect (or should have discovered it, if such Product Defect was obvious); (2) Give a detailed explanation of the product defect; (3) Obtain a Return Authorization Number; (4) Upon receipt of a claim form (which may be sent to You after You filed Your Formal Warranty Claim), fill out the claim form entirely; (5) Return the Products, shipping prepaid by You (to be refunded if you are entitled to a remedy under the Scope of this Limited Warranty), to SportEAR for verification of product defect, along with a copy of Your original sales receipts and proof

of purchase (UPC label or packaging slip) for such Products, the completed claim form, and printed Return Authorization Number on the outside of the return package (the claim form will include instructions for return).

CONTACT INFORMATION

By phone in the United States and Canada: 1-866-422-5502

By email: customersupport@sportear.com

For additional phone numbers please visit SportEAR.com for details.

FURTHER PROCEEDINGS SportEAR will determine whether a Product Defect existed. SportEAR, at its sole discretion, may direct You to obtain a repair estimate at a service center. If a repair estimate is required, You will be instructed on how to properly submit the estimate and the resulting invoice to SportEAR for payment. Any fees for repairs may be negotiated by SportEAR.

TIMING If You bring a Formal Warranty Claim and fully comply with all terms and conditions of this Limited Warranty, SportEAR will use its best efforts to provide You with a remedy within thirty (30) days after receipt of Your Formal Warranty Claim (if you reside in the United States or forty-five (45) days if You reside elsewhere), unless obstacles outside SportEAR's control delay the process.

GHOST
STRYKE

FOR EXTREME HUNTING & TACTICAL ADVANTAGE™



GHOST
STRYKE

FOR EXTREME HUNTING & TACTICAL ADVANTAGE™

