



NATtrol™ *Chlamydia trachomatis* External Run Controls
Strain: LGVII 434

Catalog #: NATCT(434)-ERCL

Catalog #: NATCT(434)-ERCM

PRODUCT DESCRIPTION:

NATtrol™ *Chlamydia trachomatis* External Run Controls (NATCT(434)-ERCL and NATCT(434)-ERCM) are formulated with purified, intact bacterial cells that have been chemically modified to render them non-infectious and refrigerator stable*. Each control pack contains 6 x 1.0 mL vials of *C. trachomatis* NATtrol™ at target concentrations listed in Table 1. These controls are supplied in a purified protein matrix that mimics the composition of a true clinical specimen.

*NATtrol™ Patents Pending

INTENDED USE:

- NATtrol™ *Chlamydia trachomatis* External Run Controls are designed to evaluate the performance of nucleic acid tests for determination of the presence of *C. trachomatis* DNA. NATCT(434)-ERCL and NATCT(434)-ERCM can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel.
- NATCT(434)-ERCL and NATCT(434)-ERCM contain intact organisms and should be run in a manner identical to that used for clinical specimens.

ETIOLOGIC STATUS/BIOHAZARD TESTING:

- NATtrol™ inactivation was carried out on the *C. trachomatis* stock used to formulate each control pack. The inactivation was verified by the absence of Chlamydial inclusion forming units (IFU).
- The purified protein matrix was manufactured from materials that were screened and found to be negative for HIV 1&2 Ab, HBsAg, HTLV I&II Ab, HCV Ab, HIV RNA, HBV DNA and HCV RNA using FDA cleared kits at the single donor level.

This product was manufactured in a facility which has a Quality Management System that is certified as being in compliance with ISO 13485.

PINATCT(434)-ERCM
PINATCT(434)-ERCL
Rev.: 12/2012
For Research Use Only

PRECAUTIONS:

- Although NATCT(434)-ERCL and NATCT(434)-ERCM contain inactivated organisms, they should be handled as if potentially infectious.
- Use Universal Precautions when handling this product.
- To avoid cross-contamination, use separate pipette tips for all reagents.

RECOMMENDED STORAGE:

- NATtrol™ *Chlamydia trachomatis* External Run Controls should be stored at 2-8°C.

INSTRUCTIONS FOR USE:

- Extract *C. trachomatis* DNA prior to use in downstream assays.

TABLE 1

Catalog Number	Target Concentration (copies/mL)*
NATCT(434)-ERCL	5,000
NATCT(434)-ERCM	50,000

*Quantitation based on in-house real time PCR assay targeting 16s RNA gene region.

DO NOT USE IN HUMANS

These products are intended for research, product development, quality assurance or manufacturing use. These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

ZeptoMetrix Corporation
878 Main Street
Buffalo, New York 14202
Phone: 716-882-0920
Fax: 716-882-0959
Customer Service: 800-274-5487
www.zeptometrix.com