



NATtrol™ EBV External Run Controls
Strain: B95-8

Catalog #: NATEBV-ERCL

Catalog #: NATEBV-ERCM

PRODUCT DESCRIPTION:

NATtrol™ EBV External Run Controls (NATEBV-ERCL and NATEBV-ERCM) are formulated with purified, intact virus particles that have been chemically modified to render them non-infectious and refrigerator stable*. Each control pack contains 6 x 1.0 mL vials of EBV NATtrol™ at concentrations listed in Table 1. These controls are supplied in a purified protein matrix that mimics the composition of a true clinical specimen.

*NATtrol™ Patents Pending

INTENDED USE:

- NATtrol™ EBV External Run Controls are designed to evaluate the performance of nucleic acid tests for determination of the presence of EBV DNA. NATEBV-ERCL and NATEBV-ERCM can also be used for validation of clinical assays, development of diagnostic tests and training of laboratory personnel.
- NATEBV-ERCL and NATEBV-ERCM contain intact virus and should be run in a manner identical to that used for clinical specimens.

ETIOLOGIC STATUS/BIOHAZARD TESTING:

- NATtrol™ inactivation was carried out on each member in the panel. The inactivation was verified by the absence of viral growth in validated tissue culture based infectivity assays.
- The purified protein matrix was manufactured from materials that were screened and found to be negative for HIV 1&2 Ab, HBsAg, HTLV I&II Ab, HCV Ab, HIV RNA, HBV DNA and HCV RNA using FDA cleared kits at the single donor level.

PRECAUTIONS:

- Although NATEBV-ERCL and NATEBV-ERCM contain inactivated virus, they should be handled as if potentially infectious.
- Use Universal Precautions when handling this product.
- To avoid cross-contamination, use separate pipette tips for all reagents.

RECOMMENDED STORAGE:

- NATtrol™ EBV External Run Controls should be stored at 2-8°C.

INSTRUCTIONS FOR USE:

- Extract EBV DNA prior to use in downstream assays.

TABLE 1:

Catalog Number	Concentration (copies/mL)
NATEBV-ERCL	5,000
NATEBV-ERCM	50,000

*Conversion factor: 1.00 IU = 9.03 copies. Based on internal testing of the 1st WHO International Standard for Epstein-Barr Virus for Nucleic acid Amplification Techniques (NIBSC code: 09/260)

DO NOT USE IN HUMANS

These products are intended for research, product development, quality assurance or manufacturing use. These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

This product was manufactured in a facility which has a Quality Management System that is certified as being in compliance with ISO 13485.

PINATEBV-ERCM
PINATEBV-ERCL
Rev.: 04/2013
For Research Use Only

ZeptoMetrix Corporation
878 Main Street
Buffalo, New York 14202
Phone: 716-882-0920
Fax: 716-882-0959
Customer Service: 800-274-5487
www.zeptometrix.com