

## NATtrol™ Influenza/RSV External Run Controls

Catalog Number: NATFLUAB-6C

Catalog Number: NATRSV-6C

Catalog Number: NATCXVA9-6C

### PRODUCT DESCRIPTION:

**NATtrol™ Influenza/RSV External Run Controls (NATFLUAB-6C, NATRSV-6C and NATCXVA9-6C)\*** are formulated with purified, intact virus particles that have been chemically modified to render them non-infectious and refrigerantly stable. Each control contains 6 x 0.5 mL vials of NATtrol™ Influenza A/B, NATtrol™ Respiratory Syncytial Virus, or NATtrol™ Coxsackievirus A9 (Negative Control). These controls are supplied in a purified protein matrix that mimics the composition of a true clinical specimen.

\*Pat.: <http://www.zeptometrix.com/patent-information/>

### INTENDED USE:

- NATtrol™ Influenza/RSV External Run Controls are full process controls designed to evaluate the performance of nucleic acid tests for determination of the presence of respiratory virus nucleic acids. NATFLUAB-6C, NATRSV-6C and NATCXVA9-6C can also be used for quality control of clinical assays and training of laboratory personnel.
- NATFLUAB-6C, NATRSV-6C and NATCXVA9-6C contain intact organisms and should be run in a manner identical to that used for clinical specimens.

### ETIOLOGIC STATUS/BIOHAZARD TESTING:

- NATtrol™ inactivation was carried out on the virus stocks used to formulate each control. The inactivation was verified by the absence of viral growth in validated tissue culture based infectivity assays.
- Purified protein matrix used in the manufacture of this product is treated with 0.09% sodium azide. It was manufactured from materials that have been tested and found non-reactive at the donor level for HIV-1/HIV-2 Antibody, HBsAg and HCV Antibody by FDA licensed donor screening test methods. All materials are also tested for HIV-1 and HCV by FDA approved Nucleic Acid Test (NAT) methods. Heat inactivated bovine based source materials used in the manufacture of this product meet applicable USDA requirements for abattoir sourced animals, traceability and country of origin. The materials were collected at USDA licensed establishments or legally imported from countries recognized by the USDA as negligible or controlled for risk for Bovine Spongiform Encephalopathy (BSE) and other exotic disease agents. Donor animals were inspected ante and post mortem at the abattoir as required by the USDA.

### PRECAUTIONS:

- Although NATtrol Influenza/RSV External Run Controls contain inactivated organisms, they should be handled as if potentially infectious.
- Use Universal Precautions when handling these products.
- To avoid cross-contamination, use separate pipette tips for all reagents.

### RECOMMENDED STORAGE:

- NATtrol™ Influenza/RSV External Run Controls should be stored at 2-8°C.

### INSTRUCTIONS FOR USE WITH Xpert® Flu/RSV XC and Xpert® Xpress Flu/RSV ASSAYS:

- Vortex NATtrol™ sample for 5-10 seconds.
- Open the cartridge lid. Using a clean 300 µL transfer pipette, transfer 300 µL (one draw) of the NATtrol™ sample to the sample chamber with large opening in the cartridge.
- Close cartridge lid and follow manufacturer's instructions.

**For all other assays follow manufacturer's Instructions.**




### EXPECTED RESULTS:

Catalog Number	Organism	Xpert® Flu/RSV XC Expected Result
NATFLUAB-6C	Influenza A/Brisbane/59/07 Influenza B/Florida/02/06	Flu A POSITIVE Flu B POSITIVE RSV NEGATIVE
NATRSV-6C	RSV B (CH93(18)-18)	Flu A NEGATIVE Flu B NEGATIVE RSV POSITIVE
NATCXVA9-6C	Coxsackievirus A9	Flu A NEGATIVE Flu B NEGATIVE RSV NEGATIVE

**DO NOT USE IN HUMANS. FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.**

These products are intended for research, product development, quality assurance or manufacturing use. These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

This product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.

REF	Catalog Number		Temperature Limitation
LOT	Lot Number		Expiration Date
	Biological Risk		

**ZeptoMetrix Corporation**  
878 Main Street,  
Buffalo, NY 14202  
800-274-5487  
[www.zeptometrix.com](http://www.zeptometrix.com)