

- 5.1 For levels of copper below 50 $\mu\text{g/L}$, either the Special Extraction Procedure, given in part 9.2 of the Atomic Absorption Methods section or the furnace technique, Method 220.2, is recommended.
- 5.2 Numerous absorption lines are available for the determination of copper. By selecting a suitable absorption wavelength, copper samples may be analyzed over a very wide range of concentration. The following lines may be used:
 327.4 nm Relative Sensitivity 2
 216.5 nm Relative Sensitivity 7
 222.5 nm Relative Sensitivity 20
- 5.3 Data to be entered into STORET must be reported as $\mu\text{g/L}$.
- 5.4 The 2,9-dimethyl-1,10-phenanthroline colorimetric method may also be used (Standard Methods, 14th Edition, p. 196).

6.0 Precision and Accuracy

- 6.1 An interlaboratory study on trace metal analyses by atomic absorption was conducted by the Quality Assurance and Laboratory Evaluation Branch of EMSL. Six synthetic concentrates containing varying levels of aluminum, cadmium, chromium, copper, iron, manganese, lead and zinc were added to natural water samples. The statistical results for copper were as follows:

Number of Labs	True Values $\mu\text{g/Liter}$	Mean Value $\mu\text{g/Liter}$	Standard Deviation $\mu\text{g/Liter}$	Accuracy as % Bias
91	302	305	56	0.9
92	332	324	56	-2.4
86	60	64	23	7.0
84	75	76	22	1.3
66	7.5	9.7	6.1	29.7