

**METHOD #:** 245.5

(Issued 1974)

**TITLE:**

Mercury In Sediment (Manual Cold Vapor Technique)

**ANALYTE:**

CAS # Hg Mercury 7439-97-6

**INSTRUMENTATION:**

AA

## 1.0 Scope and Application

- 1.1 This procedure<sup>(1)</sup> measures total mercury (organic & inorganic) in soils, sediments, bottom deposits and sludge type materials.
- 1.2 The range of the method is 0.2 to 5  $\mu\text{g/g}$ . The range may be extended above or below the normal range by increasing or decreasing sample size or through instrument and recorder control.

## 2.0 Summary of Method

- 2.1 A weighed portion of the sample is digested in aqua regia for 2 minutes at 95°C, followed by oxidation with potassium permanganate. Mercury in the digested sample is then measured by the conventional cold vapor technique.
- 2.2 An alternate digestion<sup>(2)</sup> involving the use of an autoclave is described in (8.2).

## 3.0 Sample Handling and Preservation

- 3.1 Because of the extreme sensitivity of the analytical procedure and the omnipresence of mercury, care must be taken to avoid extraneous contamination. Sampling devices and sample containers should be ascertained to be free of mercury; the sample should not be exposed to any condition in the laboratory that may result in contact or air-borne mercury contamination.
- 3.2 While the sample may be analyzed without drying, it has been found to be more convenient to analyze a dry sample. Moisture may be driven off in a drying oven at a temperature of 60°C. No mercury losses have been observed by using this drying step. The dry sample should be pulverized and thoroughly mixed before the aliquot is weighed.

## 4.0 Interferences

- 4.1 The same types of interferences that may occur in water samples are also possible with sediments, i.e., sulfides, high copper, high chlorides, etc.
- 4.2 Volatile materials which absorb at 253.7 nm will cause a positive interference. In order to remove any interfering volatile materials, the dead air space in the BOD bottle should be purged before the addition of stannous sulfate.

## 5.0 Apparatus

- 5.1 Atomic Absorption Spectrophotometer (See Note 1): Any atomic absorption unit having an open sample presentation area in which to mount the

absorption cell is suitable. Instrument settings recommended by the particular manufacturer should be followed.

NOTE 1: Instruments designed specifically for the measurement of mercury using the cold vapor technique are commercially available and may be substituted for the atomic absorption spectrophotometer.

- 5.2 Mercury Hollow Cathode Lamp: Westinghouse WL-22847, argon filled, or equivalent.
- 5.3 Recorder: Any multi-range variable speed recorder that is compatible with the UV detection system is suitable.
- 5.4 Absorption Cell: Standard spectrophotometer cells 10 cm long, having quartz end windows may be used. Suitable cells may be constructed from plexiglass tubing, 1 " O.D. X 4-1/2". The ends are ground perpendicular to the longitudinal axis and quartz windows (1" diameter X 1/16" thickness) are cemented in place. Gas inlet and outlet ports (also of plexiglass but 1/4" O.D.) are attached approximately 1/2" from each end. The cell is strapped to a burner for support and aligned in the light beam to give the maximum transmittance.

NOTE 2: Two 2" X 2" cards with one inch diameter holes may be placed over each end of the cell to assist in positioning the cell for maximum transmittance.

- 5.5 Air Pump: Any peristaltic pump capable of delivering 1 liter of air per minute may be used. A Masterflex pump with electronic speed control has been found to be satisfactory. (Regulated compressed air can be used in an open one-pass system.)
  - 5.6 Flowmeter: Capable of measuring an air flow of 1 liter per minute.
  - 5.7 Aeration Tubing: Tygon tubing is used for passage of the mercury vapor from the sample bottle to the absorption cell and return. Straight glass tubing terminating in a coarse porous frit is used for sparging air into the sample.
  - 5.8 Drying Tube: 6" X 3/4" diameter tube containing 20 g of magnesium perchlorate (See Note 3). The apparatus is assembled as shown in the accompanying diagram.
- NOTE 3: In place of the magnesium perchlorate drying tube, a small reading lamp with 60W bulb may be used to prevent condensation of moisture inside the cell. The lamp is positioned to shine on the absorption cell maintaining the air temperature in the cell about 10°C above ambient.

## 6.0 Reagents

- 6.1 Aqua Regia: Prepare immediately before use by carefully adding three volumes of conc. HCl to one volume of conc. HNO<sub>3</sub>.
  - 6.2 Sulfuric Acid, 0.5 N: Dilute 14.0 mL of conc. sulfuric acid to 1 liter.
  - 6.3 Stannous Sulfate: Add 25 g stannous sulfate to 250 mL of 0.5 N sulfuric acid (6.2). This mixture is a suspension and should be stirred continuously during use.
  - 6.4 Sodium Chloride-Hydroxylamine Sulfate Solution: Dissolve 12 g of sodium chloride and 12 g of hydroxylamine sulfate in distilled water and dilute to 100 mL.
- NOTE 4: A 10% solution of stannous chloride may be substituted for (6.3) and hydroxylamine hydrochloride may be used in place of hydroxylamine sulfate in (6.4).

- 6.5 Potassium Permanganate: 5% solution, w/v. Dissolve 5 g of potassium permanganate in 100 mL of distilled water.
- 6.6 Stock Mercury Solution: Dissolve 0.1354 g of mercuric chloride in 75 mL of distilled water. Add 10 mL of conc. nitric acid and adjust the volume to 100.0 mL. 1.0 mL = 1.0 mg Hg.
- 6.7 Working Mercury Solution: Make successive dilutions of the stock mercury solution (6.6) to obtain a working standard containing 0.1  $\mu\text{g}/\text{mL}$ . This working standard and the dilution of the stock mercury solutions should be prepared fresh daily. Acidity of the working standard should be maintained at 0.15% nitric acid. This acid should be added to the flask as needed before the addition of the aliquot.

## 7.0 Calibration

- 7.1 Transfer 0, 0.5, 1.0, 2.0, 5.0 and 10 mL aliquots of the working mercury solution (6.7) containing 0 to 1.0  $\mu\text{g}$  of mercury to a series of 300 mL BOD bottles. Add enough distilled water to each bottle to make a total volume of 10 mL. Add 5 mL of aqua regia (6.1) and heat 2 minutes in a water bath at 95°C. Allow the sample to cool and add 50 mL distilled water and 15 mL of  $\text{KMnO}_4$  solution (6.5) to each bottle and return to the water bath for 30 minutes. Cool and add 6 mL of sodium chloride-hydroxylamine sulfate solution (6.4) to reduce the excess permanganate. Add 50 mL of distilled water. Treating each bottle individually, add 5 mL of stannous sulfate solution (6.3) and immediately attach the bottle to the aeration apparatus. At this point, the sample is allowed to stand quietly without manual agitation. The circulating pump, which has previously been adjusted to rate of 1 liter per minute, is allowed to run continuously. The absorbance, as exhibited either on the spectrophotometer or the recorder, will increase and reach maximum within 30 seconds. As soon as the recorder pen levels off, approximately 1 minute, open the bypass valve and continue the aeration until the absorbance returns to its minimum value (See Note 5). Close the bypass valve, remove the fritted tubing from the BOD bottle and continue the aeration. Proceed with the standards and construct a standard curve by plotting peak height versus micrograms of mercury.

NOTE 5: Because of the toxic nature of mercury vapor precaution must be taken to avoid its inhalation. Therefore, a bypass has been included in the system to either vent the mercury vapor into an exhaust hood or pass the vapor through some absorbing media, such as:

- a) equal volumes of 0.1 N  $\text{KmnO}_4$  and 10%  $\text{H}_2\text{SO}_4$
- b) 0.25% iodine in a 3% KI solution.

A specially treated charcoal that will absorb mercury vapor is also available from Barnebey and Cheney, E. 8th Ave., and North Cassidy St., Columbus, Ohio 43219, Cat. #580-13 or #580-22.

## 8.0 Procedure

- 8.1 Weigh triplicate 0.2 g portions of dry sample and place in bottom of a BOD bottle. Add 5 mL of distilled water and 5 mL of aqua regia (6.1). Heat 2 minutes in a water bath at 95°C. Cool, add 50 mL distilled water and 15 mL potassium permanganate solution (6. 5) to each sample bottle. Mix thoroughly

and place in the water bath for 30 minutes at 95°C. Cool and add 6 mL of sodium chloride-hydroxylamine sulfate (6.4) to reduce the excess permanganate. Add 55 mL of distilled water. Treating each bottle individually, add 5 mL of stannous sulfate (6.3) and immediately attach the bottle to the aeration apparatus. Continue as described under (7.1).

- 8.2 An alternate digestion procedure employing an autoclave may also be used. In this method 5 mL of conc. H<sub>2</sub>SO<sub>4</sub> and 2 mL of conc. HNO<sub>3</sub> are added to the 0.2 g of sample. 5 mL of saturated KMnO<sub>4</sub> solution is added and the bottle covered with a piece of aluminum foil. The samples are autoclaved at 121°C and 15 lbs. for 15 minutes. Cool, make up to a volume of 100 mL with distilled water and add 6 mL of sodium chloride hydroxylamine sulfate solution (6.4) to reduce the excess permanganate. Purge the dead air space and continue as described under (7.1).

## 9.0 Calculation

- 9.1 Measure the peak height of the unknown from the chart and read the mercury value from the standard curve.
- 9.2 Calculate the mercury concentration in the sample by the formula:

$$\mu\text{g Hg/g} = \frac{\mu\text{g Hg in the aliquot}}{\text{wt of the aliquot in gms}}$$

- 9.3 Report mercury concentrations as follows: Below 0.1 μg/gm, <0.1; between 0.1 and 1 μg/gm, to the nearest 0.01 μg; between 1 and 10 μg/gm, to nearest 0.1 μg; above 10 μg/gm, to nearest μg.

## 10.0 Precision and Accuracy

- 10.1 The following standard deviations on replicate sediment samples were recorded at the indicated levels; 0.29 μg/g ± 0.02 and 0.82 ug/g ±0.03. Recovery of mercury at these levels, added as methyl mercuric chloride, was 97% and 94%, respectively.

## Bibliography

1. Bishop, J. N., "Mercury in Sediments", Ontario Water Resources Comm., Toronto, Ontario, Canada, 1971.
2. Salma, M., private communication, EPA Cal/Nev Basin Office, Alameda, California.