

volume change of the sample.

4.0 Instrument Parameters

- 4.1 Drying time and temperature: 30 sec @ 125°C
- 4.2 Charring time and temperature: 30 sec @ 1200°C
- 4.3 Atomizing time and temperature: 10 sec @ 2700°C
- 4.4 Purge Gas Atmosphere: Argon
- 4.5 Wavelength: 196.0 nm.
- 4.6 Other operating parameters should be set as specified by the particular instrument manufacturer.

5.0 Analysis Procedure

- 5.1 For the analysis procedure and the calculation see "Furnace Procedure" part 9.3 of the Atomic Absorption Methods section of this manual.

6.0 Notes

- 6.1 The above concentration values and instrument conditions are for a Perkin-Elmer HGA- 2100, based on the use of a 20 μL injection, purge gas interrupt and non pyrolytic graphite. Smaller size furnace devices or those employing faster rates of atomization can be operated using lower atomization temperatures for shorter time periods than the above recommended settings.
- 6.2 The use of background correction is recommended.
- 6.3 Selenium analysis suffers interference from chlorides ($> 800 \text{ mg/L}$) and sulfate ($> 200 \text{ mg/L}$). For the analysis of industrial effluents and samples with concentrations of sulfate from 200 to 2000 mg/L, both samples and standards should be prepared to contain 1% nickel.
- 6.4 For every sample matrix analyzed, verification is necessary to determine that method of standard addition is not required (see part 5.2.1 of the Atomic Absorption Methods section of this manual).
- 6.5 For quality control requirements and optional recommendations for use in drinking water analyses, see part 10 of the Atomic Absorption Methods section of this manual.
- 6.6 If method of standard addition is required, follow the procedure given earlier in part 8.5 of the Atomic Absorption Methods section of this manual.
- 6.7 Data to entered into STORET must be reported as $\mu\text{g/L}$.

7.0 Precision and Accuracy

- 7.1 Using a sewage treatment plant effluent containing $< 2 \mu\text{g/L}$ and spiked with a concentration of $20 \mu\text{g/L}$, a recovery of 99% was obtained.
- 7.2 Using a series of industrial waste effluents spiked at a $50 \mu\text{g/L}$ level, recoveries ranged from 94 to 112%.
- 7.3 Using a 0.1% nickel nitrate solution as a synthetic matrix with selenium concentrations of 5, 10, 20, 40, 50, and $100 \mu\text{g/L}$, relative standard deviations of 14.2, 11.6, 9.3, 7.2, 6.4 and 4.1%, respectively, were obtained at the 95% confidence level.
- 7.4 In a single laboratory (EMSL), using Cincinnati, Ohio tap water spiked at

concentrations of 5, 10, and 20 $\mu\text{g Se/L}$, the standard deviations were ± 0.6 , ± 0.4 , and ± 0.5 , respectively. Recoveries at these levels were 92%, 98%, and 100%, respectively.

8.0 Reference

- 8.1 "Determining Selenium in Water, Wastewater, Sediment and Sludge By Flameless Atomic Absorption Spectroscopy", Martin, T.D., Kopp, J. F. and Ediger, R.D. Atomic Absorption Newsletter 14, 109 (1975).