Installation and Operating Instructions

Model SM5-MC

Sound Monitor 418MHz RF Transmitter

Transmitter Features:

SM5-MC transmitters are battery operated. SM5-MC is coded at the manufacturer with a unique address ID. Transmitter has a test button for programming and follow on operational testing. When activated the RED LED will illuminate and the SM5-MC will transmit a sound alert to the Silent Call Medallion series receiver. Low Battery YELLOW LED will flash every 10 seconds when batteries are in need of replacement.

Transmitter Mounting and Location:


The SM5-MC Medallion Series Transmitter is designed for indoor use only. The transmitter case has provisions for mounting to a nail or screw on the wall. Velcro strip is also provided for mounting on desktop or wall. Avoid mounting transmitter near large metal objects, other transmitting devices in close proximity or high current appliances that will degrade the transmitter / receiver operation. Test transmitter mounting location before mounting to ensure proper system operation. Adjust antenna for best results. Transmission range with a Silent Call Medallion Series receiver is approximately 500 feet

Installation:

1. Remove your Sound Monitor from the packaging.
2. Place the Sound Monitor as close to the sound source as possible that you want to monitor (but not within a baby's reach). For example, place it on a dresser next to a baby's bed or next to a telephone, alarm clock or apartment doorbell speaker.
3. The Sound Monitor can be mounted to any flat surface with the supplied Velcro strips or you can hang it on a wall using the access hole on the back of the case. This step is not necessary if you will be moving the Sound Monitor from location to location.
4. Avoid placing the Sound Monitor on a metallic surface. Metal will interfere with the signal transmission and shorten the signal range.

Operation:

1. The Sound Monitor operates on 2 AA Alkaline batteries (included).
2. There is a Yellow LED (low battery indicator light) on the face of the transmitter. When this light comes on you know that the batteries need replacing.
3. Slide the on/off switch on the side of the Sound Monitor to the ON position.
4. Any continuous loud sound will cause the Sound Monitor to activate and the Red LED on the face of the transmitter will come on. The transmitter then sends a wireless signal, to any Silent Call Medallion Series Receiver activating the receiver.
5. Transmission range with a Medallion Series Receiver is about 500 feet.
6. The Sound Monitor is designed so that brief sounds such as door slams, regardless of how loud they are, will not activate the Sound Monitor.
7. The Sound Monitor leaves the factory set to most sensitive position, which should be adequate for almost all situations. If for some reason the sensitivity needs to be changed, there is an adjustment dial located behind the battery door. Rotating the adjustment dial clockwise decreases the sensitivity.
8. A new feature of the Sound Monitor is the ability to set the Sound Monitor to detect a specific sound. Remove the battery cover and at the bottom left hand corner there is a 3 position switch. Switch # 1 ON is for Sound Alert, Switch #2 ON is for Telephone Alert and Switch # 3 ON is for Doorbell Alert.

3 position switch

Replacement Batteries:

Use two AA 1.5V Alkaline Batteries. AA Batteries may be purchased from local retail stores. Transmitter is shipped with 2 AA Alkaline Batteries. Batteries may also be purchased Online at www.silentcall.com.

Installation and Operating Instructions

Model SM5-MC

Sound Monitor 418MHz RF Transmitter

Regulatory Information Notice:

This device complies with Industry Canada license-exempt RSS standard(s).

Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation of the device. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- *Reorient or relocate the receiving antenna.*
- *Increase the separation between the equipment and the receiver*
- *Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.*
- *Consult the dealer or an experienced radio/television technician for help*

Unauthorized changes or modifications could void the user's authority to operate the equipment.

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Technical Support

For technical support on this or any other Silent Call product, please feel free to contact us. You can reach us by phone at 800-572-5227 (voice or TTY) or by Email at support@silentcall.com

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Limited Warranty

Your transmitter is warranted to be free of defects in material and workmanship for five years from the date of initial purchase. During that time, the unit will be repaired or replaced free of charge when shipped prepaid to Silent Call Communications. This warranty is void if the defect is caused by customer abuse or neglect.

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Return for Serving To:

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Silent Call Communications

5095 Williams Lake Road, Waterford Michigan 48329

(800) 572-5227 v/tty (248) 673-7360 fax

Website: www.silentcall.com Email: silentcall@silentcall.com