


VIBRATION PAD

VIB12V


Parts & Function



Setup

This vibration pad only works with the Multi-Alert system receiver (Chime Unit).

When using the vibration pad, move the Vibration On/Off switch to the "ON" position



#1 DC jack: If using the vibration pad, you must plug in the 12V power supply into the DC 12V Input first.

#2 Vibration pad on/off: Slide this switch to turn the vibration function on/off.

#3 Vibration Out: Insert the vibration pad plug into the pad out jack.

#4 12V ready indicator: This indicator will light up red, when the 12V adaptor is connected correctly.

#5 Pad on indicator: If the vibration pad On/Off switch is moved to the "ON" position, the LED indicator will light up red to indicate the Vibration pad is "ON"

Notes.

The Vibration Pad will not work without the power supply being connected.

When the vibration pad is activated, it will run for 15 seconds.

Both the 12V ready indicator & Pad On indicator will light up Red when they are in use.

To cancel any alert on the remote receiver press the reset button. The reset button is located on the back of the receiver behind the flasher.

FCC Statement:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- * Reorient or relocate the receiving antenna.
- * Increase the separation between the equipment and receiver.
- * Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: Modifications to this product will void the user's authority to operate this equipment.

For Customer Support Please Contact:

Silent Call / MyAlert

Sales@Silentcall.com / 1-800-572-5227

scdirectsales.com for DIY instructional videos and User Guides

