

## Bio Feedback Elbow Support



**Bio Feedback Elbow Support is designed to relieve a wide range of joint and muscle pain.**

### **Sprains, strains, and pulled muscles:**

Provides firm support to injured elbows without restricting movement, helping to prevent further damage.


Adjustable compression over the elbow helps reduce pain, swelling and tenderness.

### **Weak or arthritic joints:**

Provides firm support to elbow, reducing load on joint and alleviating aches, pain and stiffness.

Neoprene fabric helps retain heat to provide gentle warming relief from aches and pain.

### **How to wear your Bio Feedback Elbow Support correctly**



1. Put hand through the inner fabric sleeve
2. Position the support over elbow. Fasten and adjust all straps to desired compression level.

### **Directions:**

For sprains, strains and pulled muscles, use the support in conjunction with the following techniques to provide symptomatic relief:

P. Protection: Protect the injury from further damage by wearing the support.

R. Rest: Avoid activity that causes pain in the first 48-72 hours after injury.

I. Ice: Cool the injured area as soon as possible for 15-20 minutes every 2-3 hours to limit pain and inflammation.

C. Compression: Use the support to apply compression to the affected area, reducing swelling.

E. Elevation: Keep the injured area raised and supported to help reduce swelling.

For weak or arthritic joint, wear the support as needed to provide daytime relief from aches, pain and stiffness.

The inner fabric sleeve is provided for additional comfort. If preferred, the support can be used without the inner sleeve.

### **Hazards and Cautions:**

1. The indications outlined may not be suitable for your condition, before use always seek medical advice from your physician
2. Only use on intact skin. Do not use if you have a known allergy to any of the specified materials
3. Avoid if you have Neoprene and Nylon sensitivity
4. The product should be applied as illustrated.
5. Ensure product is not uncomfortably tight and blood flow is not restricted. Avoid wearing for prolonged periods (e.g.) whilst sleeping
6. If a rash develops, pain is prolonged or conditions worsen, discontinue use and consult a physician
7. If your symptoms persist or worsen, seek medical advice.

### **Composition of Materials**

55% neoprene, 20% nylon, 20% polyester, 5% spandex

Material composition infused with copper.

### **Care Instructions**

Machine wash cold. Delicate cycle with like colours. Tumble dry low. Do not Iron. Do not Bleach.

Made in China for Bioenergiser Ltd.

Copyright protected. Version: 1.1.15

[www.globalproductsourcing.co.uk](http://www.globalproductsourcing.co.uk) to purchase the full range of supports and other wellbeing products.