

Introduction of Stormfiber-II

Stormfiber-II is a first rate air cooling system developed for Pre, post and During the laser treatment without interfering the laser beam emission.

Stormfiber-II is designed for reducing pains, redness, swelling and thermal damages caused by laser treatment, Intense Light (IPL) treatment and Radio frequency (RF) treatment.



Stormfiber-II is designed for reducing pains, redness, swelling and thermal damages caused by laser treatment, Intense Light (IPL) treatment and Radio frequency (RF) treatment




Working Theory

Stormfiber-II Makes the best performances, pre, post and during the laser treatments of the following treatments:

- Laser hair removal by diode laser, IPL, Co2 Fractional laser, Long pulse Nd:yag laser or alexsandrite laser.
- ② Radio frequency High frequency skin rejuvenation treatments.
- ③ Temporary anesthesia before injections
- 4 Feature of Stormfiber-II system
- ⑤ Powerful cooling air temperature (-30C)







Working Theory

- 6 High power compressor makes the system work for long time continuously
- Thigh speed fans makes air emitting strong for better cooling result, 6 degree cooling fan speeds for various treatment.
- Self defrosting system for best cooling performance
- Self water drainage system make safe and cleaning performance
- Fuse supporter make the procedure comfortable







Specifica tions

Laser source:	Stormfiber-II air cooling system	
Controlling screen:	2.4 inch multi color screen	
Totol power:	750W	
Refrigeration Power: Temperature range:	1500w -4°C~-30°C(adjustable)(-4/-8/-12/-16/- 20/-24/-28/-30)	
Fan speed::	1~8 (adjustable)	
Power requirement:	110V,50hz,6A, 220V,50hz, 3A	
Net weight:	50kg	
Package Weight:	74kg	
Treatment Tube:	2.5 meter	
Power Source:	100V~240V/ 50-60Hz	
Fuse supporter:	3 joint supporter	
Dimension of system:	610*380*750mm	
Package size:	820*500*930mm	
Safety Class:	EN60601-1, Type B MDD/MPG,Class I	





Appearance introduction

super-low temperature refrigeration APC003-750W



Compressors imported from Japan



Super wind speed blower



Independent noise reduction system



-30°C super-low temperature refrigeration



200g super light cold air tube



Stepless speed regulation, temperature adjustment knob



One-key startstop button

Operator station introduction



- ① Stepless speed regulation, temperature adjustment knob
 - ② Temperature indication $(-4^{\circ}C\sim-30^{\circ}C)$
 - ③ Strat-Stop button



01

02

STORMFIBER-II

O1 HIGHLY Compressor Imported from Japan

DUPONT Environment protecting refrigerant

System interface Introduction



Wind speed: OFFLevel 0


temperature display: -4°C~-30°C Adjustable

Wind speed level: (-4/-8/-12/-16/-20/-24/-28/-30)

Low temperature test

Test Conditions:

- 1. Ten minutes after starting up
- 2. Test Tool-Taiwan Weather Meter



Noise test

Test Conditions:

- 1. Ten minutes after starting up
- 2. Test Tool—HS5633A Sound Level Meter

The voice of a person during normal communication is 40-60 dB

Before turning on



Noise index: 33.6dB

Working state



Noise index: 58dB



Compared with other companies, our advantages:

- ① HIGHLY Compressor Imported from Japan and DUPONT environmentally friendly refrigerant R404a;
- ② Ultra-low temperature refrigeration -35 degrees
- 3 Applicable to foreign machines (same size as ZIMMER and foton);;
- 4 2.5mm length and 200g weight of the air pipe are more flexible;
- ⑤ 5.1L large capacity bellows, sustainable ultra-low temperature treatment;
- Dedicated bracket for easy operation with orientation;
- 2.4-inch true color touch screen, visual human-machine interface;
- 8 700L/min super strong wind speed, 8 gears adjustable, suitable for various treatments;
- The whole machine is equipped with silencer to reduce treatment noise and make treatment more comfortable;
- Super strong platform, placed with a variety of instruments for use.

Best Partner









Laser hair removal

 Reduce skin pain and the risk of burns during hair removal

PICO laser

 When removing freckles, moles and tattoos, it can effectively reduce skin pain and reduce the risk of blistering

CO2 Fractional Laser

 When removing scars, moles, and skin rejuvenation, it can accelerate the heat dissipation of the skin and reduce pain

Handle PICO laser

 When removing freckles, moles and tattoos, it can effectively reduce skin pain and reduce the risk of blistering

Name	Round Adapte	Middle square Adapter	Large square Adapter	Diode laser Adapter	SHR IPL Adapter
Image					
Function Descriptio n	For reducing the skin temperature of small treatment area, such as eyebrow, underarm, forehead.	For reducing the skin temperature of middle square area, especially for hair removal treatment, such Arm, under arm, leg	For reducing the skin temperature of large treatment area , such as thigh, belly. especially for hair removal treatment.	Cooperating with P-808 Diode laser system for hair removal treatment	Cooperating with SHR IPL series for hair removal and skin rejuvenation treatment.

808 connector

SHR connector





Best partner

After inserted the connector assembly, The I-Cool could works with 808nm diode laser hair removal machine, co2 fractional laser machine, SHR machine, IPL machine and Nd yag laser machine... Suitable for all machines if patients feel pain and needs cooling procedure.

	Stormfiber-II	Contact Cooling	Cryogen Spray	ICE Pack
Cooling Temperature	-20°C	-5°C	-3°C ~ -10°C	0°C ~ -10°C
Cooling Durability	Long	Mid	Short	Mid~short
Availability to use pre, post and during treatment	High	Less	Poor	Less
Patient's comfortability	High	Less	Poor	Less
Cooling efficacy	High	Less	Poor	Less
Operator's preparation	Simple	Simple	Simple	Complicated
Clean degree	No contact	Contact	No contact	Contact

THANKS

Welcome to consult more!

www.apexxmed.con

Address: 2219 SW. 59th Terrace, west park fl 33023

Phone: 954-756-3325

Email: sales@apexxmed.com

