

For Research Use Only



INSTRUCTION

MANUAL

rYeast Total RNA Mini Kit

IB47410 (4 Preparation Sample Kit)

IB47411 (100 Preparation Kit)

IB47412 (300 Preparation Kit)

Advantages

Sample:	a variety of yeast and other fungus species
Yield:	up to 30 µg of RNA (5 x 10 ⁷ <i>Saccharomyces cerevisiae</i> : 20 µg)
Convenient:	includes Sorbitol Buffer to reduce sample preparation time
Format:	certified DNase and RNase-free spin columns
Time:	within 20 minutes
Elution Volume:	50-100 µl
Kit Storage:	dry at room temperature (15-25°C)

Introduction

The rYeast Total RNA Mini Kit was designed for total RNA purification from yeast and a wide variety of other fungus species. Sorbitol Buffer is included with the kit to reduce sample preparation time and minimize hands on time. Detergents and chaotropic salt are used to lyse cells and inactivate RNase while RNA is bound by the glass fiber matrix of the RNA spin column. Once any contaminants have been removed, using the Wash Buffer (containing ethanol), the purified total RNA is eluted by RNase-free Water. High quality total RNA can be purified within 20 minutes and is ready for use in RT-PCR, Northern Blotting, Primer Extension, mRNA Selection and cDNA Synthesis.

Quality Control

The quality of the rYeast Total RNA Mini Kit is tested on a lot-to-lot basis by isolating RNA from *Saccharomyces cerevisiae* (5×10⁷) harvested by centrifugation at 5,000 x g for 10 minutes. A 5 µl aliquot of purified RNA from a 50 µl eluate is analyzed by electrophoresis on a 0.8% agarose gel.

Kit Contents

Component	IB47410	IB47411	IB47412
Sorbitol Buffer	4.5 ml	90 ml	225 ml
RB Buffer	2 ml	60 ml	130 ml
DNase I1 (2U/µl)	20 µl	550 µl	550 µl x 3
DNase I Reaction Buffer	200 µl	5 ml	15 ml
W1 Buffer	2 ml	50 ml	130 ml
Wash Buffer2 (Add Ethanol)	1.5 ml (6 ml)	25 ml + 12.5 ml (100 ml) (50 ml)	50 ml x 2 (200 ml) x 2
RNase-free Water	1 ml	15 ml	30 ml
RB Columns	4	100	300
2 ml Collection Tubes	8	8	600

¹ DNase I is shipped at room temperature and stored at -20°C for extended periods after receiving the kit.

² Add absolute ethanol (see the bottle label for volume) to Wash Buffer then mix by shaking for a few seconds. Check the box on the bottle. Be sure and close the bottle tightly after each use to avoid ethanol evaporation.

Steps to prevent RNase contamination


1. During operation, always wear a lab coat, disposable gloves, protective goggles and (anti-fog) procedure mask.
2. Disposable plasticware and automatic pipettes should be sterile (RNase-free) and used only for RNA procedures.
3. Non-disposable glassware or plasticware should also be sterile (RNase-free).




**IMPORTANT
BEFORE USE!**

During the procedure, always wear a lab coat, disposable gloves, and protective goggles.


Quick Protocol Diagram




1. Cell lysis of yeast and other fungus species



2. RNA binding to membrane while contaminants remain suspended



3. Wash (removal of contaminants while RNA remains bound to membrane)



4. Elution of pure total RNA which is ready for subsequent reactions

96 Well Blood gDNA Kit Protocol

Please read the entire instruction manual prior to starting the Protocol Procedure.



**IMPORTANT
BEFORE USE!**

Add absolute ethanol (see the bottle label for volume) to Wash Buffer then mix by shaking for a few seconds. Check the box on the bottle. Be sure and close the bottle tightly after each use to avoid ethanol evaporation.

DNA Removal Options: For DNA-free RNA perform either option 1 (following RNA Binding) or option 2 (following RNA Elution).

Additional Requirements: lyticase or zymolase, absolute ethanol, ddH₂O (RNase-free and DNase-free) to prepare 70% ethanol, microcentrifuge tubes (RNase-free), pipette tips (RNase-free), β -mercaptoethanol, EGTA (for DNA Digestion In Solution).

Yeast/Fungus Protocol Procedure

1. Sample Preparation

A. Yeast/Fungus on Agar Plate

Use an inoculating loop to transfer a small piece of yeast/fungus (up to 5×10^7) from an agar plate to a 1.5 ml microcentrifuge tube containing 600 μ l of Sorbitol Buffer.

B. Yeast/Fungus in Broth

Transfer yeast/fungus cells (up to 5×10^7) in broth to a 1.5 ml microcentrifuge tube. Centrifuge for 10 minutes at $5,000 \times g$ then discard the supernatant. Re-suspend the cells in 600 μ l of Sorbitol Buffer.

Add 200 U of lyticase or zymolase then mix well. Incubate at 30°C for 30 minutes. Centrifuge the mixture for 10 minutes at $2,000 \times g$ to form a spheroplast pellet. Discard the supernatant.

2. Cell Lysis

Add 300 μ l of RB Buffer and 3 μ l β -mercaptoethanol (or 6 μ l of freshly prepared 2M Dithiothreitol in RNase Free Water) to the spheroplast pellet then vortex to mix. Incubate at room temperature for 5 minutes then centrifuge at $14\text{-}16,000 \times g$ for 2 minutes. Transfer the supernatant to a new 1.5 ml microcentrifuge tube (RNase-free).

3. RNA Binding

Add 500 μ l of 70% ethanol to the lysate and pipette immediately. Place a RB Column in a 2 ml Collection Tube. Transfer 500 μ l of the mixture to the RB Column. Centrifuge at $14\text{-}16,000 \times g$ for 1 minute then discard the flow-through. Transfer the remaining mixture to the same RB Column and centrifuge at $14\text{-}16,000 \times g$ for 1 minute. Discard the flow-through and place the RB Column in a new 2 ml Collection Tube.

Optional Step 1: In Column DNase I Digestion

DNA contamination is significantly reduced following In Column DNase I Digestion. However, traces of residual DNA may be detected in very sensitive applications. In this situation, please perform Optional Step 2: DNA Digestion In Solution instead to efficiently remove trace amounts of DNA. Standard DNase buffers are incompatible with In Column DNase I Digestion and may affect RNA integrity and reduce yield.

1. Add 400 μ l of Wash Buffer (make sure ethanol was added) to the RB Column then centrifuge at $14\text{-}16,000 \times g$ for 30 seconds.
2. Discard the flow-through and place the RB Column back in the 2 ml Collection Tube.
3. Prepare DNase I solution in a 1.5 ml microcentrifuge tube (RNase-free) as follows:

DNase I	5 μ l (2 U/ μ l)
DNase I Reaction Buffer	45 μ l
Total Volume	50 μ l

4. Gently pipette DNase I solution to mix (DO NOT vortex) then add DNase I solution (50 μ l) into the CENTER of the RB column matrix.
5. Incubate the column for 15 minutes at room temperature ($20\text{-}30^\circ\text{C}$) then proceed with the RNA Wash step.

4. RNA Wash

Add 400 μ l of W1 Buffer to the RB Column then centrifuge at $14\text{-}16,000 \times g$ for 1 minute. Discard the flow-through and place the RB Column back in the 2 ml Collection Tube. Add 600 μ l of Wash Buffer (make sure ethanol was added) into the RB Column. Centrifuge at $14\text{-}16,000 \times g$ for 1 minute. Discard the flow-through then place the RB Column back in the 2 ml Collection Tube. Add 600 μ l of Wash Buffer (make sure ethanol was added) into the RB Column. Centrifuge at $14\text{-}16,000 \times g$ for 1 minute. Discard the flow-through then place the RB Column back in the 2 ml Collection Tube. Centrifuge at $14\text{-}16,000 \times g$ for 3 minutes to dry the column matrix.

5. RNA Elution

Place the dried RB Column in a clean 1.5 ml microcentrifuge tube (RNase-free). Add 50 μ l of RNase-free Water into the CENTER of the column matrix. Let stand for at least 3 minutes to ensure the RNase-free Water is absorbed by the matrix. Centrifuge at $14\text{-}16,000 \times g$ for 1 minute to elute the purified RNA.

Optional Step 2: DNA Digestion In Solution

1. Prepare DNase I reaction in a 1.5 ml microcentrifuge tube (RNase-free) as follows:

RNA in RNase-free Water	1-40 μ l
DNase I	0.5 μ l/ μ g RNA
DNase I Reaction Buffer	5 μ l
RNase-free Water	Add to final volume = 50 μ l
Total Volume	50 μ l

2. Gently pipette the DNase I reaction solution to mix (DO NOT vortex) then incubate the microcentrifuge tube at 37°C for 15-30 minutes.

3. Stop the reaction by adding 1 μ l of 20 mM EGTA (pH=8.0) then incubate the microcentrifuge tube at 65°C for 10 minutes.

NOTE: DNase I Reaction Buffer may cause aberrant migration or smearing of RNA on gels. If analyzing RNA by gel electrophoresis, repurify the RNA sample by using the RNA Pure Kit instead of stopping the reaction with EGTA.

rYeast Total RNA Mini Kit Functional Test Data




Figure 1.

Total RNA was extracted using the rYeast Total RNA Mini Kit.

Saccharomyces cerevisiae (5×10^7) was harvested by centrifugation at 5,000 x g for 10 minutes. A 5 μ l aliquot of purified RNA from a 50 μ l eluate was analyzed by electrophoresis on a 0.8% agarose gel.

M = 1 Kb DNA Ladder

Test	RNA Conc.	260/280	260/230	Yield
1	391.0 μ g/ml	2.19	2.48	19.6 μ g
2	2 389.9 μ g/ml	2.19	2.51	19.5 μ g
3	3 387.9 μ g/ml	2.20	2.51	19.4 μ g

Troubleshooting

Low Yield

Clogged Column.

Reduce the amount of starting material or separate it into multiple tubes. Centrifugation temperature must be between 20°C to 25°C. Yeast cells were not completely homogenized. Make sure lyticase or zymolase was added to Sorbitol Buffer immediately prior to use.

Residual Ethanol Contamination.

Following the wash step, dry the RB Column with additional centrifugation at 14-16,000 x g for 5 minutes.

RNA Degradation.

The harvested sample should be stabilized immediately prior to use. Disposable plasticware & automatic pipettes should be sterile (RNase-free) & used only for RNA procedures. Non-disposable glassware or plasticware should be sterile (RNase-free).

Related DNA/RNA Extraction Products

Total RNA Extraction		
Product	Package Size	Catalogue Number
Total RNA Mini Kit (Blood/Cell)	50/100/300 preps	IB47321/322/323
Total RNA Mini Kit (Tissue)	50/100/300 preps	IB47301/302/303
Total RNA Mini Kit (Plant)	50/100/300 preps	IB47341/342/343
rBAC Mini RNA Bacteria Kit	100/300 preps	IB47421/422
rYeast Total RNA Mini Kit	100/300 preps	IB47411/412
miRNA Isolation Kit	100 preps	IB47371
IBI Isolate	100/200 rxns	IB47601/302
Plasmid DNA Extraction		
Product	Package Size	Catalogue Number
I-Blue Mini Plasmid Kit	100/300 preps	IB47171/172
I-Blue Mini Plasmid Endo Free Kit	100 preps	IB47176
I-Blue Midi Plasmid Kit	25 preps	IB47181
I-Blue Midi Plasmid Kit	25 preps	IB47191
Midi Fast Ion Plasmid Kit	25 preps	IB47111
Midi Fast Ion Plasmid Kit (Endotoxin Free)	25 preps	IB47113
Maxi Fast Ion Plasmid Kit	10/25 preps	IB47121/122
Maxi Fast Ion Plasmid Kit (Endotoxin Free)	10/25 preps	IB47124/125
96 Well Plasmid Kit	4/10 x 96 preps	IB47151/152
Post Reaction DNA Extraction		
Product	Package Size	Catalogue Number
Gel/PCR DNA Fragments Extraction Kit	100/300 preps	IB47020/030
Small DNA Fragments Extraction Kit	100/300 preps	IB47061/062
96 Well PCR Cleanup Kit	4/10 x 96 preps	IB47040/050
Genomic DNA Extraction		
Product	Package Size	Catalogue Number
Genomic DNA Mini Kit (Blood/Cultured Cell)	100/300 preps	IB47201/202
Genomic DNA Mini Kit (Tissue)	50/300 preps	IB47221/222
gMAX Mini Kit (Blood/Tissue)	100/300 preps	IB47281/282
Genomic DNA Mini Kit (Plant)	100 preps	IB47230
gSWAB Mini Genomic DNA Kit	100/300 preps	IB47276/277
gBAC Mini DNA Bacteria Kit	100/300 preps	IB47291/292
gYEAST Genomic DNA Kit	100/300 preps	IB47266/267
96 Well Blood Genomic DNA Extraction Kit	4/10 x 96 preps	IB47251/252
gPURE Cell DNA Isolation Kit	100/1000 rxns	IB47431/432
DNA/RNA/Protein Extraction		
Product	Package Size	Catalogue Number
DNA/RNA/Protein Extraction Kit	50/100 preps	IB47701/702

We have more to offer!

Check out our wide variety of products to meet all of your research needs

www.ibisci.com





SCIENTIFIC

7445 Chavenelle Road • Dubuque, IA 52002
800-253-4942 • (563) 690-0484 • info@ibisci.com

IBISCI.com

