# RAPID SARS-COV-2 ANTIGEN TEST CARD

## INSTRUCTION GUIDE FOR ANTERIOR NASAL SWAB SPECIMENS

For private use/home use/self-testing

REF	IN40C5-2	For 1 Test/Box
REF	IN40C5-4	For 5 Tests/Box
REF	IN40C5-6	For 20 Tests/Box

Please follow the instruction leaflet carefully.

**INTENDED USE**  
Rapid SARS-CoV-2 Antigen Test Card is an immunochromatography based one step in vitro test. It is designed for the rapid qualitative determination of SARS-CoV-2 virus antigen in anterior nasal swabs from individuals suspected of COVID-19 within the first seven days of symptom onset. Rapid SARS-CoV-2 Antigen Test Card shall not be used as sole basis to diagnose or exclude SARS-CoV-2 infection. Children under 14 years of age should be assisted by adults.

**SUMMARY**  
The novel coronavirus belongs to the  $\beta$  genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients caused by the novel coronavirus are the main source of infection, asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

**MATERIALS PROVIDED**

Components	For 1 Test/Box	For 5 Tests/Box	For 20 Tests/Box
Rapid SARS-CoV-2 Antigen Test Card (sealed foil pouch)	1	5	20
Sterile swab	1	5	20
Extraction tube	1	5	20
Sample extraction buffer	1	5	20
Instructions for use (this leaflet)	1	1	1
Tube stand	1 (packaging)	1	1

**PERFORMANCE (SENSITIVITY AND SPECIFICITY)**

Rapid SARS-CoV-2 Antigen Test Card was compared to the confirmed clinical diagnosis. The study involved 156 samples.

Sensitivity	97.72%
Specificity	99.20%
Accuracy	98.72%

A feasibility study demonstrated that:  
- 99.10% of non-professionals carried out the test without requiring assistance  
- 97.87% of the different types of results were interpreted correctly

**INTERFERENCES**

None of the following substances at the tested concentration showed any interference with the test:

Whole Blood: 1%	Alkaloi: 10%	Mucin: 2%
Phenylephrin: 0.15%	Trombramycin: 0.0004%	Oxymetazolin: 15%
Menthol: 0.15%	Cromolyn: 15%	Benzocaine: 0.15%
Fluticasone Propionate: 5%	Mupirocin: 0.25%	Zicam Nasal Spray: 5%
Osseltamivir Phosphate: 0.5%	sodium chloride: 5%	Human Anti-mouse Antibody (HAMA): 60 ng/mL
Biotin: 1200 ng/mL		

**IMPORTANT INFORMATION BEFORE THE EXECUTION**

- Read this instruction guide carefully.
- Do not use the product beyond the expiration date.
- Do not use the product if the pouch is damaged or the seal is broken.
- Store the test device at 4 to 30 °C in the original sealed pouch. Do Not Freeze.
- The product should be used at room temperature (15 °C to 30 °C). If the product has been stored in a cool area (less than 15 °C), leave it at normal room temperature for 30 minutes before using.
- Handle all specimens as potentially infectious.
- Inadequate or inappropriate specimen collection, storage, and transport may yield inaccurate test results.
- Use the swabs included in the test kit to ensure optimal performance of the test.
- Proper specimen collection is the most important step in the procedure. Make sure to collect enough specimen material (nasal secretion) with the swab, especially for anterior nasal sampling.
- Blow the nose several times before collecting specimen.
- The specimens should be tested as soon as possible after collection.
- Apply the drops of test specimen well (S).
- Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.
- Children under 14 years of age should be assisted by an adult.
- LIMITATIONS

- The test is to be used exclusively for the qualitative detection of SARS-CoV-2 viral antigen in anterior nasal swab specimens. The exact concentration of SARS-CoV-2 viral antigen cannot be determined as part of this test.
- Proper specimen collection is critical. Failure to follow the procedure may result in inaccurate test results. Improper collection, storage, or even freezing and thawing of the specimen can lead to inaccurate test results.
- If the viral load of the specimen is below the detection limit of the test, the test may produce a negative result.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be made by the physician after evaluation of all clinical and laboratory results.
- A negative result does not exclude viral infection except for SARS-CoV-2 and should be confirmed by molecular diagnostic methods if COVID-19 is suspected.
- A positive result does not exclude coinfection with other pathogens.
- The SARS-CoV-2 rapid antigen test can detect both viable and non-viable SARS-CoV-2 material. The performance of the SARS-CoV-2 rapid test is dependent on viral load and may not correlate with other diagnostic methods performed on the same specimen.
- Users should test specimens as soon as possible after specimen collection and within two hours of specimen collection.
- Sensitivity for nasal or oropharyngeal swabs may be lower than nasopharyngeal swabs. It is recommended to use the nasopharyngeal swab specimens by healthcare professionals.
- Monoclonal antibodies may fail to detect, or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope region.
- The amount of antigen in a sample may decrease as the duration of illness increases. Specimens collected after day 5-7 of illness are more likely to be tested negative compared to a RT-PCR assay.
- The kit was validated with the assorted swabs. Use of alternative swabs may result in false negative results.
- The validity of Rapid SARS-CoV-2 Antigen Test Card has not been proven for identification/confirmation of tissue culture isolated and should not be used in this capacity.
- Cross-reactivity of the Test Device was evaluated by testing viruses and other microorganisms. The final test concentrations of viruses and other microorganisms are documented in the Cross-Reactivity-Study. The therein following viruses and other microorganisms except the Human SARS-coronavirus have no effect on the test results of the Test Device. Positive test results do not rule out co-infections with other pathogens. Positive results may occur in cases of infection with SARS-CoV.

**PREPARATION**






- Clear, clean and dry a flat surface.
- Check the test kit contents. Make sure that nothing is damaged or broken.
- Minut at hand.
- Blow your nose several times before collecting specimen.
- Wash hands.

**DISPOSAL**  
The test kit may be disposed of with normal household waste in accordance with the applicable local regulations.

**PROCEDURE**

This test is suitable for people of all ages. The recommended operators are aging from 14 to 90. Children under 14 years of age should be tested by an adult. Do not continue the test if the child feels any pain.

Rotate the lid of sample extraction buffer bottle.  
**Caution: Open it away from your face and be careful not to spill any of the liquid.**

**INTERPRETATION OF RESULTS**  
**Positive:**  
If two colored bands appear with one colored band in the Control Zone (C) and another in the Test Zone (T) within 15-20 minutes, the test result is positive.  
**Caution: No matter how faint the colored band is in the Test Zone (T), the result should be considered as positive.**  
**Negative:**  
If one colored band appears in the Control Zone (C) and no colored band appears in the Test Zone (T) within 15-20 minutes, the test result is negative.  
**Invalid:**  
If no color line appears in the control area (C) within 15-20 minutes, the test is invalid. Repeat the test with a new test card.

**QUALITY CONTROL**  
The control line is an integrated reagent and is used to control the procedure. The control line appears when the test has been performed correctly and the reagents are reactive.

**FREQUENTLY ASKED QUESTIONS (FAQ)**  
**How does the detection work?**  
The N protein of the SARS-CoV-2 virus reacts with the stripe-like coating of the test line and, if present, results in a color change, i.e. a red line appears. Therefore, if the sample does not contain any viral proteins or antigens, there will be no red test line (T).  
**When should I test myself?**  
You can test yourself whether you have symptoms or not. Studies show that earlier testing within the first 4 days of illness typically means a higher viral load, which is easier to detect. Since the test result is a snapshot valid for that point in time, testing should be repeated as recommended by local authorities.  
**What can affect my test result? What should I pay attention to?**  
Be sure to blow your nose multiple times before collecting the specimen.  
Be sure to visibly collect sample material (nasal secretions).  
Perform the test immediately after taking the sample.  
Follow the instructions for use carefully.  
Apply the drops of extraction solution only to the sample well (S).  
Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.  
**The test strip is clearly discolored or smudged? What is the reason for this?**  
The test strip is discolored or smudged, as the liquid absorption of the test strip is naturally limited. If the control line does not appear or the test strip is badly smudged or discolored, making it unreadable, please repeat the test according to the instructions.  
**I have taken the test, but I don't see a control line (C). What should I do?**  
Your test result is invalid. Observe the answer to question 4 and repeat the test according to the instructions for use.  
**I am unsure about reading the result. What should I do?**  
For the result to be positive, 2 straight horizontal lines must be clearly visible with the full width of the cassette. If you are still unsure about the results, contact the nearest health facility according to the recommendations of your local authorities.  
**My result is positive. What should I do?**  
If your result is positive and the test kit clearly indicates the control line as well as the test line, you should contact the nearest medical facility as recommended by your local authorities. Your test result may be double-checked and the authority or facility will explain the appropriate next steps.  
**My result is negative. What should I do?**  
If your result is only negative shows the control line, this may mean that you are negative or that the viral load is too low to be detected. If you experience symptoms (headache, fever, fatigue, loss of sense of smell or taste, etc.), please consult your primary care physician, or the nearest health care facility as recommended by your local authorities.  
If you are not sure, you can repeat the test.  
**How can I dispose of the test kit?**  
The test kit may be disposed of with normal household waste in accordance with the applicable local regulations.

**ACCESSORIES:**

Accessory	Manufacturer	EC-Representative	CE-Mark
Swab A	Jiangsu Changfeng Medical Industry Co., Ltd. Touqiao Town, Guangling District Yangzhou 225109 Jiangsu P.R. China	Lins Service & Consulting GmbH Obere Seegasse 34/2, 69124 Heidelberg Germany	acc. 93/42/EEC 0197
Swab B	Goodwood Medical Care Ltd. 1-2 Floor, 3-919 Yonzheng Street Jinchou District Dalian 116100 Liaoning China	CMC Medical Devices & Drugs S.L. C/ Horacio Lengua No18, CP 29006, Malaga, Spain	acc. 93/42/EEC 0197

**EXPLANATION FOR SYMBOLS**

	In Vitro Diagnostics Use		See instructions for Use		Expiry Date
	Tests per Kit		Keep dry		Batch Number
	Authorized Representative		Keep away from sunlight		Manufacturer
	Do not reuse		Do not use if package is damaged		Store between 4-30 °C
	CE Mark		Catalogue Number		Warning, please refer to the instruction
	H317: Warning! Liquid component may cause an allergic skin reaction.				

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REF	1N40C5-2	For 1 Test/Box
REF	1N40C5-4	For 5 Tests/Box
REF	1N40C5-6	For 20 Tests/Box

Please follow the instruction leaflet carefully.  
INTENDED USE

Rapid SARS-CoV-2 Antigen Test Card is an immunochromatography based one step in vitro test. It is designed for the rapid qualitative determination of SARS-CoV-2 virus antigen in anterior nasal swabs from individuals suspected of COVID-19 within the first seven days of symptom onset. Rapid SARS-CoV-2 Antigen Test Card shall not be used as sole basis to diagnose or exclude SARS-CoV-2 infection. Children under 14 years of age should be assisted by a caregiver.  
SUMMARY

The novel coronavirus belongs to the  $\beta$  genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients caused by the novel coronavirus are the main source of infection, asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

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### PERFORMANCES (SENSITIVITY AND SPECIFICITY)

Rapid SARS-CoV-2 Antigen Test Card was compared to the confirmed clinical diagnosis. The Study involved 156 samples.

Sensitivity	99.77%
Specificity	99.20%
Accuracy	98.72%

A feasibility study demonstrated that:

- 99.10% of non-professionals carried out the test without requiring assistance
- 97.87% of the different types of results were interpreted correctly

### INTERFERENCES

None of the following substances at the tested concentration showed any interference with the test.

Whole Blood: 1%	Alkaloi: 10%	Mucin: 2%
Phenylephrin: 0.15%	Oxymetazolin: 0,0004%	Oxymetazolin: 15%
Menthol: 0,15%	Cromolyn: 15%	Benzocain: 0,15%
Fluticason Propionate: 5%	Mupirocin: 0,25%	Zicam Nasal Spray: 5%
Osceltamivir Phosphate: 0,5%	sodium chloride: 5%	Human Anti-mouse Antibody (HAMA): 60 ng/mL
Biotin: 1200 ng/mL		

### IMPORTANT INFORMATION BEFORE THE EXAMINATION

- Read this instruction guide carefully.
- Do not use the product beyond the expiration date.
- Do not use the product if the pouch is damaged or the seal is broken.
- Store the test device at 4 to 30 °C in the original sealed pouch. Do Not Freeze.
- The product should be used at room temperature (15°C to 30°C). If the product has been stored in a cool area (less than 15°C), leave it at normal room temperature for 30 minutes before using.
- Handle all specimens as potentially infectious.
- Inadequate or inappropriate specimen collection, storage, and transport may yield inaccurate test results.
- Use the swabs included in the test kit to ensure optimal performance of the test.
- Correct specimen collection is the most important step in the procedure. Make sure to collect enough specimen material (nasal secretion) with the swab, especially for anterior nasal sampling.
- Blow the nose several times before collecting specimen.
- The specimens should be tested as soon as possible after collection.
- Apply the drops of test specimen only to the specimen well (S).
- Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.
- Children under 14 years of age should be assisted by an adult.

### LIMITATIONS

- The test is to be used exclusively for the qualitative detection of SARS-CoV-2 viral antigen in anterior nasal swab specimens. The exact concentration of SARS-CoV-2 viral antigen cannot be determined as part of this test.
- Proper specimen collection is critical. Failure to follow the procedure may result in inaccurate test results. Improper collection, storage, or even freezing and thawing of the specimen can lead to inaccurate test results.
- If the viral load of the specimen is below the detection limit of the test, the test may produce a negative result.
- As with all diagnostic tests, a definitive clinical diagnosis should not be based on the result of a single test, but should be made by the physician after evaluation of all clinical and laboratory results.
- A negative result does not exclude viral infection except for SARS-CoV-2 and should be confirmed by molecular diagnostic methods if COVID-19 is suspected.
- A positive result does not exclude coinfection with other pathogens.
- The SARS-CoV-2 rapid antigen test can detect both viable and non-viable SARS-CoV-2 material. The performance of the SARS-CoV-2 rapid test is dependent on viral load and may not correlate with other diagnostic methods performed on the same specimen.
- Users should test specimens as soon as possible after specimen collection and within two hours of specimen collection.
- Sensitivity for nasal or oropharyngeal swabs may be lower than nasopharyngeal swabs. It is recommended to use the nasopharyngeal swab specimens by healthcare professionals.
- Monoclonal antibodies may fail to detect, or detect with less sensitivity, SARS-CoV-2 viruses that have undergone minor amino acid changes in the target epitope region.
- The amount of antigen in a sample may decrease as the duration of illness increases. Specimens collected after day 5-7 of illness are more likely to be tested negative compared to an RT-PCR assay.
- The kit was validated with the assorted swabs. Use of alternative swabs may result in false negative results.
- The validity of Rapid SARS-CoV-2 Antigen Test Card has not been proven for identification/confirmation of tissue culture isolates and should not be used in this capacity.
- Cross-reactivity of the Test Device was evaluated by testing viruses and other microorganisms. The final test concentrations of viruses and other microorganisms are documented in the Cross-Reactivity-Study. The therein following viruses and other microorganisms except the Human SARS-coronavirus have no effect on the test results of the Test Device. Positive test results do not rule out co-infections with other pathogens. Positive results may occur in cases of infection with SARS-CoV.

### PREPARATION

- Clear, clean and dry a flat surface.
- Check the test kit contents. Make sure that nothing is damaged or broken.
- Wash your hands.
- Blow your nose several times before collecting specimen.
- Wash hands.

**DISPOSAL.** The test may be disposed of with normal household waste in accordance with the applicable local regulations.

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- 
- 
- 
- 
- 
- 
- 
- 
- 
- 
- 

### INTERPRETATION OF RESULTS

**Positive:**  
If two colored bands appear with one colored band in the Control Zone (C) and another in the Test Zone (T) within 15-20 minutes, the test result is positive.  
**Caution: No matter how faint the colored band is in the Test Zone (T), the result should be considered as positive.**

**Negative:**  
If one colored band appears in the Control Zone (C) and no colored band appears in the Test Zone (T) within 15-20 minutes, the test result is negative.

**Invalid:**  
If no color line appears in the control area (C) within 15-20 minutes, the test is invalid. Repeat the test with a new test card.

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### FREQUENTLY ASKED QUESTIONS (FAQ)

- How does the detection work?**  
The N protein of the SARS-CoV-2 virus reacts with the stripe-like coating of the test line and, if present, results in a color change, i.e. a red line appears. Therefore, if the sample does not contain any viral proteins or antigens, there will be no red test line (T).
- When should I test myself?**  
You can test yourself whether you have symptoms or not. Studies show that earlier testing within the first 4 days of illness typically means a higher viral load, which is easier to detect. Since the test result is a snapshot valid for that point in time, testing should be repeated as recommended by local authorities.
- What can affect my test result? What should I pay attention to?**  
Be sure to blow your nose multiple times before collecting the specimen.  
Be sure to visibly collect sample material (nasal secretions).  
Perform the test immediately after taking the sample.  
Follow the instructions for use carefully.  
Apply the drops of extraction solution only to the sample well (S).  
Too many or too few drops of extraction solution can lead to an invalid or incorrect test result.
- The test strip is clearly discolored or smudged? What is the reason for this?**  
The test strip is naturally limited. If the control line does not appear or the test strip is badly smudged or discolored, making it unreadable, please repeat the test according to the instructions.
- I have taken the test, but I don't see a control line (C). What should I do?**  
Your test result is invalid. Observe the answer to question 4 and repeat the test according to the instructions for use.
- I am unsure about reading the result. What should I do?**  
For the result to be positive, 2 straight horizontal lines must be clearly visible with the full width of the cassette. If you are still unsure about the results, contact the nearest health facility according to the recommendations of your local authorities.
- My result is positive. What should I do?**  
If your result is positive and the test kit clearly indicates the control line as well as the test line, you should contact the nearest medical facility as recommended by your local authorities. Your test result may be double-checked and the authority or facility will explain the appropriate next steps.
- My result is negative. What should I do?**  
If the test kit only shows the control line, this may mean that you are negative or that the viral load is too low to be detected. If you experience symptoms (headache, fever, fatigue, loss of sense of smell or taste, etc.), please consult your primary care physician, or the nearest health care facility as recommended by your local authorities.  
If you are not sure, you can repeat the test.
- How can I dispose of the test card?**  
The test kit may be disposed of with normal household waste in accordance with the applicable local regulations.

### ACCESSORIES:

Accessory	Manufacturer	EC-Representative	CE-Mark
Swab A	Jiangsu Changfeng Medical Industry Co., Ltd. Touqiao Town, Guangling District Yangzhou 225109 Jiangsu P.R. China	Lins Service & Consulting GmbH Obere Seegrasse 34/2,69124 Heidelberg Germany	acc. 93/42/EEC 0197
Swab B	Goodwood Medical Care Ltd. 1-2 Floor, 3-919 Yonzheng Street Jinzhou District Dalian 116100 Liaoning China	C/MC Medical Devices & Drugs S.L C/ Horacio Lengo No18, CP 29006, Malaga, Spain	0197 acc. 93/42/EEC

**EXPLANATION FOR SYMBOLS**

IVD	In Vitro Diagnostics Use	See Instructions for Use	Expiry Date
Tests per Kit	Keep dry	LOT	Batch Number
Authorized Representative	Keep away from sunlight	Manufacturer	
Do not reuse	Do not use if package is damaged	Store between 4-30 °C	
CE Mark	CE Mark	REF	Catalogue Number
Warning	H317: Warning! Liquid component may cause an allergic skin reaction.		

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